# medicalprint®

## shell

3D PRINTING MATERIAL

MED RESIN  
biocompatible

CLASS IIA

### ITE Shells

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**INTENDED USE:**

Plastic material for the production of ITE shells in 3D printing processes

**INDICATION**

Hearing aid shells for connection to a hearing system, foil earmolds and in-ear monitoring

**PATIENT TARGET GROUP**

Persons for whom an earmold is to be created.

**INTENDED USERS**

Earmold laboratory, ENT specialist, hearing aid acoustician

**SUITABLE FOR THE FOLLOWING DLP-PRINTER**

Asiga 385 nm	MiiCraft 385 nm	Rapidshape 385 nm	W2P 385 nm
MAX UV Pro2 UV Pico2 UV Pico2 HD UV	MiiCraft Series	HA Series	Solflex Series

Printers may only be operated using material parameters authorised by DETAX!

**PROCESSING**

- The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (observe process description on page 8).
- After storage, the material in the bottle should be shaken intensively and homogenized with a bottle roller before use.
- Maximum curing depth upon direct post-exposure: blue-opaque: 1.5 mm, red-opaque: 1.2 mm, black: 1.2 mm, white: 2 mm, skin 2.0: 1.3 mm, beige 2.0: 1.3 mm

\*In case of solid earmolds and exposure on both sides, the material thickness can be up to 3 mm

(Example: **medicalprint® shell** blue-opaque – with a curing depth of 1.5 mm).

► Lacquer surface or polish surface mechanically.

► Processing temperature  $23^{\circ}\text{C} \pm 2^{\circ}\text{C}$ .

**SAFETY INFORMATION**

- Only for the specified intended use by trained specialists.
- Avoid direct contact with the liquid material and the components before post-curing, in particular in pregnant / breastfeeding women. Irritating to eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- After contact with eyes rinse thoroughly with water immediately and consult a doctor.

► After contact with skin wash immediately with water and soap.

► Biocompatibility is only guaranteed with complete polymerisation.

► The processing instructions and precautionary measures must be observed. Otherwise, irreparable damage to the hearing organ or eardrum may result.

► Refer to the relevant safety data sheet for hazard and safety information.

**NOTES**

- DETAX shall not be held liable for any damage caused by misuse.
- Always keep container tightly sealed, immediately close the container carefully after each use.
- Read and understand the safety data sheet!

**For users and/or patients:**

Any serious incidents occurring in relation to this product should be reported immediately to [incident@detax.de](mailto:incident@detax.de) and to the competent authority of the Member State in which the user and/or patient is established.

**STORAGE**

**medicalprint® shell** is to be stored dry (at  $15^{\circ}\text{C} - 28^{\circ}\text{C}$ ) and protected from light. Minimal influence of light can already induce polymerisation.

**CONTRAINDICATION**

Contains (meth)acrylics and phosphine oxide.

Some ingredients of **medicalprint® shell** may cause allergic reactions in predisposed persons. In such cases refrain from using the product. **medicalprint® shell** only insert in completely polymerised state.

**ADVERSE EFFECTS**

Product may cause allergic reactions.

**DISPOSAL**

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

**Symbols:**

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**Storage:****Application:**

At  $23^{\circ}\text{C} \pm 2^{\circ}\text{C} / 73^{\circ}\text{F} \pm 4^{\circ}\text{F}$ ,  
 $50 \pm 5\%$  rel. humidity

# Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

## Construction process

Generation of a Print Job complying with machine and material parameters

## Post-processing

If possible, post-processing should commence immediately following the construction process. After raising the platform, a drip time of approx. 10 minutes is recommended.

## Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity ≥ 98%) for 3 min. in an ultrasonic bath.

## Cleaning

Then thoroughly clean the openings, cavities and gap areas, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

## Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity ≥ 98%) for 3 min. in an ultrasonic bath. Check the openings and additional bore holes for residues. Then blow off with compressed air.

## Post-exposure

Post-exposure is performed with the xenon photoflash unit Otoflash G171 with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.

## Surface processing

Varnishing (e.g. with luxaprint® shellac) or mechanical finishing\* (polishing) of the surface

\*Mechanical finishing of the earmould and the related work steps may possibly impair the biocompatibility of the earmould.

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GmbH & Co. KG