mould

3D PRINTING MATERIAL

MED RESIN
biocompatible

CLASS IIA

Earmoulds & ITE Shells

high initial hardness
max. scratch resistance
non brittle

DETAX

DE	Gebrauchsanweisung	3
EN	Instructions for use	6
FR	Mode d'emploi	9
ES	Modo de empleo	12
IT	Istruzioni per l'uso.....	15
DA	Brugsanvisning.....	18
EL	Οδηγίες χρήσης.....	21
HR	Upute za uporabu	24
LV	Lietošanas instrukcija.....	27
LT	Naudojimo instrukcija	30
NL	Gebruiksaanwijzing	33
NO	Bruksanvisning	36
PL	Instrukcja stosowania	39
PT	Manual de instruções.....	42
RO	Instructiuni de utilizare	45
SV	Bruksanvisning	48
SK	Návod na použitie	51
SL	Navodila za uporabo	54
CS	Návod k použití	57
HU	Használati útmutató.....	60
RU	Инструкция по применению	63
TR	Kullanım kılavuzu	66

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INTENDED USE:

Resin material for the production of earmolds and ITE shells in 3D printing processes

INDICATION

Earmoulds for ear defenders and hearing aid earmoulds for connection to a hearing system, In-Ear-Monitoring

PATIENT TARGET GROUP

Persons for whom an earmold is to be created.

INTENDED USERS

Earmold laboratory, ENT specialist, hearing aid acoustician

SUITABLE FOR THE FOLLOWING DLP-PRINTER

Asiga 385 nm	MiiCraft 385 nm	Rapidshape 385 nm	W2P 385 nm
MAX UV Pro2 UV Pico2 UV Pico2 HD UV	MiiCraft y-Serie	HA Series	Solflex Serie

Printers may only be operated using material parameters authorised by DETAX!

PROCESSING

- The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the earmoulds are completely cured (observe process description on page 8).
- After storage, the material in the bottle should be shaken intensively and homogenized with a bottle roller before use.
- Maximum curing depth* at direct post-exposure: clear:
6 mm, rose: 6 mm, rose-orange: 6 mm, light beige: 6 mm, blue: 6 mm, intensive blue: 5 mm, red: 4 mm *For large earmoulds and exposure on both sides, the material thickness can be up to 12 mm(Example luxaprint® mould clear - with a curing depth of 6 mm).
- The final colour shade is achieved after a few days.
- Processing temperature $23^{\circ}\text{C} \pm 2^{\circ}\text{C}$.
- During the first days after their fabrication don't expose the finished earmoulds to the direct sunlight.

SAFETY INFORMATION

- Only for the specified intended use by trained specialists.
- Avoid direct contact with the liquid material and the components before post-curing, in particular in pregnant / breastfeeding women. Irritating to eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.

- Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.

- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

- In case of contact with skin, wash immediately with plenty of soap and water.

- Biocompatibility is guaranteed only after complete polymerization.

- The processing instructions and precautionary measures must be observed.

- Refer to the relevant safety data sheet for hazard and safety information.

IMPORTANT NOTES

- DETAX shall not be held liable for any damage caused by misuse.

- Always keep container tightly sealed, immediately close the container carefully after each use.

- Read and understand the safety data sheet!

For users and/or patients:

Any serious incidents occurring in relation to this product should be reported immediately to incident@detax.de and to the competent authority of the Member State in which the user and/or patient is established.

STORAGE

luxaprint® mould is to be stored dry (at $15^{\circ}\text{C} - 28^{\circ}\text{C}$) and protected from light. Minimal influence of light can already induce polymerisation.

CONTRAINDICATION

Contains (meth)acrylics and phosphine oxides. Some ingredients of luxaprint® mould may cause allergic reactions in predisposed persons. In such cases refrain from using the product.

ADVERSE EFFECTS

Product may cause allergic reactions.

DISPOSAL

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Symbols:

Medical device

Storage:**Application:**

$At 23^{\circ}\text{C} \pm 2^{\circ}\text{C} / 73^{\circ}\text{F} \pm 4^{\circ}\text{F}$

Manufacturing process for ear moulds

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters, see table.

Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity ≥ 98%) for 3 min. in an ultrasonic bath.

Cleaning of the additional bore holes

Then thoroughly clean the openings as well as additional bore holes (i.e. with compressed air) and, if applicable, remove the construction components carefully from the support structure.

Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity ≥ 98%) for 3 min. in an ultrasonic bath. Prior to drying, check the openings and additional bore holes for residues.

Drying

Heat the construction components for 30 min. in an oven to approx. 40 °C to remove the solvent residues from the cleaning process.

Post-exposure

Post-exposure is performed with a xenon photoflash unit (e.g. Otoflash G171) with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.

Surface processing

Option 1

Varnish the surface with luxaprint® shellac

Option 2

Mechanical finishing of the surface*
(polishing, vibratory finishing)

*Mechanical post-processing of the earmould and the related work steps may possibly impair the biocompatibility of the ear moulds.

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GmbH & Co. KG