luxaprint® cast 2.0

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luxaprint® cast 2.0

UV-curing resin based on (meth) acrylate for the generative fabrication of casts for soft silicone earmoulds, for DLP-systems with UV / LED 385 nm

Suitable for the following DLP-Printer

Asiga 385 nm	MiiCraft 385 nm	Rapidshape 385 nm	W2P 385 nm
MAX UV Pro2 UV Pico2 UV Pico2 HD UV	MiiCraft y-Serie	HA-Serie	Solflex-Serie

Printers may only be operated using material parameters authorised by DETAX!

Processing

- Observe process description on page 2.
- After storage, the material in the bottle should be shaken intensively and homogenized with a bottle roller before use.
- Processing temperature $23^{\circ}\text{C} \pm 2^{\circ}\text{C}$.

Safety Information

- Only for the specified intended use by trained specialists.
- Contains (meth)acrylates and phosphine oxides. Some ingredients of luxaprint® 3D cast 2.0 may cause allergic reactions in predisposed persons. In such cases refrain from using the product.
- Avoid direct contact with the liquid material and the components before post-curing. Irritating to eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- In case of contact with skin, wash immediately with plenty of soap and water.
- The processing instructions and precautionary measures must be observed.
- Refer to the relevant safety data sheet for hazard and safety information.

Important notes

- DETAX shall not be held liable for any damage caused by misuse.
- Always keep container tightly sealed, immediately close the container carefully after each use.
- Product may cause allergic reactions.
- Not suitable for fabricating ear moulds.
- Read and understand the safety data sheet!

Storage

- luxaprint® 3D cast 2.0 is to be stored dry (at $15^{\circ}\text{C} - 28^{\circ}\text{C}$) and protected from light. Minimal influence of light can already induce polymerisation.

Disposal

- Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Tip: If you use the earmould silicone **detax softwear® 2.0** for filling a cast and you have applied a post-exposure with 2 x 750 flashes under inert gas conditions, the use of **Cast Separator 2.0** (separating agent) is not required.

Area of application:

Generative fabrication of casts

Ordering information:

luxaprint® 3D cast 2.0
1.000 g
green-transparent 03918

Cast Separator 2.0
500 ml 03636

detax softwear® 2.0
Automix2 System

Cartridge pack
40 Shore A

8 x 50 ml cartridges
clear-transparent 04230
rose-transparent 04231

2 x 50 ml cartridges
red opaque 04232
blue opaque 04233
green opaque 04235
yellow-fluorescent 04234
pink-fluorescent 04236
night glow 04229

Cartridge pack
60 Shore A

8 x 50 ml cartridges
clear-transparent 03906
rose-transparent 03911

2 x 50 ml cartridges
red opaque 04242
blue opaque 04238
green opaque 04240
yellow-fluorescent 04239
pink-fluorescent 04241
night glow 04237

Mixing tips rose 02770
Pack à 25 pcs.

Intra Tips white 02771
Pack à 25 pcs.

Mixing gun 02699
1 pc.

softwear coat
500 ml 03509

Processing:
at $23^{\circ}\text{C} \pm 2^{\circ}\text{C}$

Storage:

More information at
www.detax.de

Made in
Germany

06/2021
4

DETAX

GmbH & Co. KG
Carl-Zeiss-Str. 4 - 76275 Ettlingen/Germany
Telefon: 0 72 43/510-0 · Fax: 0 72 43/510-100
www.detax.de · post@detax.de

luxaprint® cast 2.0

Manufacturing process for casts

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer.

Construction process

Generation of a Print Job complying with machine and material parameters. (see table)

Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

Cleaning interior of cast

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity ≥ 98%) for 3 min. in an ultrasonic bath.

Ensure that the cast is completely filled with isopropyl alcohol.
Then air-clean with compressed air.

Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity ≥ 98%) for 3 min. in an ultrasonic bath.

Then air-clean with compressed air.

Drying

Heat the construction components for 30 min. in an oven to approx. 40°C to remove the solvent residues from the cleaning process.

Post-exposure

Post-exposure is performed with a xenon photoflash unit (e.g. Otoflash G171) with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.

Coating with Cast Separator 2.0

see instructions for Cast Separator 2.0

Filling of the cast

Then the cast can be filled with silicone, i.e. detax softwear® 2.0.

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