# luxaprint® 3D cast

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# luxaprint® 3D cast

UV-curing resin based on (meth) acrylate, free of bisphenol A, for the generative fabrication of casts for soft silicone earmoulds, for systems with UV / LED 385 nm

## Machine and material parameters

| Selection of device            | Software                      | Resin                                  | Parameters        |
|--------------------------------|-------------------------------|--|-------------------|
| Rapidshape HA90;<br>HA60; HA30 | netfabb Professional<br>V.5.0 | luxaprint® 3D cast<br>grün-transparent | luxaprint 3D cast |

Printers may only be operated using material parameters authorised by DETAX!

### Processing

- Observe process description on page 2.
- After storage, the material in the bottle should be shaken intensively and homogenized with a bottle roller before use.
- Processing temperature  $23^{\circ}\text{C} \pm 2^{\circ}\text{C}$ .

### Safety information

- Safety Information
- Only for the specified intended use by trained specialists.
- Contains (meth)acrylates and phosphine oxides. Some ingredients of luxaprint® 3D cast may cause allergic reactions in predisposed persons. In such cases refrain from using the product.
- Avoid direct contact with the liquid material and the components before post-curing.
- Irritating to eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- In case of contact with skin, wash immediately with plenty of soap and water.
- The processing instructions and precautionary measures must be observed.
- Refer to the relevant safety data sheet for hazard and safety information.

### Important notes

- DETAX does not accept liability for any damage caused by misuse of the impression material.
- Always keep container tightly sealed, immediately close the container carefully after each use.
- Product may cause allergic reactions.
- Not suitable for fabricating ear moulds.
- Read and understand the safety data sheet!

### Storage

- **luxaprint® 3D cast** UV is to be stored dry (at  $15^{\circ}\text{C}$  -  $25^{\circ}\text{C}$ ) and protected from light. Minimal influence of light can already induce polymerisation.

### Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations. Product may cause allergic reactions.

**Area of application:**  
generative cast moulds

**Processing:**  
at  $23^{\circ}\text{C} \pm 2^{\circ}\text{C}$

### Storage

  $25^{\circ}\text{C}$   
 $59^{\circ}\text{F}$

**Ordering information:**

**Cast Separator**  
500 ml **03636**

**luxaprint® 3D cast** 385  
1.000 g  
green-transparent **03520**

**detax softwear® 2.0  
Automix2 System**

**Cartridge pack**  
**40 Shore A**  
8 x 50 ml cartridges  
clear-transparent **04230**  
rose-transparent **04231**

2 x 50 ml cartridges  
red opaque **04232**  
blue opaque **04233**  
green opaque **04235**  
yellow-fluorescent **04234**  
pink-fluorescent **04236**  
night glow **04229**

**Cartridge pack**  
**60 Shore A**  
8 x 50 ml cartridges  
clear-transparent **03906**  
rose-transparent **03911**

2 x 50 ml cartridges  
red opaque **04242**  
blue opaque **04238**  
green opaque **04240**  
yellow-fluorescent **04239**  
pink-fluorescent **04241**  
night glow **04237**

**Mixing tips yellow** **02706**  
Pack à 48 pcs.

**Mixing gun** **02699**  
1 pc.

**softwear coat** **03509**  
500 ml

**More information at**  
[www.detax.de](http://www.detax.de)

**DETAX**

GmbH & Co. KG  
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany  
Telefon: + 49 (0) 7243/510-0 · Fax: + 49 (0) 7243/510-100  
[www.detax.de](http://www.detax.de) · [post@detax.de](mailto:post@detax.de)

# luxaprint® 3D cast

## Manufacturing process for casts

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer.

### Construction process

Generation of a Print Job complying with machine and material parameters, see table - netfabb Professional: *Selection of device + parameters*

### Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

### Cleaning interior of cast

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity  $\geq$  98%) for 3 min. in an ultrasonic bath.

Ensure that the cast is completely filled with isopropyl alcohol.  
Then air-clean with compressed air.

### Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity  $\geq$  98%) for 3 min. in an ultrasonic bath.

### Drying

Heat the construction components for 30 min. in an oven to approx. 40°C to remove the solvent residues from the cleaning process.

### Post-exposure

Post-exposure is performed with a xenon photoflash (e.g. Otoflash G171) with 2 x 1000 flashes under inert gas conditions (nitrogen), rotate components in between.

### Coating with Cast Separator

see instructions for Cast Separator

### Filling of the cast

Then the cast can be filled with silicone, i.e. detax softwear® 2.0.