# Implant Mask scan

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# Implant Mask scan

EN  
Instructions for use

Flexible gingival mask for dental laboratory implant models, A-silicone based 1:1, light bodied


Fig. 1


Fig. 2

Fig. 3

Fig. 4

Fig. 5


Fig. 6


Fig. 7


Fig. 8

Fig. 9


Fig. 10


Fig. 11

## 1. Preparation

a) Fabrication by direct method, application to impression

Mark the edge of the gingival areas for the impression with soft wax or similar. Then insulate this area by spraying lightly with silicone Separating Agent. Drying time 30 sec. Silicone Separating Agent is suitable for A- and C-silicones and polyether impression material.

For the application of **Implant Mask scan** with the Automix2 cartridge, always use the blue, short mixer. Inject directly around the impression or transfer posts and leave the tip of the mixer in the material during dispensing (Fig. 1). Dispense with steady pressure. Ensure that the **Implant Mask scan** has the desired dimensions.

Continue model fabrication only after the **Implant Mask scan** has set completely (Fig. 2). Any separating medium residues on the laboratory posts can be removed with isopropyl alcohol. Insulation of the mask against plaster or resins is not necessary.

b) Fabrication by indirect method, on the model

To fabricate the **Implant Mask scan** first make a silicone mould of addition- or condensation curing putty (e.g. **blue eco** or **compact lab putty**), which covers all of the areas of the unsawn model to be reproduced (Fig. 3). An available impression can be used as mould after removing the tray. Remove any areas in the index or impression that are in the way outside the gingival mask region so that it can be repositioned easily on the model.

Reduce the regions of the model that are to be reproduced by **Implant Mask scan** by milling (Fig. 4). Take care to remove enough material to ensure the later thickness of the mask. Sawing and working the dies can take place before or after mask fabrication depending on the situation. If sawing is performed beforehand, block out the saw cuts.

Drill one or more injection openings (palatal/lingual) in the silicone. Drill air vent channels through the mould with a round bur starting from the highest points of the mask (Fig. 5), clean the mould and insulate it on the inside by spraying lightly with silicone Separating Agent. Drying time 30 sec. Put the mould back on the model, check that it is seated exactly and fix it.

Inject material quickly through the injection opening of the mould with the Automix2 cartridge (Fig. 6). Inspection: air vent channels must be filled with silicone. Demould the mask only after complete setting. Setting can take place in the pressure pot under pressure (max. 2.5 bar), not in water and without the addition of heat.

## 2. Mixing and applying

The cartridge is extruded with the Automix2 system gun. (Fig. 7) Insert cartridge into Automix2-System mixing gun. Remove and dispose of closure cap. For exact flow control extrude slightly material until uniform flow from both orifices is achieved (Fig. 8). Mount the mixing cannula by using the guides on the cannula and the cartridge and lock in place by rotating in the opposite direction (Fig. 9). Feed material applying even pressure. Prior to use, dispense a small amount of material and visually check that the base and catalyst paste are mixed omogeneously (Fig. 10). Only then dosing can be performed individually. After each use leave the used mixing cannula in place to seal the cartridge.

## 3. Processing

When the mask has been fabricated on the model, remove the mould carefully from the model, trim any casting residues from the air vents with a scalpel and remove the mask. Remove any flash with a scalpel or sharp scissors.

Further processing of the removed mask is possible with suitable rotating tools (Fig. 11) (Molloplast® cutters, grinding sleeves with 15.000 - 25.000 r.p.m., tungsten carbide drills and arbors).

## Important working hints

- **Implant Mask scan** does not undergo any changes in dimensions.
- Latex gloves and latex-contaminated surfaces can affect the setting of **Implant Mask scan**.
- The exact fit of the impression posts and model implants must be checked prior to fabrication.
- For demo and presentation models **Lustrol** gloss varnish can be applied for a high-gloss surface.
- Please follow the instructions of the safety data sheet!

Only for dental use by qualified personnel.

aprox. 2 min.  
Working time


aprox. 6 min.  
Setting time

## Indications for use:

- Uses gingival mask for:
- Implant constructions
- Combination technique, taper and telescopic technique
- Master and saw models
- For unpowdered digital 3D recording for all CAD/CAM technologies

## Technical data:

- **Mixing volume:** 50 ml (cartridge)
- **Mixing ratio:** 1:1
- **Colour code:** base: brown catalyst: pink
- **Mixing time:** omitted (Automix2 system)
- **Working time:** approx. 2 Min.\*
- **Setting time:** approx. 6 Min.\*
- **Recovery from deformation:** approx. 99 %
- **Linear dimensional change:** approx. 0.3 %
- **Application:** at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- **Storage:**

\* from initiation of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard these times.

## Ordering information:

**Implant Mask scan**  
Cartridge packing 03552  
cartridges of 2 x 50 ml  
10 mixing cannulas,  
blue, short  
15 ml Separating Agent

Mixing cannulas 02709  
blue, short, 48 pcs.

Standard packing 03374  
cartridges of 2 x 10 ml  
6 mixing cannulas,  
light green  
15 ml Separating Agent

Mixing cannulas 02605  
light green 1:1, 25 pcs.

**Implant Mask**  
Standard packing 02529  
cartridges of 2 x 10 ml  
6 mixing cannulas,  
light green  
15 ml Separating Agent

**Esthetic Mask**  
Standard packing 02340  
base, tube of 50 ml  
catalyst, tube of 50 ml  
15 ml Separating Agent  
2 disposable syringes

**Esthetic Mask automix**  
Standard packing 03113  
cartridges of 2 x 50 ml  
10 mixing cannulas, yellow  
15 ml Separating Agent

Mixing gun 02699  
Automix2, 1 pc.

Mixing cannulas 02706  
yellow, 48 pcs.

Intra-Oral-Tips 02712  
yellow, 96 pcs.

Separating Agent 02690  
separating liquid  
for silicones, vaporizer  
of 15 ml

Lustrol 03008  
gloss varnish bottle of 6 ml  
catalyst, bottle of 6 ml

Grinding sleeves 02615  
pointed, 10 pcs.

Grinding sleeve holder 02616  
pointed, 1 pcs.

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