Implant Mask / scan

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Implant Mask / scan

EN
Instructions for use

Flexible gingival mask for dental laboratory implant models, A-silicone based 1:1, light bodied

Indications for use:

- Uses Gingival mask for:
- Implant constructions
- Combination technique
- Master and saw models

Implant Mask scan:

- for unpowdered digital 3D recording for all CAD/CAM technologies

Technical data:

- **Mixing volume:** 10 ml (cartridges)
- **Mixing ratio:** 1:1
- **Colour code:** base: dark red catalyst: white
- **Mixing time:** omitted (mini-mix system)
- **Working time:** approx. 2 Min.*
- **Setting time:** approx. 5 - 6 Min.*
- **Recovery from deformation:** approx. 99 %
- **Linear dimensional change:** approx. 0,3 %
- **Application:** at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- **Storage:**

* from initiation of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard these times.

Ordering information:

Implant Mask
Standard packing 02529
Cartridges of 2 x 10 ml
6 Mixing cannulas,
light green
15 ml Separating Agent

Implant Mask scan
Standard packing 03374
Cartridges of 2 x 10 ml
6 Mixing cannulas,
light green
15 ml Separating Agent

Mixing cannulas 02605
light green 1:1, 25 pcs.

Separating Agent 02690
separating liquid
for silicones, vaporizer
of 15 ml

Lustrol 03008
gloss varnish bottle of
6 ml
catalyst, bottle of 6 ml

Grinding sleeves 02615
pointed, 10 units

Grinding sleeve holder 02616
pointed, 1 unit

Esthetic Mask
Standard packing 02340
base, tube of 50 ml
catalyst, tube of 50 ml
15 ml Separating Agent
2 disposable syringes

Esthetic Mask automix
Standard packing 03113
cartridges of 2 x 50 ml
10 mixing cannulas,
yellow
15 ml Separating Agent

Mixing gun 02699
Automix2, 1 pc.

Mixing cannulas 02706
yellow, 48 pcs.

Intra-Oral-Tips 02712
yellow, 96 pcs.

1. Preparation

a) Fabrication by direct method, application to impression

Mark the edge of the gingival areas for the impression with soft wax or similar. Then insulate this area by spraying lightly with silicone **Separating Agent**. Drying time 30 sec. Silicone **Separating Agent** is suitable for A- and C-silicones and polyether impression material.


Fig. 1

For the application of **Implant Mask / scan** with the mini-mix cartridge, always use a mixer with a fine tip. Inject directly around the impression or transfer posts and leave the tip of the mixer in the material during dispensing (Fig. 1). Dispense with steady pressure. Ensure that the **Implant Mask / scan** has the desired dimensions.


Fig. 2

Continue model fabrication only after the **Implant Mask / scan** has set completely (Fig. 2). Any separating medium residues on the laboratory posts can be removed with isopropyl alcohol. Insulation of the mask against plaster or resins is not necessary.

b) Fabrication by indirect method, on the model

To fabricate the **Implant Mask / scan** first make a silicone mould of addition- or condensation curing putty (e.g. **blue eco** or **compact lab putty**), which covers all of the areas of the unsawn model to be reproduced (Fig. 3). An available impression can be used as mould after removing the tray. Remove any areas in the index or impression that are in the way outside the gingival mask region so that it can be repositioned easily on the model.

Fig. 3

Reduce the regions of the model that are to be reproduced by **Implant Mask / scan** by milling (Fig. 4). Take care to remove enough material to ensure the later thickness of the mask. Sawing and working the dies can take place before or after mask fabrication depending on the situation. If sawing is performed beforehand, block out the saw cuts.

Fig. 4

Drill one or more injection openings (palatal/lingual) in the silicone. Drill air vent channels through the mould with a round bur starting from the highest points of the mask (Fig. 5), clean the mould and insulate it on the inside by spraying lightly with silicone **Separating Agent**. Drying time 30 sec. Put the mould back on the model, check that it is seated exactly and fix it.

Fig. 5

Inject material quickly through the injection opening of the mould with the mini-mix cartridge (Fig. 6). Inspection: air vent channels must be filled with silicone. Demould the mask only after complete setting. Setting can take place in the pressure pot under pressure (max. 2.5 bar), not in water and without the addition of heat.

Fig. 6

2. Mixing and applying

The material is extruded with the mini-mix system. Place the plunger in the body of the cartridge. Remove the cartridge cap by twisting off. Before attaching the mixing cannula extrude a small amount of material until silicone is dispensed evenly from the two openings (Fig. 7).


Fig. 7

To attach the mixing cannula, note the guides on the cannula and cartridge. Lock by turning them in opposite directions (Fig. 8). Feed material applying even pressure. Prior to use, dispense a small amount of material and visually check that the base and catalyst paste are mixed homogeneously (Fig. 9). Only then dosing can be performed individually.

Fig. 8

3. Processing

When the mask has been fabricated on the model, remove the mould carefully from the model, trim any casting residues from the air vents with a scalpel and remove the mask. Remove any flash with a scalpel or sharp scissors.


Fig. 9

Further processing of the removed mask is possible with suitable rotating tools (Fig. 10) (**Molloplast®** cutters, grinding sleeves with 15,000 - 25,000 r.p.m., tungsten carbide drills and arbors).

Fig. 10

Important working hints

- **Implant Mask / scan** does not undergo any changes in dimensions.
- Latex gloves and latex-contaminated surfaces can affect the setting of **Implant Mask / scan**.
- The exact fit of the impression posts and model implants must be checked prior to fabrication.
- Where movement is impeded or blocked, discard the cartridge. Do not apply force.
- For demo and presentation models **Lustrol** gloss varnish can be applied for a high-gloss surface.
- Please follow the instructions of the safety data sheets!

Only for dental use by qualified personnel.

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