

# AUDIO 3D RESINS

## HIGH PERFORMANCE POLYMERS



» PRODUCT CATALOGUE «  
2023

**DETAX**  
HIGHEND MEDICAL MATERIALS



# CONTENT

Product overview 3D Premium Resins.....	3
Product overview Accessories.....	4
Applications.....	5
Matrix.....	6
LUXAPRINT® MOULD.....	8
LUXAPRINT® SHELL.....	10
LUXAPRINT® FLEX / COAT.....	12
LUXAPRINT® CAST 2.0.....	14
LUXAPRINT® COCOON.....	16
MEDICALPRINT® MOULD.....	18
MEDICALPRINT® SHELL.....	20
MEDICALPRINT® NOBREAK.....	22
LUXAPRINT® shellac.....	24
LUXAPRINT® shellac color.....	26
LUXAPRINT® softseal Primer.....	28
Cast Separator 2.0.....	30
Cast Separator PU.....	30
Validated 3D printers.....	32
Validated curing devices.....	34
Good to know.....	35
3D Workflow.....	36
Certification.....	38
Product clips.....	39

# PRODUCT OVERVIEW

## 3D PREMIUM RESINS

### MEDICALPRINT® NOBREAK

In-Ear-Monitoring  
Earmoulds  
Hearing Protection  
ITE Shells



### MEDICALPRINT® MOULD

In-Ear-Monitoring  
Earmoulds  
Hearing Protection  
ITE Shells



### MEDICALPRINT® SHELL

ITE Shells  
In-Ear-Monitoring  
Foil-Earmoulds  
Hearing Protection (active)



### LUXAPRINT® MOULD

In-Ear-Monitoring  
Earmoulds  
Hearing Protection  
ITE Shells



### LUXAPRINT® SHELL

ITE Shells  
In-Ear-Monitoring  
Foil-Earmoulds  
Hearing Protection (active)



### LUXAPRINT® FLEX

Swim Plugs  
Hearing Protection  
Earmoulds  
In-Ear-Monitoring



### LUXAPRINT® CAST 2.0

CASTINGS/Casting Mould



### LUXAPRINT® COCOON

CASTINGS/Special Casting Mould



# PRODUCT OVERVIEW

## ACCESSORIES

### luxaprint® shellac

High gloss sealing:  
ITE Shells  
In-Ear-Monitoring  
Earmould  
Hearing Protection



### luxaprint® shellac color

Coloured sealing:  
ITE Shells  
In-Ear-Monitoring  
Earmould  
Hearing Protection



### luxaprint® flex coat

Soft surface sealing:  
In-Ear-Monitoring  
Hearing Protection  
Earmoulds  
Swim Plugs



### luxaprint® softseal Primer

Bonding of acrylates and silicones  
(SoftTip)



### Cast Separator 2.0

Separating agent for  
CASTINGS/Casting Mould



### Cast Separator PU

Separating agent for  
CASTINGS/Casting Mould



# APPLICATIONS

### ITE SHELLS

luxaprint® shell  
luxaprint® mould  
medicalprint® shell  
medicalprint® mould  
medicalprint® nobreak




### EARMOULD

luxaprint® mould  
luxaprint® flex  
medicalprint® mould  
medicalprint® nobreak

### HEARING PROTECTION

luxaprint® mould  
luxaprint® flex  
medicalprint® mould  
medicalprint® nobreak




### FOIL-EARMOULD

luxaprint® shell  
luxaprint® mould  
medicalprint® shell  
medicalprint® mould  
medicalprint® nobreak

### HEARING PROTECTION (ACTIVE)

luxaprint® shell  
luxaprint® mould  
medicalprint® shell  
medicalprint® mould  
medicalprint® nobreak




### IN-EAR-MONITORING

luxaprint® shell  
luxaprint® mould  
medicalprint® shell  
medicalprint® mould  
medicalprint® nobreak

### SWIM PLUG

luxaprint® cast 2.0  
luxaprint® flex








### CASTINGS/ CASTING MOULD

luxaprint® cast 2.0  
luxaprint® cocoon

# MATRIX

LUXAPRINT®

Material type	Application	Colour	Colour intensity low high	Product features	Medical devices Class MDR
<b>mould</b> 	• In-Ear-Monitoring • Earmoulds • Hearing Protection • ITE Shells	blue, clear, light beige, intensive blue, rose, rose-orange, red	transparent	• low viscosity • highest precision • maximum initial hardness • very high surface hardness • biocompatible	Ila
<b>shell</b> 	• ITE Shells • In-Ear-Monitoring • Foil-Earmoulds • Hearing Protection (active)	beige, intensive blue, intensive red, black, white	transparent opaque	• low sedimentation tendency • highest precision • maximum initial hardness • very high surface hardness • optimal depth curing • biocompatible	Ila
<b>flex</b> 	• In-Ear-Monitoring, • Hearing Protection • Earmoulds • Swim Plugs	clear	transparent	• high impact resistance with memory effect • fast elastic recovery • 90 Shore A at 23 °C • 70 Shore A at 37 °C • biocompatible	Ila
<b>cast 2.0</b> 	• Casting Mould for silicone earmoulds / polyurethane	green-transparent	transparent	• easy to remove • very low viscosity	TEC resin
<b>COCOON</b> 	• Special Casting Mould for silicone earmoulds	clear-transparent	transparent	• easy to remove < 2 sec • flexible and stable in form • filling level control • usage without separator	TEC resin




## LUXAPRINT®

- › Bisphenol A-free, MMA free
- › highest initial hardness
- › minimal initial discolouration
- › highest surface hardness
- › drying process




## MEDICALPRINT®

- › MMA free
- › high initial hardness
- › no initial discolouration
- › highest impact resistance
- › no drying process

MEDICALPRINT®

Material type	Application	Colour	Colour intensity low high	Product features	Medical devices Class MDR
<b>mould</b> 	• Earmoulds • Hearing Protection • In-Ear-Monitoring • ITE Shells	brilliant-clear, rose, rose-orange	transparent	• low viscosity • high impact resistance • high initial transparency • maximal precision • biocompatible	Ila
<b>shell</b> 	• ITE Shells • Foil-Earmoulds • Hearing Protection (active) • In-Ear-Monitoring	beige, blue-opaque, red-opaque, black, skin, white	opaque	• low sedimentation tendency • high impact resistance • break resistance • maximal precision • optimal depth curing • biocompatible	Ila
<b>nobreak</b> 	• Earmoulds • Hearing Protection • In-Ear-Monitoring • ITE Shells	beige, clear, rosé	transparent	• highest breaking strength • max. impact resistance • high initial transparency • easy processing	in process

LUXAPRINT® LACQUERS

Material type	Application	Colour	Colour intensity low high	Product features	Medical devices Class MDR
<b>shellac</b> 	High gloss sealing: • ITE Shells • In-Ear-Monitoring • Earmoulds • Hearing Protection	transparent	transparent	• easy to clean surface • „blue ray“ yellowing protection • high surface hardness • strong adhesion • very low viscosity • biocompatible	Ila
<b>shellac color</b> 	Coloured sealing: • ITE Shells • In-Ear-Monitoring • Earmoulds • Hearing Protection	blue, brown, yellow, green, orange, red, black, violet	transparent	• colours freely mixable • long-term colour stable • extended colour stability • scratch proof • very low viscosity • strong adhesion • biocompatible	Ila
<b>flex coat</b> 	Soft surface sealing for hard & soft Earmoulds: • In-Ear-Monitoring, • Hearing Protection • Earmoulds • Swim Plugs	transparent	transparent	• permanently soft • strong bonding • „blue ray“ yellowing protection • low viscosity • perfect for dipping • biocompatible	Ila

 MDR certified	 Medical Product Class I	 Medical Product Class IIa	 Technical Product	 Eco-Bag 3/5 kg
 36 Months Shelf life	 Premium quality	 THF-MA free	 Bisphenol A free	 MMA free

**THF-MA-free** does not contain any tetrahydrofurfuryl methacrylate [Reproductive Toxicity, Cat. 1B] **BPA-free** does not contain any raw material based on bisphenol A [Reproductive Toxicity, Cat. 1B] **MMA-free** does not contain any methyl methacrylate **MDR** Medical Device Regulation EU



# LUXAPRINT® MOULD

IN-EAR-MONITORING  
EARMOULDS  
HEARING PROTECTION  
ITE SHELLS

Light-curing formulation for 3D printing of hard earmoulds & hearing protection.

**Colours:** blue, clear, light beige, intensive blue, rose, rose-orange, red  
**Wavelength:** 385 nm  
**Medical Device Class IIa**

- High initial hardness
- Max construction precision
- Optimum depth curing
- Low material consumption




Property	Standard	Unit measurement	Result
Flexural strength	DIN EN ISO 178 <sup>1)</sup>	MPa	> 75
Flexural modulus	DIN EN ISO 178 <sup>1)</sup>	MPa	> 1750
Tensile strength	DIN EN ISO 527-1 <sup>2)</sup>	MPa	> 47
Elongation	DIN EN ISO 527-1 <sup>2)</sup>	-	> 9 %
Hardness	-	Shore D	> 84
Biocompatibility	DIN EN ISO 10993-1 <sup>3)</sup>	-	complies

<sup>1)</sup> Plastics: Determination of flexural properties (in accordance with the norm at room temperature)  
<sup>2)</sup> Plastics: Determination of tensile properties (in accordance with the norm at room temperature)  
<sup>3)</sup> Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process


light beige*		1.000 g	03973	
blue	500 g	03716	1.000 g	03610
clear	500 g	03717	1.000 g	03608
intensive blue		1.000 g	03728	5 kg 03325
rose	500 g	03718	1.000 g	03609
rose-orange*	500 g	03718	1.000 g	03946
red	500 g	03715	1.000 g	03611
				5 kg 03324
				5 kg 03323

\* upon request



Wide range of classic colours of supreme transparency, especially for hearing protection and hearing aid earmoulds.

Maximum production process reliability thanks to high green part/initial hardness and strength.



Very high mechanical flexural and tensile strength without brittleness.



# LUXAPRINT® SHELL

ITE SHELLS  
IN-EAR-MONITORING  
FOIL-EARMOULDS  
HEARING PROTECTION (ACTIVE)

Light-curing formulation for 3D printing of hard ITE Shells.

**Colours:** beige, intensive blue, intensive red, black, white

**Wavelength:** 385 nm

**Medical Device Class IIa**

- Low viscosity
- High initial hardness
- Max construction precision
- Very high surface hardness



Property	Standard	Unit measurement	Result
Flexural strength	DIN EN ISO 178 <sup>1)</sup>	MPa	> 70
Flexural modulus	DIN EN ISO 178 <sup>1)</sup>	MPa	> 1800
Tensile strength	DIN EN ISO 527-1 <sup>2)</sup>	MPa	> 47
Elongation	DIN EN ISO 527-1 <sup>2)</sup>	-	> 9 %
Hardness	-	Shore D	> 82
Biocompatibility	DIN EN ISO 10993-1 <sup>3)</sup>	-	complies


<sup>1)</sup> Plastics: Determination of flexural properties (in accordance with the norm at room temperature)  
<sup>2)</sup> Plastics: Determination of tensile properties (in accordance with the norm at room temperature)  
<sup>3)</sup> Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process

beige	500 g	03538	1.000 g	03513
intensive blue	500 g	03720	1.000 g	03591
intensive red			1.000 g	03601
black	500 g	03844	1.000 g	03843
white	500 g	03842	1.000 g	03841




Classic or modern: The color palette offers the right colour for every application. Especially suitable for high-end monitoring.

Very high mechanical flexural and fracture strength without brittleness, even in thin-walled ITE shells.



Targeted colour setting for very high opacity with simultaneous optimum curing.





# LUXAPRINT® FLEX

SWIM PLUGS  
HEARING PROTECTION  
EARMOULDs  
IN-EAR-MONITORING

**luxaprint® flex:** Light-curing formulation for 3D printing of soft, massive earmoulds.  
**luxaprint® flex coat:** UV curing one component lacquer for sealing of 3D printed hard and soft earmoulds.

**Colour:** clear

**Wavelength:** 385 nm

**Medical Device Class IIa**

- Permanently soft
- Memory effect
- High impact resistance
- Fast elastic recovery



Property	Standard	Unit measurement	Result	
Tensile strength	DIN EN ISO 527-1 <sup>1)</sup>	MPa	> 8	> 8
Elongation	DIN EN ISO 527-1 <sup>1)</sup>	-	> 60 %	> 60 %
Tear strength	DIN ISO 34-1 <sup>2)</sup>	N/mm	> 35 without luxaprint® flex coat	> 45 with luxaprint® flex coat
Hardness	-	Shore A	> 90 at room temperature	> 70 at body temperature
Biocompatibility	DIN EN ISO 10993-1 <sup>3)</sup>	-	complies	


<sup>1)</sup>Plastics: Determination of flexural properties (in accordance with the norm at room temperature)

<sup>2)</sup>Rubber, vulcanized or thermoplastic: Determination of tear strength (in accordance with the norm at room temperature)

<sup>3)</sup>Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process



04246	luxaprint® flex	500 g
04245	luxaprint® flex	1.000 g
04247	luxaprint® flex coat*	100 ml

\* not MMA-free



Thanks to the memory effect, the earmould returns to its original shape after deformation, retaining the dimensions.

Earmoulds printed with luxaprint® flex guarantee a natural feel: rigid at room temperature, flexible and supple at body temperature.



The permanently soft sealing with luxaprint® flex coat improves tear resistance and impact strength and protects the flexible earmold from dirt and cerumen adhesion.



# LUXAPRINT® CAST 2.0

## CASTINGS/CASTING MOULD

Light-curing formulation for 3D printing of hard cast forms, manufacturing of silicone earmoulds or polyurethane earmoulds.

**Colour:** green-transparent

**Wavelength:** 385 nm

**Technical Product**

- Easy to break
- Predefined brittleness
- Low viscosity
- „easy peel off“ with Cast Separator 2.0



Property	Standard	Unit measurement	Result
Flexural strength	DIN EN ISO 178 <sup>1)</sup>	MPa	> 75
Flexural modulus	DIN EN ISO 178 <sup>1)</sup>	MPa	> 2300
Hardness	-	Shore D	> 85

<sup>1)</sup> Plastics: Determination of flexural properties (in accordance with the norm at room temperature)

03918	green-transparent	1.000 g
03327	green-transparent	5 kg



Highest process reliability and production precision with minimum wall thickness.



The high stability of the material enables production of very thin-walled cast forms. The transparency enables visual control of the filling process.



Effortless break-up of the cast form due to defined brittleness and easy peel-off when using the cast separator, especially for more complicated structures.





# LUXAPRINT® COCOON

## CASTINGS/CASTING MOULD

Light-curing formulation for 3D printing of soft elastic, transparent cast forms, manufacturing of silicone earmoulds.

**Colour:** clear-transparent

**Wavelength:** 385 nm

**Technical Product**

- Flexible and stable in form
- Easy to remove
- No sticking to the form
- Common VPS silicones



Property	Standard	Unit measurement	Result
Tensile strength	DIN EN ISO 527-1 <sup>1)</sup>	MPa	> 7
Tear strength	DIN ISO 34-1 <sup>2)</sup>	N/mm	> 30
Hardness	-	Shore A	> 90

<sup>1)</sup> Plastics: Determination of flexural properties (in accordance with the norm at room temperature)


<sup>2)</sup> Rubber, vulcanized or thermoplastic: Determination of tear strength (in accordance with the norm at room temperature)

03031	clear-transparent	1.000 g
03326	clear-transparent	5 kg



Fast peel-off – quick removal of the blank: No adhesion of silicone to the cast form; no additional separator necessary.

The transparent cast form made of luxaprint® cocoon allows visual control of the filling process; air pockets are avoided.



Due to the predetermined tear line and the flexible material properties, the earmould can be demolded without tools. This prevents formation of sharp-edged fragments that could damage the silicone earmould. The nominal tear line can be easily created directly in the design software (3shape, Cyfex).



# MEDICALPRINT® MOULD

IN-EAR-MONITORING  
EARMOULDS  
HEARING PROTECTION  
ITE SHELLS

Light-curing formulation for 3D printing of earmoulds, In-Ear-Monitoring & ITE shells.

**Colours:** brilliant-clear, rose, rose-orange

**Wavelength:** 385 nm

**Medical Device Class IIa**

- High initial transparency
- High impact resistance
- Max. construction precision
- No drying process required



Property	Standard	Unit measurement	Result
Flexural strength	DIN EN ISO 178 <sup>1)</sup>	MPa	> 75
Flexural modulus	DIN EN ISO 178 <sup>1)</sup>	MPa	> 1900
Tensile strength	DIN EN ISO 527-1 <sup>2)</sup>	MPa	> 50
Elongation	DIN EN ISO 527-1 <sup>2)</sup>	-	> 4 %
Hardness	-	Shore D	> 80
Biocompatibility	DIN EN ISO 10993-1 <sup>3)</sup>	-	complies

<sup>1)</sup> Plastics: Determination of flexural properties (in accordance with the norm at room temperature)

<sup>2)</sup> Plastics: Determination of tensile properties (in accordance with the norm at room temperature)

<sup>3)</sup> Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process



brilliant-clear 2.0	1.000 g	04224	5 kg	03316
rose 2.0	1.000 g	02449	5 kg	03314
rose-orange 2.0*	1.000 g	02904	5 kg	03304

\* upon request



Brilliant transparency directly after printing. Long-term colour stability, without yellowing effects due to post exposure.

The high impact strength and good mechanical flexural and fracture strengths ensure long service life of the earmoulds.



Particularly suitable for delicate foil earmolds and RIC earmoulds, with long supports.





# MEDICALPRINT® SHELL

## ITE SHELLS

## IN-EAR-MONITORING

## FOIL-EARMOULDS

## HEARING PROTECTION (ACTIVE)

Light-curing formulation for 3D printing of hearing aid shells, In-Ear-Monitoring, foil earmoulds..

**Colours:** beige, blue-opaque, red-opaque, black, skin, white

**Wavelength:** 385 nm

**Medical Device Class IIa**

- Very high construction precision
- High impact resistance
- Reduced sedimentation tendency
- No drying process required



Property	Standard	Unit measurement	Result
Flexural strength	DIN EN ISO 178 <sup>1)</sup>	MPa	> 75
Flexural modulus	DIN EN ISO 178 <sup>1)</sup>	MPa	> 1900
Tensile strength	DIN EN ISO 527-1 <sup>2)</sup>	MPa	> 50
Elongation	DIN EN ISO 527-1 <sup>2)</sup>	-	> 4 %
Hardness	-	Shore D	> 80
Biocompatibility	DIN EN ISO 10993-1 <sup>3)</sup>	-	complies


<sup>1)</sup>Plastics: Determination of flexural properties (in accordance with the norm at room temperature)

<sup>2)</sup>Plastics: Determination of tensile properties (in accordance with the norm at room temperature)

<sup>3)</sup>Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process


03016	beige 2.0	1.000 g
04164	blue-opaque*	1.000 g
04165	red-opaque*	1.000 g
03043	black*	1.000 g
02192	skin 2.0	1.000 g
04073	white*	1.000 g

\* not THF-MA free



Wide range of classic and modern colours, long-term stability, particularly for thin-walled earmolds.

Very high production precision and surface hardness, outstanding mechanical flexural and breaking strength.



The intense colours are particularly suitable for hearing protection earmoulds and can be optimally used in the area of in-ear monitoring for visual side marking.





# MEDICALPRINT® NOBREAK

IN-EAR-MONITORING  
EARMOULDS  
HEARING PROTECTION  
ITE SHELLS

Light-curing formulation for 3D printing of earmoulds, In-Ear-Monitoring & ITE shells.

Colours: beige, clear, rose  
Wavelength: 385 nm

- Highest breaking strength
- Max. impact resistance
- High initial transparency
- Easy processing




Property	Standard	Unit measurement	Result
Tensile strength	DIN EN ISO 527-1 <sup>1)</sup>	MPa	> 45*
Elongation	DIN EN ISO 527-1 <sup>1)</sup>	-	> 7 %*
Hardness	-	Shore D	> 70*

<sup>1)</sup> Plastics: Determination of flexural properties (in accordance with the norm at room temperature)

04429	medicalprint® nobreak	1.000 g
04430	medicalprint® nobreak	5 kg

\* Q1/2023



Classic color selection for the fabrication of ITE shells and earmolds.

No cracks, breaks or chipping due to extremely high impact resistance.

Maximum space efficiency thanks to extremely thin-walled shells and foil earmolds.

# LUXAPRINT® shellac

HIGH GLOSS SEALING:

ITE SHELLS

IN-EAR-MONITORING

EARMOULD

HEARING PROTECTION

Light curing lacquer for permanent surface sealing of 3D printed earmoulds made of (meth)acrylat.

**Colour:** transparent

**Medical Device Class IIa**


- Easy to clean surface
- High surface hardness
- Strong bonding
- Very low viscosity
- „blue ray“ - yellowing protection



04006	transparent	50 ml
03594	transparent	100 ml
03595	transparent	300 ml



Brilliant finish without mechanical repolishing.



Easy to process (dip/brush). Creates a homogeneous, scratch-resistant and easy-to-clean surface (easy clean) and reduces cerumen adhesion.



Can be used in combination with luxaprint® shellac color for mixing individual colours or for marking lettering.





# LUXAPRINT® shellac color

## COLOURED SEALING:

ITE SHELLS  
IN-EAR-MONITORING  
EARMOULD  
HEARING PROTECTION


Light curing lacquer for permanent coloured sealing of 3D printed earmoulds made of (meth)acrylat.

**Colours:** blue, brown, yellow, green, orange, red, black, violet  
**Medical Device Class IIa**

- Colours freely mixable
- Long-term colour stable
- Extended colour stability
- Scratch-proof
- Very low viscosity



blue	50 ml	03701	100 ml	03683
brown	50 ml	03995	100 ml	03994
yellow	50 ml	03702	100 ml	03684
green	50 ml	03705	100 ml	03678
orange	50 ml	03703	100 ml	03685
red	50 ml	03700	100 ml	03682
black	50 ml	03857	100 ml	03856
violet	50 ml	03704	100 ml	03686



**FREELY  
MIXABLE**

Transparent earmoulds can be simultaneously sealed with a high-gloss finish and permanently colour-coated. The surface is protected and has long-term colour stability.

**Skin tones unlimited:** All colours are individually mixable. As a result, any desired skin tone can be achieved. There are no limits to the colour wishes.



Easy to process (dip/brush). Creates a homogeneous, scratch-resistant, easy-to-clean surface and reduces cerumen adhesion.





# LUXAPRINT® softseal Primer

## BONDING OF (METH-)ACRYLATES AND SILICONES (SOFTTIP)

Primer for perfect binding of (Meth-)acrylates (luxaprint®, medicalprint®) and silicones (detax softwear® 2.0, earflex®).

**Colour:** transparent  
**Medical Device Class IIa**


- Especially for combination earmoulds (hard / soft)
- Excellent and permanent bonding
- Firm fit, even during chewing movements
- Easy handling




04161

Bottle with brush

15 ml



Reliable adhesive composite of (meth)acrylate and silicone. Easy surface conditioning by means of integrated brush.



The cast form printed with luxaprint® cast 2.0 brings the silicone into the desired shape for the soft-tip application.



Combination earmoulds (hard/soft) improve fit and wearing comfort. The soft-tip application adapts harmoniously to the chewing movements, among others.



# Cast Separator 2.0

## SEPARATING AGENT FOR CASTINGS/CASTING MOULD

Separating agent for an effective isolation of 3D printed cast forms against VPS or A-silicone.

**Colour:** transparent  
**Technical Product**



03636

Cast Separator 2.0

500 ml

# Cast Separator PU

## SEPARATING AGENT FOR CASTINGS/CASTING MOULD

Separating agent for an effective isolation of 3D printed cast forms against thermoplastic polyurethane.

**Colour:** transparent  
**Technical Product**



04142


Cast Separator PU

500 ml



Provides a durable and safe insulating layer for the additively manufactured cast forms.

Effortless and easy demolding of the silicone and PU blank.























The ready-to-use immersion solution enables quick and easy application.



### Qualification

- ☒ Done  
☐ In process

												
<b>luxaprint® mould</b> In-Ear-Monitoring Earmoulds Hearing Protection ITE Shells 	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<b>luxaprint® shell</b> ITE Shells In-Ear-Monitoring Foil-Earmoulds Hearing Protection (active) 	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<b>luxaprint® flex</b> Swim Plugs Hearing Protection Earmoulds In-Ear-Monitoring 	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<b>luxaprint® cast 2.0</b> CASTINGS / Casting Mould 	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<b>luxaprint® cocoon</b> CASTINGS / Special Casting Mould 	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<b>medicalprint® mould</b> In-Ear-Monitoring Earmoulds Hearing Protection ITE Shells 	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<b>medicalprint® shell</b> ITE Shells In-Ear-Monitoring Foil-Earmoulds Hearing Protection (active) 	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<b>medicalprint® nobreak</b> In-Ear-Monitoring Earmoulds Hearing Protection ITE Shells 	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>



**CERTIFIED  
VALIDATED  
RELIABLE**

Only the combination of high-performance resins with scientifically based expert knowledge from all areas of the digital workflow leads to cumulated expertise, to truly innovative products, and thus to an unlimited choice of materials. DETAX 3D materials are validated for all standard DLP printers. Our validation portfolio is continuously being expanded with new materials and qualified printers. To this end, our experts check and document complete process sequences in accordance with the relevant standards and regulatory requirements. This ensures permanently reproducible results and constant product quality.

Digital Workflow requires profound material competence and a close cooperation with the technology partners in order to perfectly match individual elements of the process chain. For transparency and process reliability, all instructions for use comprise an overview of validated printers, certified finishing equipment (post-exposure, cleaning, etc.) and detailed flowcharts of the manufacturing process.

Our expert team will support you with useful tips.

edition: 05.10.22





# PROCESS VALIDATION

## CURING DEVICES

### Qualification

- ☒ Done
- ☐ In process
- ☐ On request

	NK Optik Ofioflash GT71 N2	Rapid Shape RS Cure	Dentalfarm Photopol	Meccatroni- core BB Cure	Scheu Imprime Cure	Phrozen Phrozen Cure	Sprint Ray Pro Cure	Formlabs Form Cure	Dreve PCU LED N2
<b>luxaprint® mould</b>									
In-Ear-Monitoring Earbuds Hearing Protection ITE Shells									
rose-orange, clear, rose, light beige red, blue, int.-blue									
<b>luxaprint® shell</b>									
ITE Shells In-Ear-Monitoring Foal-Ear-Monitoring Hearing Protection (active)									
beige, white black, int.-blue, int.-red									
<b>luxaprint® flex</b>									
Swim Plugs Hearing Protection Earbuds In-Ear-Monitoring									
clear									
<b>luxaprint® cast 2.0</b>									
CASTINGS / Casting Mould									
green-transparent									
<b>luxaprint® cocoon</b>									
CASTINGS / Special Casting Mould									
transparent									

<b>medicalprint® mould</b>									
In-Ear-Monitoring Earbuds Hearing Protection ITE Shells									
brilliant-clear 2.0, rose 2.0, rose-orange 2.0, clear, rose									
<b>medicalprint® shell</b>									
beige 2.0, skin 2.0 ITE Shells In-Ear-Monitoring Foal-Ear-Monitoring Hearing Protection (active)									
black, blue-opaque, red-opaque white									
<b>medicalprint® nobreak</b>									
In-Ear-Monitoring Earbuds Hearing Protection ITE Shells									

edition: 06.10.22

# GOOD TO KNOW ...

## BOTTLE ROLLER

By using a bottle roller, optimum mixing of the material is achieved, thus preventing possible segregation. The Eco Bags can be homogenized with an appropriate attachment.

## CLEANING

Best cleaning results of the production jobs are achieved when the pre- and post-cleaning are carried out in separate tanks in the ultrasonic unit. After cleaning with isopropanol, it is recommended to clean the bores/openings with compressed air.

## POST-CURING UNIT

The post-curing units recommended in the instructions for use ensure optimum through-hardening and surface curing, thus a biocompatible end product, and ensure high color brilliance and transparency, without discoloration.

## DETAX EXPERTS@



# 3D WORKFLOW

After completion of the design (CAD), the slicing software prepares the objects for printing. The slicing process creates the individual layers to be exposed. The software serves as an intermediary between the 3D model and the 3D printer.

After printing, the non-polymerized material on the surface must be removed so as to leave no residue before the final post-exposure. Drain the print job off in the printer, then carry out a 2-stage secondary cleaning with isopropanol in an ultrasonic device. Cleaning can also be carried out in suitable separate devices.

Finally, the surface is finished as required, e.g. mechanically polished. Perfect fit, optimal product properties and reliable reproduction are the results of a validated and certified process.

## SCANNING




Digitization of the patient's initial situation is the basis for the digital manufacturing process. It is done by scanning the ear directly, or the impression. From the data thus generated, a three-dimensional surface structure is generated – in most cases in the form of an STL file –, which is then transferred to a design software.

## SLICING



## PRINTING



For a precise print job, the setting parameters of the corresponding material in the printer are necessary. These data are used not only to control the exposure process for the material, but also to determine the corresponding movement mechanics of the printers. Coordination of these processes is the prerequisite for successful DLP/LCD printing of challenging structures.

## CLEANING



The properties of the final product depend, among other things, on the finishing process. Correct post-exposure is very important for biocompatibility. To ensure that the printed structures are fully cured, post-exposure in devices with LED lamps or xenon flashlight in an inert gas atmosphere is recommended.

## CURING



## FINISHING



# CERTIFICATION



All luxaprint® and medicalprint® Class IIa printing resins have been MDR-certified since October 2020. Thus, DETAX 3D materials are among the first in the industry with MDR certification.



All DETAX 3D Premium printing materials have a shelf life of 36 months and can be used during this period without any loss of quality in the printing process.




## #HELLO ECOBAG



In addition to the 1-kg standard bottles, many luxaprint® and medicalprint® materials will in future be offered in practical 3- or 5-kg Eco Bags. The bags are perfect for frequent users and are handy to use: The 2 handles (top and bottom) make it easy to fill the printer tray. Highly pigmented materials can be homogenized with a roller mixer (with appropriate attachment). The empty bag can be rolled up to a tiny ball, thus taking up much less waste volume and significantly reducing plastic waste.



## PRODUCTCLIPS „HOW TO“







DETAX GmbH Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany  
Tel. +49 72 43/510-0 · Fax: +49 72 43/510-100 · [www.detax.com](http://www.detax.com) · [post@detax.com](mailto:post@detax.com)