# Detaseal® hydroflow putty 5:1

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# Detaseal<sup>®</sup> hydroflow putty 5:1

**GB**  
Instructions for use

VPS precision impression silicone, kneadable consistency

Fig. 1

## 1. Preparation of the impression tray

Select a tray type suited for desired application (standard-size or custom, in metal or hard plastics). To avoid teeth contact with the tray, place stops from wax or light-curing resin into the tray at non-relevant areas of the impression. For optimal adhesion of impression material, we recommend applying a thin layer of **trayloc<sup>®</sup> A** to tray (Fig. 1). Let dry for approx. 3 minutes.


Fig. 2

## 2. Mixing and dispensing

Disengage and remove the tamper-proof seal on the cartridge by raising the locking lever (Fig. 2). Insert the cartridge into the dispenser. Start the device and dispense material until material flows uniformly. Discard the first inhomogeneous material, then attach the dynamic mixer (yellow). Fix the mixer by the applying and turning the yellow rotating ring clockwise (Fig. 3). Feed material applying even pressure. Prior to use, dispense a small amount of material and visually check that the base and catalyst paste are mixed homogeneously (Fig. 4). Then fill the tray with the desired quantity of **Detaseal<sup>®</sup> hydroflow putty**, holding the tray at an angle and pressing against it slightly. While dispensing, leave the mixer in the material (Fig. 5). Leave the filled mixer on the cartridge holder as a closure. Before the next application remove the used mixer by undoing the rotating ring and check the cartridge openings for blockages. Attach a new mixer and proceed in the normal way. When starting up, discard the initial material emerging.


Fig. 3

## 3. Procedure

Firstly, apply a sufficient amount of **Detaseal<sup>®</sup> hydroflow putty** into the impression tray, giving the teeth area a distinctly trough-like shape (Fig. 6). Then apply a layer of **Detaseal<sup>®</sup> hydroflow lite** or **-Xlite** into the trough so formed for the mixing technique (Fig. 7). At the same time coat the relevant areas intraorally with **Detaseal<sup>®</sup> hydroflow lite** or **-Xlite**. Intra-Oral-Tips or application syringes (Fig. 8) can be used for this purpose. Insert the filled tray into the patient's mouth and position it using light pressure. Then seat the impression without pressure until it has fully set (Fig. 9).


Fig. 4

## 4. Disinfecting

After removing tray from the mouth, rinse impression under running, lukewarm water. If desired, impression can then be disinfected by immersing into 2% glutaraldehyde or other usual solutions for 15 minutes.

Fig. 5

## 5. Pouring the model

We recommend waiting at least 30 minutes before pouring the model. The impression will retain its dimensional stability for up to 14 days. All class III and IV dental plasters and standard modelling resins can be used.

Fig. 6

## 6. Cleaning the tray

After setting, material can be removed with a blunt instrument. Remaining film of **trayloc<sup>®</sup> A** is dissolved by immersion in acetone or similar cleaning agents (be sure to use only in well-ventilated rooms). Tray can then be cleaned and disinfected in the usual manner.

Fig. 7

## Important working hints

- Where movement is impeded or blocked, discard the cartridge. Do not apply force.
- Substances such as retraction liquids can adversely affect the setting reaction and should therefore be thoroughly rinsed off before taking the impression.
- Do not combine with condensation curing silicones.
- Contact with latex examination gloves or other latex material may impair the setting reaction of **Detaseal<sup>®</sup> hydroflow putty**.
- Do not leave any residues of the impression material in the patient's mouth, and prevent inadvertent pressing of material into areas outside the oral cavity through mouth-antrum connections or malformations of the palate.
- Cured impression materials are chemically inert – spots on clothing should be avoided.
- For the one-step putty wash technique we recommend **Detaseal<sup>®</sup> hydroflow lite** or **-Xlite** as precision wash material for single-stage impression-taking.
- To ensure trouble-free use we only recommend using dynamic mixers from DETAX following the mixer manufacturer's instructions
- **Detaseal<sup>®</sup> hydroflow putty** attains its typically kneadable consistency after approx. 30 seconds (tray application).
- Please follow the instructions of the safety data sheet!

Fig. 8


Further information:

Silicone based materials are proven a million times. On condition of a proper application, **undesired effects** are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test **before the application** of the material.

Fig. 9

Only for dental use by qualified personnel.

## Indications for use:

- Impressions with sandwich and double-mix technique
- Situation impressions

## Technical data:

DIN EN ISO 4823 - Type 0

- **Mixed volume:** maxi cartridge: 380 ml base + catalyst
- **Mixing ratio:** 5:1
- **Product colour:** base: green catalyst: white
- **Working time:** approx. 1 min. 30 sec.\*
- **Setting time:** approx. 4 min.\*
- **Time in mouth:** approx. 2 min. 30 sec.\*
- **Strain in compression:** approx. 3,0 %
- **Recovery from deformation:** approx. 99,8 %
- **Linear dimensional change:** approx. 0,2 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity.
- **Storage:** at ambient temperature

\* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard these times.

## Ordering information:

**Detaseal<sup>®</sup> hydroflow putty**  
Jumbo packing **03283**

maxi cartridge of 380 ml, base + catalyst  
10 dynamic mixers, yellow  
1 rotating ring, yellow

**Dynamic mixers** **02833**  
yellow, 40 pcs.

**Rotating ring** **02836**  
yellow, 1 pc.

**Standard pack** **02727**  
2 jars of 250 ml  
2 measuring scoops

**Detaseal<sup>®</sup> hydroflow heavy**  
Jumbo packing **02729**

maxi cartridge of 380 ml, base + catalyst  
10 dynamic mixers, yellow  
1 rotating ring, yellow

**Detaseal<sup>®</sup> hydroflow mono**  
Jumbo packing **03373**

maxi cartridge of 380 ml, base + catalyst  
10 dynamic mixers, yellow  
1 rotating ring, yellow

**Standard packing** **02754**  
2 cartridges of 50 ml  
6 mixing cannulas, green

**Multi packing** **02757**  
4 cartridges of 50 ml  
12 mixing cannulas, green

**Mixing cannulas** **02705**  
green, 48 pcs.

**Detaseal<sup>®</sup> hydroflow lite**  
Standard packing

2 cartridges of 50 ml  
12 mixing cannulas, yellow  
**regular set** **02730**  
**fast set** **02733**

**Detaseal<sup>®</sup> hydroflow Xlite**  
Standard packing

2 cartridges of 50 ml  
12 mixing cannulas, yellow  
**regular set** **02741**  
**fast set** **02744**

**trayloc<sup>®</sup> A** **03098**  
17 ml bottle with brush  
15 ml spray bottle **03392**

**DETAX**

GmbH & Co. KG  
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany  
Telefon: 072 43/510-0 · Fax: 072 43/510-100  
www.detax.de · post@detax.de

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personnel).

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