Detaseal® hydroflow putty

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15 ml botella de aerosol

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Detaseal® hydroflow putty

VPS precision impression silicone,
kneadable consistency – fast set

Fig. 1

1. Tray preparation

Select a tray type suited for desired application (standard-size or custom, in metal or hard plastics). To avoid teeth contact with the tray, place stops from wax or light-curing resin into the tray at non-relevant areas of the impression. For optimal adhesion of impression material, we recommend applying a thin layer of **trayloc® A** to tray (Fig. 1). Let dry for approx. 3 minutes.

Fig. 2

2. Mixing and dispensing

Dispense in a mixing ratio of 1:1, using the enclosed measuring scoops. Mix base and catalyst thoroughly by manual kneading until a uniform colour is achieved (Fig. 2).

Fig. 3

3. Using

First apply the mixed **Detaseal® hydroflow putty** onto the impression tray, giving the teeth area a distinctly trough-like shape (Fig. 3). Then apply a layer of **Detaseal® hydroflow lite** or **-Xlite** into the shaped depression (Fig. 4). Insert the filled tray into patient's mouth and position it within 3 – 4 seconds, under light pressure. Then seat the impression without pressure, until **Detaseal® hydroflow** has set.


Fig. 4

4. Disinfecting

After removing tray from the mouth, rinse impression under running, lukewarm water. If desired, impression can then be disinfected by immersing into 2% glutaraldehyde or other usual solutions for 15 minutes.

Fig. 5

5. Pouring the model

We recommend waiting at least 30 minutes before pouring the model. The impression will retain its dimensional stability for up to 14 days. All class III and IV dental plasters and standard modelling resins can be used.

6. Cleaning the tray

After setting, material can be removed with a blunt instrument. Remaining film of **trayloc® A** is dissolved by immersion in acetone or isopropanol (be sure to use only in well-ventilated rooms). Tray can then be cleaned and disinfected in the usual manner.


Important working hints

- Make sure that only base and catalyst components from the same lot number are used.
- Substances such as retraction liquids can adversely affect the setting reaction and should therefore be thoroughly rinsed off before taking the impression.
- Do not combine with condensation curing silicones.
- Contact with latex examination gloves or other latex material may impair the setting reaction of **Detaseal® hydroflow putty**.
- Do not leave any residues of the impression material in the patient's mouth, and prevent inadvertent pressing of material into areas outside the oral cavity through mouth-antrum connections or malformations of the palate.
- Cured impression materials are chemically inert – spots on clothing should be avoided.
- We recommend use of **Detaseal® hydroflow lite** or **-Xlite** as correction material in the double mix technique for the one time impression.
- In order to enhance easy mouth removal in case of one-step impressions and to avoid blocking in strong undercuts, fill undercuts beforehand with **Detaseal® antilock**.
- Please follow the instructions of the safety data sheet!

Further information:

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

Only for dental use by qualified personnel.

Indications for use:

- Preliminary impressions for the double impression technique
- Impressions with sandwich and double-mix technique
- Situation impressions

Technical data:

DIN EN ISO 4823 - Type 0

- **Mixed volume:** 2 x 250 ml jars
- **Mixing ratio:** 1:1
- **Product colour:** base: green catalyst: white
- **Mixing time:** approx. 30 sec.
- **Working time:** approx. 1 min. 30 sec.*
- **Setting time:** approx. 3 min. 30 sec.*
- **Time in mouth:** approx. 2 min.
- **Strain in compression:** approx. 3 %
- **Recovery from deformation:** approx. 99,8 %
- **Linear dimensional change:** approx. 0,2 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity.
- **Storage:**

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperature accelerate, decreased temperatures retard a. m. times.

Ordering information:

Detaseal® hydroflow putty
Standard pack **02727**
2 jars of 250 ml
2 measuring scoops

Multi packing **02728**
4 jars of 250 ml
2 measuring scoops

Detaseal® hydroflow lite
Standard packing
2 cartridges of 50 ml
12 mixing cannulas, yellow
regular set **02730**
fast set **02733**

Detaseal® hydroflow Xlite
Standard packing
2 cartridges of 50 ml
12 mixing cannulas, yellow
regular set **02741**
fast set **02744**

Detaseal® antilock
Standard pack **03535**
1 x 50 ml cartridge
6 Mixing cannulas, yellow
6 Intra-oral tips, yellow

trayloc® A
17 ml bottle with brush **03098**
15 ml spray bottle **03392**

Mixing cannulas
yellow, 48 pcs. **02706**

DETAX

GmbH & Co. KG
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
Telefon: 07243/510-0 · Fax: 07243/510-100
www.detax.de · post@detax.de

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personnel).

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