
# **FREEPRINT<sup>®</sup>**

## **3D PRINTING MATERIAL**

### **TEMPORARY CROWNS & BRIDGES**

# **DETAX**

<b>DE</b>	Gebrauchsanweisung .....	3	<b>NO</b>	Bruksanvisning .....	36
<b>EN</b>	Instructions for use .....	6	<b>PL</b>	Instrukcja stosowania .....	39
<b>FR</b>	Mode d'emploi .....	9	<b>PT</b>	Manual de instruções.....	42
<b>ES</b>	Modo de empleo .....	12	<b>RO</b>	Instrucțiuni de utilizare .....	45
<b>IT</b>	Istruzioni per l'uso .....	15	<b>SV</b>	Bruksanvisning .....	48
<b>DA</b>	Brugsanvisning.....	18	<b>SK</b>	Návod na použitie.....	51
<b>EL</b>	Οδηγίες χρήσης.....	21	<b>SL</b>	Navodila za uporabo .....	54
<b>HR</b>	Upute za uporabu .....	24	<b>CS</b>	Návod k použití .....	57
<b>LV</b>	Lietošanas instrukcija.....	27	<b>HU</b>	Használati útmutató.....	60
<b>LT</b>	Naudojimo instrukcija .....	30	<b>RU</b>	Инструкция по применению .....	63
<b>NL</b>	Gebruiksaanwijzing .....	33	<b>TR</b>	Kullanım kılavuzu .....	66

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Temporary temps and bridges

## **PATIENT TARGET GROUP**

Persons being treated in the context of a dental procedure.

## **INTENDED USERS**

Dentist, dental technician

## **SUITABLE FOR THE FOLLOWING DLP PRINTERS/CLEANING/POST-EXPOSURE**

see "Annex 1" (enclosed separately)

## **PROCESSING**

- ▶ The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (observe process description on page 8).
- ▶ Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- ▶ Do not use heat-based methods for disinfection or sterilisation. This could possibly deform the workpiece.
- ▶ **FREEPRINT® TEMP** is applied for the production of temporary c & b in the anterior and posterior area and can be used for single temps and bridges up to 4 elements containing one pontic. The occlusal minimum thickness amounts to 1.5 mm, circularly up to 0.8 mm. The connector cross-section for pontics in the anterior area is 12 mm<sup>2</sup> minimum, in the posterior area 15 mm<sup>2</sup>. Any non-compliance of these requirements can lead to an undesirable outcome.
- ▶ A layer thickness of 50 µm is suitable for the printing process.
- ▶ **FREEPRINT® TEMP** can be shaded individually by means of our smartrepair® System.
- ▶ Polish the surface mechanically. Prepolish by means of rotating brushes and prepolishing paste, highgloss by means of buffing wheels and highshine polishing paste for resins.
- ▶ The temporaries can be cemented by means of usual temporary cements, such as tempolink® (eugenol-free). Temporary cements containing eugenol must be removed completely after the removal of the provisional as otherwise the setting of subsequently used permanent cements could be affected.
- ▶ Processing temperature 23 °C ± 2 °C.

## **SAFETY INFORMATION**

- ▶ Only for the specified intended use by trained specialists.
- ▶ Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant / breastfeeding women. Irritating to airways, eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- ▶ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ▶ After contact with skin wash immediately with water and soap.
- ▶ Biocompatibility is only guaranteed with complete polymerisation.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

## **NOTES**

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly sealed, immediately close the container carefully after each use.

- ▶ Please follow the instructions on the safety data sheet!

#### For users and/or patients:

Any serious incidents occurring in relation to this product should be reported immediately to [incident@detax.com](mailto:incident@detax.com) and to the competent authority of the Member State in which the user and/or patient is established.

### **STORAGE**

- ▶ Store **FREEPRINT® TEMP** in a dry place (at 15 °C – 28 °C) and away from light. Even slight exposure to light can trigger polymerization.
- ▶ To protect against contamination, cover the material in the tray with the lid or a glass plate.

### **CONTRAINDICATION**

Contains (meth)acrylics and phosphine oxides.

Some ingredients of **FREEPRINT® TEMP** may cause allergic reactions in predisposed persons. In such cases refrain from using the product. **FREEPRINT® TEMP** only insert intraorally in completely polymerised state.

### **ADVERSE EFFECTS**

Product may cause allergic reactions.

### **DISPOSAL**

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

#### Storage:

#### Application:

At 23 °C ± 2 °C

# Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

## Construction process

Generation of a Print Job complying with machine and material parameters

## Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

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+49 7243/510-0 · +49 7243/510-100

[www.detax.com](http://www.detax.com) · [post@detax.com](mailto:post@detax.com)

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personnel).

IM05

Made in  
Germany

07/2023

5

# Annex 1

## DIGITAL LIGHT PROCESSING (DLP) PRINTER, OPERATION SOFTWARE AND PARAMETER

Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
Asiga Max 	385 nm	6.1 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint temp 385_5
Asiga Pico 2 	385 nm	20 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint temp 385_5
Asiga PRO 2 	385 nm	5.7 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint temp 385_5
Asiga PRO 4K 	385 nm	7.0 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint temp 385_5
Ivoclar PrograPrint PR5 	388 nm	16 mW/cm <sup>2</sup>	PrograPrint CAM 1.1.10.1	Detax FREEPRINT temp
Microlay Versus 	385 nm	4.3 mW/cm <sup>2</sup>	Microform 1.0.3.7	DETAX Freeprint Temp 385 50 microns v5.3
Micraft Prime Series 	385 nm	6.2 mW/cm <sup>2</sup>	Utility 6.3.0	Detax_Freeprint temp 385_50_PH / 100_PH
Mircraft Ultra Series 	385 nm	5.7 mW/cm <sup>2</sup>	Utility 6.3.0	Detax_Freeprint temp UV_50 / 100
Rapidshape D10+/D20+/D30+ 	385 nm	2.0 mW/cm <sup>2</sup>	Netfabb 2020	DETAX Freeprint-temp 385

## DIGITAL LIGHT PROCESSING (DLP) PRINTER, OPERATION SOFTWARE AND PARAMETER

Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
Rapidshape D70+/D90+ 	385 nm	5.0 mW/cm <sup>2</sup>	Netfabb 2020	DETX Freeprint-temp 385
Way2Production SolFlex Series 	385 nm	8.0 mW/cm <sup>2</sup>	Netfabb 2020	Freeprint temp 385 04-05-2021

\*The set of parameters includes all relevant material- and printer specific information

## CLEANING EQUIPMENT

Cleaning unit Manufacturer, Model	Cleaning process
Ivoclar PrograPrint Clean	Clean the parts with isopropyl alcohol (purity $\geq$ 98 %) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq$ 98 %) for 3 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Rapidshape RS Wash	Use the following settings: DETX Freeprint-temp uv Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Ultrasonic bath Bandelin Sonorex	Clean the parts with isopropyl alcohol (purity $\geq$ 98 %) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq$ 98 %) for 3 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.

## LIGHT CURING EQUIPMENT

Light Curing unit Manufacturer, Model	Curing process
Dentalfarm Photopol	2x3 min, progressiv + N <sub>2</sub> , turnaround components after 3 min
Ivoclar PrograPrint Cure	Post curing A: Wavelength = 405 nm; Intensity = 100 %; Duration = 120 s Post curing B: Wavelength = 460 nm; Intensity = 100 %; Duration = 120 s
NK Optik Otofash G171	2 x 2000 flashes under inert gas, turnaround components after 2000 flashes
NK Optik Otofash 250/500	4000 flashes under inert gas @15 Hz
Rapidshape RS Cure	Use the following settings: DETX Freeprint-temp uv

**DETX**

DETX GmbH

Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany

+49 7243/510-0 · +49 7243/510-100

www.detax.com · post@detax.com

Made in  
Germany

09/2023  
3

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DETAX GmbH