# FREEPRINT®

## 3D PRINTING MATERIAL

INDIVIDUAL  
FUNCTIONAL TRY-INS

# DETAX

<b>DE</b>	Gebrauchsanweisung .....	3	<b>PT</b>	Manual de instruções.....	42
<b>EN</b>	Instructions for use .....	6	<b>RO</b>	Instrucțiuni de utilizare .....	45
<b>FR</b>	Mode d'emploi .....	9	<b>SV</b>	Bruksanvisning .....	48
<b>ES</b>	Modo de empleo.....	12	<b>SK</b>	Návod na použitie.....	51
<b>IT</b>	Istruzioni per l'uso.....	15	<b>SL</b>	Navodila za uporabo .....	54
<b>DA</b>	Bruksanvisning.....	18	<b>CS</b>	Návod k použití .....	57
<b>EL</b>	Οδηγίες χρήσης.....	21	<b>HU</b>	Használati útmutató.....	60
<b>HR</b>	Upute za uporabu .....	24	<b>RU</b>	Инструкция по применению .....	63
<b>LV</b>	Lietošanas instrukcija.....	27	<b>TR</b>	Kullanım kılavuzu .....	66
<b>LT</b>	Naudojimo instrukcija .....	30	<b>FI</b>	Käyttöohje.....	69
<b>NL</b>	Gebruiksaanwijzing .....	33	<b>ET</b>	Kasutusjuhend .....	72
<b>NO</b>	Bruksanvisning .....	36	<b>BG</b>	Инструкция за употреба.....	75
<b>PL</b>	Instrukcja stosowania .....	39			

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Resin for dental 3D printing

## **INDICATION**

Individual functional try-ins of complete dentures

## **PATIENT TARGET GROUP**

Persons being treated in the context of a dental procedure.

## **INTENDED USERS**

Dentist, dental technician

## **SUITABLE FOR THE FOLLOWING DLP PRINTERS**

"Annex 1" (enclosed separately)

## **PROCESSING**

- ▶ The properties of the final product depend, among other things, on the finishing process. Correct post-exposure is important for biocompatibility. Therefore, it must be ensured that the exposure unit is in proper condition and the molded parts are completely cured (refer to the process description).
- ▶ Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- ▶ A layer thickness of 100 µm is suitable for the printing process.
- ▶ Polish the surface mechanically if necessary. Pre-polishing is done with rotating brushes and pre-polishing paste, high-gloss polishing with buffing and high-gloss polishing agents for resin.
- ▶ Processing temperature 23 °C ± 2 °C.

## **SAFETY INFORMATION**

- ▶ Only for the specified intended use by trained specialists.
- ▶ Avoid direct contact with the liquid material and the components before post-curing, in particular in pregnant / breastfeeding women. Irritating to airways, eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- ▶ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ In case of contact with skin, wash immediately with plenty of soap and water.
- ▶ Biocompatibility is guaranteed only after complete polymerization.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

## **NOTES**

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly closed, carefully close immediately after each use.
- ▶ Read and understand the safety data sheet!

### **For users and/or patients:**

Any serious incident that has occurred in connection with this product must be reported immediately to [incident@detax.com](mailto:incident@detax.com) and to the competent public authority of the Member State in which the user and/or patient resides.

## **STORAGE**

- ▶ Store **FREEPRINT® tryin** in a dry place (at 15 °C – 28 °C) and away from light. Even slight exposure to light can trigger polymerization.
- ▶ To protect against contamination, cover the material in the tray with the lid or a glass plate.

## **CONTRAINDICATION**

Contains (meth)acrylates and phosphine oxides.

Ingredients of **FREEPRINT® tryin** may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further. Apply **FREEPRINT® tryin** intraorally only when fully polymerized.

## **SIDE EFFECTS**

Product may cause allergic reactions.

## **DISPOSAL**

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

### **Storage:**

### **Processing:**

at 23 °C ± 2 °C

# Manufacturing process

Data processing and creation of the support structure according to the specifications of the CAD software manufacturer

## Building process

Creation of a print job in compliance with the machine and material parameters.

## Finishing process

After raising the platform, a dripping time of approx. 10 minutes is recommended. If possible, the finishing should take place immediately after the building process.

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# Annex 1

## DIGITAL LIGHT PROCESSING (DLP) PRINTER, OPERATION SOFTWARE AND PARAMETER

Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
<b>Asiga Max</b> 	385 nm	6.1 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint tryin 385_5
<b>Asiga Pico 2</b> 	385 nm	20 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint tryin 385_5
<b>Asiga PRO 2</b> 	385 nm	5.7 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint tryin 385_5
<b>Asiga PRO 4K</b> 	385 nm	7.0 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint tryin 385_5
<b>Microlay Versus</b> 	385 nm	4.3 mW/cm <sup>2</sup>	Microform 1.0.3.7	DETAX Freeprint Tryin 385 100 microns v5.3
<b>Miicraft Ultra Series</b> 	385 nm	5.7 mW/cm <sup>2</sup>	Utility 6.3.0	Detax_Freeprint tryin 385_100

## DIGITAL LIGHT PROCESSING (DLP) PRINTER, OPERATION SOFTWARE AND PARAMETER

Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
<b>Rapidshape D20+/D30+/D40+</b> 	385 nm	2.0 mW/cm <sup>2</sup>	Netfabb 2020	DETAX Freeprint-tryin 385
<b>Way2Production SolFlex Series</b> 	385 nm	8.0 mW/cm <sup>2</sup>	Netfabb 2020	Freeprint tryin 385 24-05-2022

\*The set of parameters includes all relevant material- and printer specific information

## CLEANING EQUIPMENT

Cleaning unit Manufacturer, Model	Cleaning process
<b>Rapidshape RS Wash</b>	Use the following settings: DETAX Freeprint-tryin 385 Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
<b>Ultrasonic bath Bandelin Sonorex</b>	Clean the parts with isopropyl alcohol (purity $\geq$ 98 %) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq$ 98 %) for 3 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.

## LIGHT CURING EQUIPMENT

Light Curing unit Manufacturer, Model	Curing process
<b>Dentalfarm Photopool</b>	2 x 3 min, progressive + N <sub>2</sub> , turnaround components after 3 min
<b>NK Optik Otoflash G171</b>	2 x 2000 flashes under inert gas, turnaround components after 2000 flashes
<b>NK Optik Otoflash 250/500</b>	4000 flashes under inert gas @ 15 Hz
<b>Rapidshape RS cure</b>	Use the following settings: DETAX Freeprint-tryin 385

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DETAX GmbH