
# **FREEPRINT<sup>®</sup>**

## **3D PRINTING MATERIAL**

DENTAL MED RESIN

» **CROWN** «

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Composite for dental 3D printing

## INDICATION

Permanent single crowns (anterior and posterior areas), permanent inlays and onlays, veneers, long-term provisional crowns, bridges and permanent denture teeth.

## PATIENT TARGET GROUP

Persons being treated in the context of a dental procedure.

## INTENDED USERS

Dentist, dental technician

## SUITABLE FOR THE FOLLOWING DLP PRINTERS/CLEANING/POST-EXPOSURE

see "Annex 1" (enclosed separately)

## PROCESSING

- ▶ The properties of the final product depend, among other things, on the finishing process. Correct post-exposure is important for biocompatibility. Therefore, it must be ensured that the exposure unit is in proper condition and the molded parts are completely cured (refer to the process description).
- ▶ Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- ▶ Do not use thermal methods for disinfection or sterilization. This might deform the workpiece.
- ▶ **FREEPRINT® CROWN** is used for the fabrication of permanent single crowns in the anterior and posterior region, denture teeth, inlays and onlays. It is also used for long-term provisional crowns and bridges having up to 3 pontics with an intermediate pontic and veneers. The minimum occlusal thickness is 1.5 mm, circularly 1.5 mm and cervically 1.0 mm. The connection cross-section for intermediate pontics of a bridge must be at least 16 mm<sup>2</sup>. Failure to observe this information may lead to an undesirable result.
- ▶ A layer thickness of 50 µm is suitable for the printing process.
- ▶ Mechanically polish the surface. Pre-polishing is done with rotating brushes and pre-polishing paste, a high-gloss polish with buffing and high-gloss polishing agents for resin composite.
- ▶ The restorations can be used with the commonly used permanent luting cements, e.g. Variolink Esthetic DC by Ivoclar or RelyX Unicem2 by 3M Espe. Observe the manufacturer's instructions.
- ▶ **Attachment of permanent denture teeth**  
Once the denture base (see additional IFU **FREEPRINT® DENTURE**) and the denture teeth have been removed from the respective build platforms and the corresponding cleaning processes completed, leave the printed components in their green body state (not post-exposed). Place the printed denture teeth in the corresponding tooth receptacles of the denture base, check the correct seating of the teeth and remove them again. Apply liquid **FREEPRINT® DENTURE** to the receptacles with a small brush or a syringe, and connect the repositioned denture teeth by exposure in a xenon flash unit (Otoflash G171) with 400 flashes under inert gas conditions (nitrogen). Finally post-cure the complete denture in a xenon flash unit (Otoflash G171) with 2 × 1800 flashes under inert gas conditions (nitrogen).
- ▶ Processing temperature 23 °C ± 2 °C.

## SAFETY INFORMATION

- ▶ Only for the specified intended use by trained specialists.
- ▶ Avoid direct contact with the liquid material and the components before post-curing. Irritating to airways, eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (e.g. protective gloves, goggles) when handling the uncured material.

- ▶ For finishing the cured material, suitable personal protective equipment must be used (protective gloves, goggles, face mask).
- ▶ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ In case of contact with skin, wash immediately with plenty of soap and water.
- ▶ Biocompatibility is ensured only after complete polymerization.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

## **NOTES**

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly closed, carefully close immediately after each use.
- ▶ Read and understand the safety data sheet!

### **For users and/or patients:**

Any serious incident that has occurred in connection with this product must be reported immediately to [incident@detax.com](mailto:incident@detax.com) and to the competent public authority of the Member State in which the user and/or patient resides.

### **Brief Report on Safety and Clinical Performance (SSCP):**

The Summary Safety and Clinical Performance Report (SSCP) for BASIS UDI-DI ++EDET0996K (FREEPRINT® CROWN) is stored in the European Database for Medical Devices (EUDAMED) under the internet link <https://ec.europa.eu/tools/eudamed>.

Until full functionality of the EUDAMED is announced, the SSCP can be requested from [regulatory.affairs@detax.com](mailto:regulatory.affairs@detax.com).

## **STORAGE**

- ▶ Store FREEPRINT® CROWN in a dry place (at 15 °C – 28 °C) and away from light. Even slight exposure to light can trigger polymerization.
- ▶ To protect against contamination, cover the material in the tray with the lid or a glass plate.

## **CONTRAINDICATION**

Contains (meth)acrylates and phosphine oxides.

Ingredients of FREEPRINT® CROWN may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further. Apply FREEPRINT® CROWN intraorally only when fully polymerized.

## **SIDE EFFECTS**

Product may cause allergic reactions.

## **DISPOSAL**

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

# Manufacturing process

Data processing and creation of the support structure according to the specifications of the CAD software manufacturer

## Building process

Creation of a print job in compliance with the machine and material parameters, see table

## Finishing process

After raising the platform, a dripping time of approx. 10 minutes is recommended. If possible, the finishing should take place immediately after the building process.

## Cleaning

see "Annex 1, CLEANING EQUIPMENT"

## Drying

Store the component protected from light at room temperature for 30 minutes before post-exposure.

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

\*The set of parameters includes all relevant material and printer specific information

## CLEANING EQUIPMENT

Cleaning unit Manufacturer, Model	Cleaning process
Ivoclar PrograPrint Clean	Clean the parts with isopropyl alcohol (purity $\geq 98\%$ ) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$ ) for 2 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Rapidshape RS wash	Use the following settings: DETAX Freeprint-crown Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Ultrasonic bath Bandelin Sonorex	Clean the parts with isopropyl alcohol (purity $\geq 98\%$ ) for 1 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$ ) for 1 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.

## LIGHT CURING EQUIPMENT

Light Curing unit Manufacturer, Model	Curing process
Ivoclar PrograPrint Cure	Post curing A: Wavelength = 405 nm; Intensity = 45 %; Duration = 5 min Post curing B: Wavelength = 460 nm; Intensity = 45 %; Duration = 5 min
NK Optik Otofash G171	2 x 2000 flashes under inert gas, turnaround components after 2000 flashes
Rapidshape RS Cure	Use the following settings: DETAX freeprint-crown

# DETAX

DETAX GmbH

Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany

+49 7243/510-0 · +49 7243/510-100

www.detax.com · post@detax.com

Made in  
Germany

09/2023

6.1