**FREEPRINT®**

3D PRINTING MATERIAL

DENTAL TEC RESIN

**CAST**

CASTING TECHNIQUE

**DETAX**

# FREEPRINT® CAST

**DEUTSCH**

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# FREEPRINT® CAST

**ENGLISH**

Light-curing resin based on (meth)acrylate, for cast parts in dental precision casting, burn-out for DLP printers with UV-LED 385 nm

## Area of application

Fabrication of residue-free combustible cast objects for precision casting technology employing the DLP process.

## SUITABLE FOR THE FOLLOWING DLP PRINTERS/CLEANING/POST-EXPOSURE

see "Annex 1" (enclosed separately)

## Processing

- ▶ Phosphate-bonded investment material can be used for fabrication of the casting mould. A preheating temperature of approx. 800 °C, is to be maintained for approx. 30-45 minutes, depending on the size of the casting mould, regardless of the final casting temperature.
- ▶ Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- ▶ The properties of the final product depend, among other things, on post-processing. Therefore it must be ensured that the light unit is in an orderly condition and that the models are completely cured (observe process description).
- ▶ Processing temperature 23 °C ± 2 °C.

## Safety Information

- ▶ Only to be used by trained specialists for the specified area of application.
- ▶ Product contains (meth)acrylates and phosphine oxides. The ingredients of **FREEPRINT® CAST** may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further.
- ▶ Avoid direct contact with the liquid material and the components before post-curing. Irritating to eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ In case of contact with skin, wash immediately with plenty of soap and water.
- ▶ The processing instructions and precautionary measures must be observed.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

## Notes

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly closed, carefully close immediately after each use.
- ▶ Product may cause allergic reactions.
- ▶ Read and understand the safety data sheet!

## Storage

- ▶ **FREEPRINT® CAST** is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.
- ▶ To protect against contamination, cover the material in the tray with the lid or a glass plate.

## Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

## Manufacturing process

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

## Construction process

Generation of a Print Job complying with machine and material parameters

## Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

## Cleaning

see "Annex 1, CLEANING EQUIPMENT"

## Drying

Heat the construction components for 30 minutes in an oven to approx. 40 °C to remove the solvent residues from the cleaning process.

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\*The set of parameters includes all relevant material- and printer specific information

## CLEANING EQUIPMENT

Cleaning unit Manufacturer, Model	Cleaning process
Ultrasonic bath Bandelin Sonorex	Clean the parts with isopropyl alcohol (purity $\geq 98\%$ ) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$ ) for 3 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.

## LIGHT CURING EQUIPMENT

Light Curing unit Manufacturer, Model	Curing process
NK Optik Otoflash G171	2 x 2000 flashes under inert gas, turnaround components after 2000 flashes
NK Optik Otoflash 250/500	4000 flashes under inert gas @15 Hz

# DETAX

DETAX GmbH

Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany

☎ +49 7243/510-0 · 📠 +49 7243/510-100

www.detax.com · post@detax.com

Made in  
Germany

01/2024

3