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# FREEPRINT® GINGIVA

Light-curing resin based on (meth)acrylate, for the generative fabrication of flexible gingiva masks for DLP printers with UV-LED 385 nm

## SUITABLE FOR THE FOLLOWING DLP PRINTERS/CLEANING/POST-EXPOSURE

see "Annex 1" (enclosed separately)

### Processing

- ▶ The properties of the final product depend, among other things, on post-processing. Therefore it must be ensured that the light unit is in an orderly condition and that the models are completely cured (observe process description on page 2).
- ▶ Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- ▶ Processing temperature 23 °C ± 2 °C.

### Safety Information

- ▶ Only to be used by trained specialists for the specified area of application.
- ▶ Product contains (meth)acrylates and phosphine oxides. The ingredients of **FREEPRINT® GINGIVA** may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further.
- ▶ Avoid direct contact with the liquid material and the components before post-curing. Irritating to airways, eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ In case of contact with skin, wash immediately with plenty of soap and water.
- ▶ The processing instructions and precautionary measures must be observed.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

### Notes

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly closed, carefully close immediately after each use.
- ▶ Product may cause allergic reactions.
- ▶ Read and understand the safety data sheet!

### Storage

- ▶ **FREEPRINT® GINGIVA** is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.
- ▶ To protect against contamination, cover the material in the tray with the lid or a glass plate.

### Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

### Area of application:

Gingival masks for dental 3D models

### Processing:

at 23 °C ± 2 °C

### Storage

15 °C  
59 °F

28 °C  
82 °F

### Ordering information:

**FREEPRINT® GINGIVA 385**  
500 g **02820**  
1.000 g **02843**

**FREEPRINT® MODEL 2.0 385**  
500g  
caramel **02844**

1.000 g  
sand **02128**  
light grey **02099**  
grey **02177**  
caramel **02850**  
white **02148**

3 kg Eco Bag  
caramel **04328**  
grey **04329**  
light-grey **04330**  
sand **04331**  
white **04332**

5 kg Eco Bag  
caramel **04015**  
grey **04106**  
light-grey **04107**  
sand **04117**  
white **04118**

**FREEPRINT® MODEL 385**  
1.000 g  
sand **03778**  
ivory **03780**  
grey **03782**  
caramel **02935**

3 kg Eco Bag  
sand **04333**

5 kg Eco Bag  
sand **04321**

# DETAX

DETAX GmbH  
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany  
☎ +49 7243/510-0 · 📠 +49 7243/510-100  
www.detax.com · post@detax.com

Made in  
Germany  
10/2023  
5

## Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

### Construction process

Generation of a Print Job complying with machine and material parameters

### Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

### Cleaning

see "Annex 1, CLEANING EQUIPMENT"

### Main cleaning

see "Annex 1, CLEANING EQUIPMENT"

### Post-curing

see "Annex 1, CURING LIGHT EQUIPMENT"

# Annex 1

## DIGITAL LIGHT PROCESSING (DLP) PRINTER, OPERATION SOFTWARE AND PARAMETER

Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
<b>Asiga Max</b> 	385 nm	6.1 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint gingiva 385_5
<b>Asiga Pico 2</b> 	385 nm	20 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint gingiva 385_5
<b>Asiga PRO 2</b> 	385 nm	5.7 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint gingiva 385_5
<b>Asiga PRO 4K</b> 	385 nm	7.0 mW/cm <sup>2</sup>	Composer 1.2.11	Detax_Freeprint gingiva 385_5
<b>Ivoclar PrograPrint PR5</b> 	388 nm	16 mW/cm <sup>2</sup>	PrograPrint CAM 1.1.10.1	Detax FREEPRINT gingiva
<b>Microlay Versus</b> 	385 nm	4.3 mW/cm <sup>2</sup>	Microform 1.0.3.7	DETAX Freeprint Gingiva 385 100 microns v5.3
<b>Miicraft Prime Series</b> 	385 nm	6.2 mW/cm <sup>2</sup>	Utility 6.3.0	Detax_Freeprint gingiva 385_100_PH
<b>Miicraft Ultra Series</b> 	385 nm	5.7 mW/cm <sup>2</sup>	Utility 6.3.0	Detax_Freeprint gingiva 385_100
<b>Rapidshape D20+ / D30+ / D40+</b> 	385 nm	2.0 mW/cm <sup>2</sup>	Netfabb 2020	DETAX Freeprint-gingiva uv

# DIGITAL LIGHT PROCESSING (DLP) PRINTER, OPERATION SOFTWARE AND PARAMETER

Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
Way2Production SolFlex Series 	385 nm	8.0 mW/cm <sup>2</sup>	Netfab 2020	Freeprint gingiva 385 04-05-2021

\*The set of parameters includes all relevant material- and printer specific information

## CLEANING EQUIPMENT

Cleaning unit Manufacturer, Model	Cleaning process
Ivoclar PrograPrint Clean	Clean the parts with isopropyl alcohol (purity $\geq$ 98 %) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq$ 98 %) for 3 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Rapidshape RS wash	Use the following settings: DETAX Freeprint-gingiva uv Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Ultrasonic bath Bandelin Sonorex	Clean the parts with isopropyl alcohol (purity $\geq$ 98 %) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq$ 98 %) for 3 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.

## LIGHT CURING EQUIPMENT

Light Curing unit Manufacturer, Model	Curing process
Dentalfarm Photopool	2 x 6 min, 120 % + N <sub>2</sub> , turnaround components after 6 min
Formlabs FormCure	30 min
Ivoclar PrograPrint Cure	Post curing A: Wavelength = 405 nm; Intensity = 100 %; Duration = 120 s Post curing B: Wavelength = 460 nm; Intensity = 100 %; Duration = 120 s
Meccatronicore BB Cure	30 min, 85 % Power
NK Optik Otofash G171	2 x 2000 flashes under inert gas, turnaround components after 2000 flashes
NK Optik Otofash 250/500	4000 flashes under inert gas @15 Hz
Phrozen PhrozenCure	30 min
Rapidshape RS cure	Use the following settings: DETAX Freeprint-gingiva uv

DETAX GmbH  
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany  
☎ +49 7243/510-0 · 📠 +49 7243/510-100  
www.detax.com · post@detax.com

Made in  
Germany  
09/2023  
3