FREEPRINT® MODEL 2.0

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FREEPRINT® MODEL 2.0

Light-curing resin based on (meth)acrylate, for the generative fabrication of dental models for DLP printers with UV-LED 385 nm

SUITABLE FOR THE FOLLOWING DLP PRINTERS/CLEANING/POST-EXPOSURE

see "Annex 1" (enclosed separately)

Processing

- ▶ The properties of the final product depend, among other things, on post-processing. Therefore it must be ensured that the light unit is in an orderly condition and that the models are completely cured (observe process description on page 2).
- ▶ Bottle: Before use, the material should be shaken intensively and homogenized with a bottle roller.
- ▶ Eco Bag: Before use, the material should be shaken intensively and homogenized by kneading.
- ▶ Processing temperature 23 °C ± 2 °C.

Safety Information

- ▶ Only to be used by trained specialists for the specified area of application.
- ▶ Product contains (meth)acrylates and phosphine oxides. The ingredients of **FREEPRINT® MODEL 2.0** may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further.
- ▶ Avoid direct contact with the liquid material and the components before post-curing.
Irritating to airways, eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ In case of contact with skin, wash immediately with plenty of soap and water.
- ▶ The processing instructions and precautionary measures must be observed.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

Notes

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly closed, carefully close immediately after each use.
- ▶ Product may cause allergic reactions.
- ▶ Read and understand the safety data sheet!

Storage

- ▶ **FREEPRINT® MODEL 2.0** is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.
- ▶ To protect against contamination, cover the material in the tray with the lid or a glass plate.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations..

Area of application:

3D Print working and presentation models for the entire spectrum of models in dentistry

Processing:

at 23 °C ± 2 °C

Storage

15 °C
59 °F

28 °C
82 °F

Ordering information:

FREEPRINT® MODEL 2.0 385

sand	
1.000 g	02128
3 kg Eco Bag	04331
5 kg Eco Bag	04117

light-grey	
1.000 g	02099
3 kg Eco Bag	04330
5 kg Eco Bag	04107

grey	
1000 g	02177
3 kg Eco Bag	04329
5 kg Eco Bag	04106

caramel	
500 g	02844
1.000 g	02850
3 kg Eco Bag	04328
5 kg Eco Bag	04015

white	
1.000 g	02148
3 kg Eco Bag	04332
5 kg Eco Bag	04118

FREEPRINT® GINGIVA 385	
500 g	02820
1.000 g	02843

FREEPRINT® MODEL T 385	
blue	
1.000 g	02332
3 kg Eco Bag	04334
5 kg Eco Bag	04322

FREEPRINT® MODEL 385	
sand	
1.000 g	03778
3 kg Eco Bag	04333
5 kg Eco Bag	04321

ivory	
1.000 g	03780

grey	
1.000 g	03782

caramel	
1.000 g	02935

DETAX

DETAX GmbH
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
☎ +49 7243/510-0 · 📠 +49 7243/510-100
www.detax.com · post@detax.com

Made in
Germany

10/2023
6

Manufacturing process

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters

Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

Cleaning

see "Annex 1, CLEANING EQUIPMENT"

Post-exposure

see "Annex 1, CURING LIGHT EQUIPMENT"

Annex 1

DIGITAL LIGHT PROCESSING (DLP) PRINTER, OPERATION SOFTWARE AND PARAMETER

Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
Asiga Max 	385 nm	6.1 mW/cm ²	Composer 1.2.11	Detax_Freeprint model 2.0 385 COLOR_5
Asiga Pico 2 	385 nm	20 mW/cm ²	Composer 1.2.11	Detax_Freeprint model 2.0 385 COLOR_5
Asiga PRO 2 	385 nm	5.7 mW/cm ²	Composer 1.2.11	Detax_Freeprint model 2.0 385 COLOR_5
Asiga PRO 4K 	385 nm	7.0 mW/cm ²	Composer 1.2.11	Detax_Freeprint model 2.0 385 COLOR_5
Microlay Versus 	385 nm	4.3 mW/cm ²	Microform 1.0.3.7	DETAX Freeprint Model 2.0 385 COLOR 100 microns v5.3
Miicraft Prime Series 	385 nm	6.2 mW/cm ²	Utility 6.3.0	Detax_Freeprint model 2.0 385_100_PH
Miicraft Ultra Series 	385 nm	5.7 mW/cm ²	Utility 6.3.0	Detax_Freeprint model 2.0 UV_100
Rapidshape D20+/D30+/D40+ 	385 nm	2.0 mW/cm ²	Netfabb 2020	DETAX Freeprint model 2.0 385 COLOR

DIGITAL LIGHT PROCESSING (DLP) PRINTER, OPERATION SOFTWARE AND PARAMETER

Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
 Way2Production SolFlex Series	385 nm	8.0 mW/cm ²	Netfabb 2020	Freeprint model 2.0 385 04-05-2021 Freeprint model 2.0 385 white 05.10.2021

*The set of parameters includes all relevant material- and printer specific information. The placeholder **COLOR** shall be replaced by **caramel, grey, light grey, sand** or **white**

CLEANING EQUIPMENT

Cleaning unit Manufacturer, Model	Cleaning process **
Rapidshape RS Wash	Use the following settings: DETAX Freeprint model 2.0 385 COLOR Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Ultrasonic bath Bandelin Sonorex	Clean the parts with isopropyl alcohol (purity ≥ 98 %) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity ≥ 98 %) for 3 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.

The placeholder **COLOR shall be replaced by **caramel, grey, light grey, sand** or **white**

LIGHT CURING EQUIPMENT

Light Curing unit Manufacturer, Model	Curing process ***
Dentalfarm Photopool	2 x 6 min, 120 % + N ₂ , turnaround components after 6 min
Formlabs FormCure	30 min
Meccatronicore BB Cure	30 min, 85 % Power
NK Optik Otofflash G171	2 x 2000 flashes under inert gas, turnaround components after 2000 flashes
NK Optik Otofflash 250/500	4000 flashes under inert gas @15 Hz
Phrozen PhrozenCure	30 min
Rapidshape RS Cure	Use the following settings: DETAX Freeprint model 2.0 385 COLOR

***The placeholder **COLOR** shall be replaced by **caramel, grey, light grey, sand** or **white**