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FREEPRINT® model WW

Light-curing resin based on (meth)acrylate, for the generative fabrication of dental models for DLP printers with UV-LED 385 nm

Suitable for the following DLP-Printer

Asiga 385 nm	MiiCraft 385 nm	Rapidshape 385 nm	W2P 385 nm
MAX UV Pro2 UV Pico2 UV	MiiCraft y-Serie	DII-Serie	Solflex Serie

Printers may only be operated using material parameters authorised by DETAX.

Processing

- ▶ The properties of the final product depend, among other things, on post-processing. Therefore it must be ensured that the light unit is in an orderly condition and that the models are completely cured (observe process description on page 2).
- ▶ After storage, the material in the bottle should be shaken intensively and homogenized with a bottle roller before use.
- ▶ Processing temperature 23 °C ± 2 °C.

Safety Information

- ▶ Only to be used by trained specialists for the specified area of application.
- ▶ Product contains (meth)acrylates and phosphine oxides. The ingredients of FREEPRINT® model WW may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further.
- ▶ Avoid direct contact with the liquid material and the components before post-curing. Irritating to eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ In case of contact with skin, wash immediately with plenty of soap and water.
- ▶ The processing instructions and precautionary measures must be observed.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

Notes

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly closed, carefully close immediately after each use.
- ▶ Product may cause allergic reactions.
- ▶ Read and understand the safety data sheet!

Storage

FREEPRINT® model WW is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations..

Area of application:

light curing resin for the generative fabrication of dental models for the thermoforming technique“ cleaning process with water

Processing:

at 23 °C ± 2 °C

Storage

Ordering information:

FREEPRINT® model WW 385
1.000 g
blue-transparent **03105**

FREEPRINT® model 2.0 385
1.000 g
sand **02128**
light-grey **02099**
grey **02177**
caramel **02850**
white **02148**

FREEPRINT® gingiva 385
500 g **02820**
1.000 g **02843**

FREEPRINT® model T 385
1.000 g
blue **02332**

FREEPRINT® model 385
1.000 g
sand **03778**
ivory **03780**
grey **03782**

FREEPRINT® model basic 385
1.000 g
light-grey **02111**

More information at
www.detax.de

Manufacturing process

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters

Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with water for 3 min. in an ultrasonic bath.

Cleaning

Then thoroughly clean the openings, cavities and gap areas, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

Main cleaning process

The main cleaning process is performed in a separate vessel with fresh water for 3 min. in an ultrasonic bath. Prior to post-exposure, check the openings and additional bore holes for residues. Then blow off with compressed air.

Post-exposure

Post-exposure is performed with a xenon photoflash unit (e.g. Otofash G171) with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.