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Light-curing resin based on (meth)acrylate, for the generative fabrication of dental models for DLP printers with UV-LED 385 nm

## Suitable for the following DLP-Printer

Asiga 385 nm	MiiCraft 385 nm	Rapidshape 385 nm	W2P 385 nm
MAX UV Pro2 UV Pico2 UV	MiiCraft y-Serie	DII-Serie	Solflex Serie

Printers may only be operated using material parameters authorised by DETAX.

## Processing

- ▶ The properties of the final product depend, among other things, on post-processing. Therefore it must be ensured that the light unit is in an orderly condition and that the models are completely cured (observe process description on page 2).
- ▶ After storage, the material in the bottle should be shaken intensively and homogenized with a bottle roller before use.
- ▶ Processing temperature 23 °C ± 2 °C.

## Safety Information

- ▶ Only to be used by trained specialists for the specified area of application.
- ▶ Product contains (meth)acrylates and phosphine oxides. The ingredients of **FREEPRINT® model basic** may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further.
- ▶ Avoid direct contact with the liquid material and the components before post-curing. Irritating to eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ In case of contact with skin, wash immediately with plenty of soap and water.
- ▶ The processing instructions and precautionary measures must be observed.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

## Notes

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly closed, carefully close immediately after each use.
- ▶ Product may cause allergic reactions.
- ▶ Read and understand the safety data sheet!!

## Storage

FREEPRINT® model basic is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

## Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

## Area of application:

3D Print working and presentation models for the entire spectrum of models in dentistry

## Processing:

at 23 °C ± 2 °C

## Storage

15 °C  
59 °F

28 °C  
82 °F

## Ordering information:

<b>FREEPRINT® model basic 385</b>	
1.000 g	
light grey	02111
beige	02068
<b>FREEPRINT® model 385</b>	
1.000 g	
sand	03778
ivory	03780
grey	03782
caramel	02850
<b>FREEPRINT® model 405</b>	
1.000 g	
sand	03065
ivory	03779
grey	03781
<b>FREEPRINT® model 2.0 385</b>	
1.000 g	
sand	02128
light-grey	02099
grey	02177
caramel	02850
white	02148
500 g	
caramel	02844
<b>FREEPRINT® model T 385</b>	
1.000 g	
blue	02332
<b>FREEPRINT® gingiva 385</b>	
500 g	02820
1.000 g	02843

## Manufacturing process

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

## Construction process

Generation of a Print Job complying with machine and material parameters

## Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

## Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity  $\geq 98\%$ ) for 3 min. in an ultrasonic bath.

## Cleaning

Then thoroughly clean the openings, cavities and gap areas, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

## Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity  $\geq 98\%$ ) for 3 min. in an ultrasonic bath. Prior to drying, check the openings and additional bore holes for residues.

## Drying

Heat the construction components for 30 min. in an oven to approx. 40 °C to remove the solvent residues from the cleaning process.

## Post-exposure

Post-exposure is performed with a xenon photoflash unit (e.g. Otofash G171) with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.