
FREEPRINT®

3D PRINTING MATERIAL

MODELS

DETAX

FREEPRINT® MODEL

DEUTSCH

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FREEPRINT® MODEL

ENGLISH

Light-curing resin based on (meth)acrylate, for the generative fabrication of dental models for DLP printers with UV-LED 385 nm

Area of application

3D Print working and presentation models for the entire spectrum of models in dentistry

SUITABLE FOR THE FOLLOWING DLP PRINTERS/CLEANING/POST-EXPOSURE

see "Annex 1" (enclosed separately)

Processing

- The properties of the final product depend, among other things, on post-processing. Therefore it must be ensured that the light unit is in an orderly condition and that the models are completely cured (observe process description on page 2).
- Bottle: Before use, the material should be shaken intensively and homogenized with a bottle roller.
- Eco Bag: Before use, the material should be shaken intensively and homogenized by kneading.
- Processing temperature 23 °C ± 2 °C.

Safety Information

- Only to be used by trained specialists for the specified area of application.
- Product contains (meth)acrylates and phosphine oxides. The ingredients of **FREEPRINT® MODEL** may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further.
- Avoid direct contact with the liquid material and the components before post-curing. Irritating to airways, eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- In case of contact with skin, wash immediately with plenty of soap and water.
- The processing instructions and precautionary measures must be observed.
- Refer to the relevant safety data sheet for hazard and safety information.

Notes

- DETAX does not accept liability for any damage caused by misuse.
- Always keep container tightly closed, carefully close immediately after each use.
- Product may cause allergic reactions.
- Read and understand the safety data sheet!

Storage

- FREEPRINT® MODEL** is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.
- To protect against contamination, cover the material in the tray with the lid or a glass plate.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Manufacturing process

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters

Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

Cleaning

see "Annex 1, CLEANING EQUIPMENT"

Drying

Heat the construction components for 30 minutes in an oven to approx. 40 °C to remove the solvent residues from the cleaning process.

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*The set of parameters includes all relevant material- and printer specific information. The placeholder **COLOR** shall be replaced by **grey, ivory or sand**

CLEANING EQUIPMENT

Cleaning unit Manufacturer, Model	Cleaning process**
Rapidshape RS Wash	Use the following settings: DETAX Freeprint-model-COLOR uv Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Ultrasonic bath Bandelin Sonorex	Clean the parts with isopropyl alcohol (purity $\geq 98\%$) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 3 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.

The placeholder **COLOR shall be replaced by **grey, ivory or sand**

LIGHT CURING EQUIPMENT

Light Curing unit Manufacturer, Model	Curing process***
Dentalfarm Photopol	2 x 6 min, 120% + N ₂ , turnaround components after 6 min
Formlabs FormCure	30 min
Meccatronicore BB Cure	30 min, 85 % Power
NK Optik Otoflash G171	2 x 2000 flashes under inert gas, turnaround components after 2000 flashes
NK Optik Otoflash 250/500	4000 flashes under inert gas @15 Hz
Phrozen PhrozenCure	30 min
Rapidshape RS Cure	Use the following settings: DETAX Freeprint-model-COLOR uv

***The placeholder **COLOR** shall be replaced by **grey, ivory or sand**

DETX

DETX GmbH
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
☎ +49 7243/510-0 · 📧 +49 7243/510-100
www.detax.com · post@detax.com

Made in
Germany
01/2024
3