

DENTAL 2023



MEDICAL MATERIALS

DETAX
HIGHEND MEDICAL MATERIALS

TRUST & INTEGRITY
INTERCULTURAL COMPETENCE
ENTHUSIASM FOR POLYMER CHEMISTRY

WERK STOFF 4.0

ADVANCED, INNOVATIVE MATERIALS ARE
THE DRIVING FORCE OF PRODUCT DEVELOPMENT
IN THE FIELD OF MEDICAL DEVICES

| | | | |
|---|--|---|--|
| <p>We develop biocompatible silicones and high-performance polymers for medical applications and 3D printing.</p> |  |  | <p>Our Business relationships are build upon appreciation, integrity and inter-cultural expertise.</p> |
|  | <p>Our enthusiasm and innovative strength lies in every one of the products you will get to know over the next pages</p> | <p>State-of-the-art warehouse logistics and lean IT processes support the global flow of goods and consistent availability.</p> |  |
|  |  |  | <p>The unique physical properties of polymer compounds offer nearly unlimited possibilities for tomorrow's products.</p> |
|  | <p>Modern materials are the crucial driver for the development of novel medical devices.</p> |  |  |

CONTENTS

6 - 7 COMMUNICATION & SERVICE

9 USER MATRIX A-SILICONES

20 USER MATRIX C-SILICONES

50 USER MATRIX 3D RESINS

11 - 72 PRODUCTS


73 INDEX



IMPRESSION TAKING

8

A-SILICONE
BITE REGISTRATION
C-SILICONE



SOFT RELINING

31


COLD CURING MATERIAL
HEAT CURING MATERIAL
SILICONE GLOSS VARNISH
ACCESSORIES



PREVENTION & REGENERATION

43

SEALING OF ROOT DENTINE
FISSURE SEALING
MOUTH RINSING SOLUTION



3D RESINS

49

TOTAL DENTURES
CROWNS & BRIDGES
ORTHODONTIC BASE COMPONENTS
BRACKET TRANSFER TEMPLATES
SPLINTS
TRAYS & PLATES
CASTINGS
MODEL PRODUCTION
GINGIVAL MASK

26

ACCESSORIES IMPRESSION TAKING

AUTOMIX
INLAY TRAYS
TRAY ADHESIVE
ACCESSORIES



35

TEMPORARY TREATMENT

C & B COMPOSITE
FIXATION
IMPLANT CEMENT
REPAIR
SEALING



47

LIQUIDS, SOLUTIONS

SPRAY SOLUTION



61

LABORATORY

LIGHT CURING POLYMERS
PUTTY
ISULATION
GINGIVAL MASK
RELINING
RESOURCES



SERVICES

Do you have any questions regarding our products, a shipment or another matter? You can reach our service-team Monday to Friday between 8 a.m. and 5 p.m.

Do you have technical questions about a certain product, an application technique or chemistry? Our technicians are happy to assist you.

At our Website www.detax.com, you will find even more detailed product information, international events, contact persons, news and many more useful information.
Now NEW: Live Chat with experts at www.detax.com

Do you need documents, such as studies, specific information, MSDS, declarations of conformity or instructions for use? All documentations are available for you.



CUSTOMER SERVICE

service@detax.com
Phone +49 7243 - 510-103
Phone +49 7243 - 510-104
Fax +49 7243 - 510-100

PRODUCT APPLICATION

application@detax.com
3Dapplication@detax.com
Phone +49 7243 - 510-146
Fax +49 7243 - 510-100

REGULATORY AFFAIRS

regulatory.affairs@detax.com
Phone +49 7243 - 510-268
Fax +49 7243 - 510-243

MDR

medi.guide@detax.com
Phone +49 7243 - 510-208
Fax +49 7243 - 510-222

EXPERT LIVE CHAT

www.detax.com



CERTIFICATION



We support our customers in all matters pertaining to medical devices. Our specialists have long-standing experience in the interpretation and implementation of international certification standards and individual country requirements. We will give you competent advice on required registrations for the market access of DETAX brands, OEM or custom-made materials. We will supply your quality management with the required documentation (e.g. ISO certificates, declarations of conformity, product labelling, safety data sheets or technical data sheets, etc.).



Certified according to the guidelines for medical devices (already since 1996!) as well as current QMS standards. Certification according to Annex IX - Chapter I, Regulation (EU) 2017/745, and Directive 93/42/EEC; DIN EN ISO 13485:2016 (also Taiwan) and MDSAP für Canada, Brazil, Australia, Japan, USA as well as GOST R for Russia, GOST B for Belarus. The requirements for the biocompatibility of our medical devices are based on ISO 10993-1. The tests required for this standard are performed exclusively in accredited laboratories according to EN ISO/IEC 17025. Registered with the Union Data B for Safety in the Supply Chain (RAKCD) as known consignor since February 2012. All company processes are accompanied by occupational protection management.

IMPRESSION TAKING

A-SILICONE

| | |
|---|-------|
| Detaseal® hydroflow lite | p. 11 |
| Detaseal® hydroflow Xlite | p. 11 |
| Detaseal® hydroflow putty | p. 12 |
| Detaseal® hydroflow soft ² putty | p. 12 |
| Detaseal® hydroflow heavy | p. 13 |
| Detaseal® hydroflow mono | p. 13 |
| Detaseal® antilock | p. 14 |
| Detaseal® function | p. 15 |
| monoprint supra | p. 18 |
| onetime perfect | p. 18 |
| silaplast Futur A | p. 19 |
| silasoft® Futur A | p. 19 |

BITE REGISTRATION












| | |
|------------------|-------|
| Detaseal® bite | p. 15 |
| greenbite apple | p. 16 |
| greenbite colour | p. 16 |
| milkbite | p. 17 |
| glassbite | p. 17 |

C-SILICONE





























| | |
|-------------------|-------|
| silaplast Futur | p. 21 |
| Katalysatoren | p. 21 |
| silasoft® direct | p. 22 |
| silasoft® Special | p. 22 |
| silasoft® Normal | p. 23 |
| silone® | p. 23 |
| sta-seal | p. 23 |
| sta-seal f | p. 24 |
| hydro C | p. 25 |
| exaplast | p. 25 |























MATRIX A-SILICONES

| TRAY | | Viscosity | | Flowability | | Setting characteristics | Working time min. approx. | Time in the mouth min. approx. | Impression technique | combinable with | Features |
|-----------------------------------|---|------------------------|------|-------------|------|-------------------------|------------------------------|-----------------------------------|---|--|--|
| | | low | high | low | high | slow | | | | | |
| Detaseal® putty |  | very heavy bodied ▾ | | putty ▾ | | fast set ▾ | 1:30 | 2:00 | Correction & double mix technique | Detaseal® lite regular/fast set Detaseal® Xlite regular/fast set | Short time in the mouth |
| Detaseal® soft ² putty |  | very heavy bodied ▾ | | putty ▾ | | fast set ▾ | 1:15 | 1:30 | Correction & double mix technique | Detaseal® lite regular/fast set Detaseal® Xlite regular/fast set | Soft initial consistency, reduced final hardness, less time in the mouth |
| Detaseal® putty 5:1 |  | very heavy bodied ▾ | | putty ▾ | | regular set ▾ | 1:30 | 2:30 | Correction & double mix technique | Detaseal® lite regular/fast set Detaseal® Xlite regular/fast set | For dynamic mixing unit, hydrophilic |
| Detaseal® heavy 5:1 |  | heavy bodied ▾ | | stable ▾ | | fast set ▾ | 1:30 | 2:00 | Double mix technique | Detaseal® lite regular/fast set Detaseal® Xlite regular/fast set | For dynamic mixing unit, hydrophilic activity, thixotropic |
| Detaseal® mono |  | medium bodied ▾ | | stable ▾ | | fast set ▾ | 1:30 | 2:00 | Monophase & double mix technique | Detaseal® lite regular/fast set Detaseal® Xlite regular/fast set | Universal cartridge material, hydrophilic activity |
| Detaseal® mono 5:1 |  | medium bodied ▾ | | stable ▾ | | fast set ▾ | 1:30 | 2:00 | Monophase & double mix technique | Detaseal® lite regular/fast set, Detaseal® Xlite regular/fast set | Universal cartridge material, for dynamic mixing unit, hydrophilic activity |
| monoprint supra |  | medium bodied ▾ | | stable ▾ | | regular set ▾ | 2:00 | 2:00 | Monophase & double mix technique | Detaseal® lite regular/fast set Detaseal® Xlite regular/fast set | Transfer and impressions, implant- & combined technique, hydrophilic activity, thixotropic |
| onetime putty |  | very heavy bodied ▾ | | putty ▾ | | fast set ▾ | 1:30 | 2:00 | Double mix technique | onetime perfect wash | System for the one step technique, snap set, hydrophilic |
| onetime heavy 5:1 |  | heavy bodied ▾ | | stable ▾ | | fast set ▾ | 1:30 | 2:00 | Double mix technique | onetime perfect wash | For dynamic mixing unit, thixotropic |
| silaplast Futur A |  | very heavy bodied ▾ | | putty ▾ | | regular set ▾ | 1:30 | 3:00 | Correction & double mix technique | silasoft® Futur A | VPS variant of the classic, silaplast, soft and smooth initial consistency |
| FUNCTIONAL IMPRESSION | | | | | | | | | | | |
| Detaseal® function |  | heavy bodied ▾ | | putty ▾ | | slow set ▾ | 2:30 | 3:30 | Functional impressions for the complete prosthetics, functional border moulding | Detaseal® Xlite regular/fast set | Smooth & supple, putty, with defined plastic phase (2 min. 30 sec.) for functional movements |

CORRECTION / WASH

| | | Viscosity | | Flowability | | Setting characteristics | Working time min. approx. | Time in the mouth min. approx. | Impression technique | combinable with | Features |
|------------------------------------|---|---|------|---|------|---|------------------------------|-----------------------------------|--|---|---|
| | | low | high | low | high | slow | fast | | | | |
| Detaseal® lite regular set |  |  | |  | |  | 2:00 | 2:00 | Correction & double mix technique | Detaseal® putty, Detaseal® heavy | Hydrophilic activity, thixotropic, excellent wetting of humid areas |
| Detaseal® lite fast set |  |  | |  | |  | 1:00 | 1:30 | Correction & double mix technique | Detaseal® soft² putty, Detaseal® putty, Detaseal® heavy | Hydrophilic activity, thixotropic, excellent wetting of humid areas |
| Detaseal® Xlite regular set |  |  | |  | |  | 2:00 | 2:00 | Correction & double mix technique | Detaseal® putty, Detaseal® heavy | Hydrophilic activity, very good flow properties, excellent wetting of humid areas |
| Detaseal® Xlite fast set |  |  | |  | |  | 1:00 | 1:30 | Correction & double mix technique | Detaseal® soft² putty, Detaseal® putty, Detaseal® heavy | Hydrophilic activity, very good flow properties, excellent wetting of humid areas |
| onetime wash |  |  | |  | |  | 1:30 | 2:00 | Double mix technique | onetime perfect putty, onetime perfect Jumbo, monoprint supra | System for the one step technique, snap set |
| Detaseal® antilock |  |  | |  | |  | 1:00 | 1:20 | Double mix & correction technique, single phase impression | | Spacer Silicone supersoft, 10 Shore A, thin flowing, highly elastic, easy to demould, tear resistant, fast-curing |
| silasoft® Futur A |  |  | |  | |  | 1:00 | 1:20 | Correction & double mix technique | silaplast Futur A | VPS variant of the classic, excellent flow properties |

BITE REGISTRATION

| | | | | | | | | | | | |
|-------------------------|---|---|--|---|--|---|------|------|--------------------------------|--|--|
| Detaseal® bite |  |  | |  | |  | 0:45 | 1:30 | Bite key & intraoral fixations | | 90 Shore A |
| greenbite apple |  |  | |  | |  | 0:30 | 1:00 | Bite key & intraoral fixations | | 45 Shore D, apple flavoured |
| greenbite colour |  |  | |  | |  | 0:30 | 0:50 | Bite key & intraoral fixations | | 45 Shore D, with colour change |
| milkbite |  |  | |  | |  | 0:30 | 0:50 | Bite key & intraoral fixations | | 90 Shore A, light creamy consistency, scannable |
| glassbite |  |  | |  | |  | 0:30 | 1:00 | Bite key & intraoral fixations | | 80 Shore A, transparent, suitable as matrix material during light curing of composites |

Detaseal® hydroflow lite



Hydroflow precision impression silicone, easy-flow consistency, ultimate hydrophilicity, highest possible wetting capacity immediately after starting mixing (contact angle < 10°). Exact reproduction of the finest details, particularly in moist regions. „Perfect Flow“ – thixotropic setting with outstanding flow capacity, thermosensitive setting characteristics (SnapSet) for short oral dwelling time, maximum elastic properties after stretching, highly elastic when demoulding, extremely tear resistant. Medical device Cl. IIa, colour: Base pink, Catalyst white 1:1. In 2 setting types:

normal setting – regular set

fast setting – fast set



| Art.-No. | Packing | Unit |
|----------|---|----------------------|
| 02730 | Standard packing regular set base + catalyst mixing cannulas, yellow | 2 x 50 ml 12 pcs. |
| 02732 | Multi packing regular set base + catalyst mixing cannulas, yellow | 4 x 50 ml 24 pcs. |
| 02733 | Standard packing fast set base + catalyst mixing cannulas, yellow | 2 x 50 ml 12 pcs. |
| 02740 | Multi packing fast set base + catalyst mixing cannulas, yellow | 4 x 50 ml 24 pcs. |
| 02706 | mixing cannulas, yellow | 48 pcs. |
| 02712 | Intra-Oral-Tips, yellow | 96 pcs. |


Detaseal® hydroflow Xlite



Hydroflow precision impression silicone, extra light-flowing consistency. Precise detail reproduction, particularly in moist impression regions with ultimate hydrophilicity. Highest possible wetting capacity immediately after starting mixing (contact angle < 10°). „Perfect Flow“ – ideal flow properties, excellent SnapSet, maximum elastic properties – no stretching deformation when demoulding, extremely tear-resistant. Medical device Cl. IIa, colour: Base blue, Catalyst white 1:1. Thermoreactive setting characteristics:

normal setting – regular set

fast setting – fast set



| Art.-No. | Packing | Unit |
|----------|---|----------------------|
| 02741 | Standard packing regular set base + catalyst mixing cannulas, yellow | 2 x 50 ml 12 pcs. |
| 02743 | Multi packing regular set base + catalyst mixing cannulas, yellow | 4 x 50 ml 24 pcs. |
| 02744 | Standard packing fast set base + catalyst mixing cannulas, yellow | 2 x 50 ml 12 pcs. |
| 02745 | Multi packing fast set base + catalyst mixing cannulas, yellow | 4 x 50 ml 24 pcs. |
| 02706 | mixing cannulas, yellow | 48 pcs. |
| 02712 | Intra-Oral-Tips, yellow | 96 pcs. |

IMPRESSION TAKING A-SILICONES

Detaseal® hydroflow putty



Hydroflow impression silicone, kneadable consistency. Smooth initial consistency for perfect mixing results. Ideal for preliminary impressions in the impression technique fot corrections and as a tray material in the double-mix technique. Harmonized consistencies of tray and correction materials, no bubbling or flash. Easy to remove due to: high, elastic recovery from deformation with medium final hardness and maximum tensile strength. Time-saving thanks to short mixing time and time in the mouth (only 2 minutes!), neutral tasting, now also for 5:1 mixers; reaches its high, typically kneadable consistency immediately after filling the tray. Medical device Cl. IIa, colour: Base apple green, Catalyst white 1:1/jumbo 5:1



| Art.-No. | Packing | Unit |
|----------|---|-------------------------------|
| 02727 | Standard packing putty jar of base + catalyst measuring scoops | 2 x 250 ml 2 pcs. |
| 02728 | Multi packing putty jar of base + catalyst measuring scoops | 4 x 250 ml 2 pcs. |
| 03283 | Jumbo packing putty base + catalyst dynamic mixers, yellow rotating ring, yellow | 380 ml 10 pcs. 1 pc. |
| 03385 | Double Jumbo packing putty base + catalyst dynamic mixers, yellow rotating ring, yellow | 2 x 380 ml 2 pcs. 1 pc. |




Detaseal® hydroflow soft² putty



Hydroflow silicone impression material, kneadable consistency. Fluffy-soft initial consistency with simultaneously reduced final hardness (60 Shore A). For all precision impression techniques in combination with Detaseal® hydroflow lite/Xlite. Mixing & demoulding with unprecedented ease. Time savings due to a shorter curing time (2:45 min.) with a sufficient processing time of 1:15 min. Maximum patient comfort with an especially short intraoral curing time (1:30 min.!). Easy removal and pleasant mango aroma. Medical device Cl. IIa, colour: Base bright yellow, Catalyst white 1:1



| Art.-No. | Packing | Unit |
|----------|---|----------------------------|
| 02880 | Standard packing soft² putty jar of base jar of catalyst measuring scoops | 500 ml 500 ml 2 pcs. |



Detaseal® hydroflow heavy



Hydroflow impression silicone, heavy bodied consistency, hydrophilic, stable, easy to remove because of medium final hardness and high elastic recovery from deformation. Fast set, time in mouth 2 minutes. Ideal as a tray material for double-mix / sandwich technique in combination with Detaseal® hydroflow lite/Xlite. Harmonized consistencies, dimensional and volume stability, 5:1 hard cartridges for use with dynamic mixers. Medical device Cl. IIa, colour: Base apple green, Catalyst white 5:1



| Art.-No. | Packing | Unit |
|----------|---|-------------------------------|
| 02729 | Jumbo packing heavy base + catalyst dynamic mixers, yellow rotating ring, yellow | 380 ml 10 pcs. 1 pc. |
| 03388 | Double Jumbo packing heavy base + catalyst dynamic mixers, yellow rotating ring, yellow | 2 x 380 ml 2 pcs. 1 pc. |
| 02833 | Dynamic mixers, yellow | 40 pcs. |
| 02836 | Rotating ring, yellow | 1 pc. |

Detaseal® hydroflow mono



Hydroflow mono-phase impression silicone, medium bodied consistency for single step impression-taking or mixed technique, with distinct snap set. Constant precision, especially in humid areas, due to hydrophilic activation of the surface. Highly thixotropic, flows under pressure. Very tear-resistant with high, elastic recovery from deformation (99.8%), easy to remove. Fast set (total 3:30 min.), dimensional and volume stability, scannable. New: now also in 5:1 hard cartridges for use with dynamic mixers. Medical device Cl. IIa, colour: Base bright yellow, Catalyst white 1:1/jumbo 5:1




| Art.-No. | Packing | Unit |
|----------|--|-------------------------------|
| 02754 | Standard packing base + catalyst mixing cannulas, green | 2 x 50 ml 6 pcs. |
| 02757 | Multi packing base + catalyst mixing cannulas, yellow | 4 x 50 ml 12 pcs. |
| 03373 | Jumbo packing mono base + catalyst dynamic mixers, yellow rotating ring, yellow | 380 ml 10 pcs. 1 pc. |
| 03386 | Double Jumbo packing mono base + catalyst dynamic mixers, yellow rotating ring, yellow | 2 x 380 ml 2 pcs. 1 pc. |

Detaseal® antilock



VPS Spacer-Silicone, thin-flowing, soft consistency. Facilitates reliable demoulding and removal from the mouth, also in case of pronounced retention areas in the impression zone. For versatile use in undercut regions (wide open interdental spaces, pontic bridges or bar constructions, loose teeth with reduced stability, etc.). Short dwell time in the mouth, maximum recovery following expansion, highly elastic during demoulding, tear-resistant. Medical device Cl. IIa, colour: Base pink, catalyst white 1:1



| Art.-No. | Packing | Unit |
|----------|-------------------------|--------|
| 03535 | Standard packing | |
| | base + catalyst | 50 ml |
| | mixing cannulas, yellow | 6 pcs. |
| | Intra-Oral-Tips, yellow | 6 pcs. |

Spacer-Silicone
to block out undercuts

Detaseal® function



Precision impression material for functional impressions and functional border moulding. Addition curing, putty. Ideal smooth and supple consistency, easy to mould, well retaining its shape and volume. With defined plastic phase and prolonged setting time in the mouth. Optimum mixing check due to contrasting colours. Medical device Cl. IIa, colour: Base red, Catalyst yellow 1:1



| Art.-No. | Packing | Unit |
|----------|---|-------------------------|
| 02521 | Standard packing dispenser of base dispenser of catalyst moulding syringe | 80 ml 80 ml 1 pc. |
| 03098 | trayloc® A, bottle with brush | 17 ml |
| 03392 | trayloc® A, vaporizer of | 15 ml |
| 02802 | moulding syringes/application syringes | 10 pcs. |
| 03069 | easy application syringe, short | 1 pc. |

Detaseal® bite



Bite registration material, medium bodied consistency, stays on the impression or preparation. Direct application onto teeth, does not flow into interdental spaces. Easy to reseal due to neon colour, allows easy trimming (cut and fraise), high dimensional stability, fast setting, very high final hardness, approx. 90 Shore A. Medical device Cl. IIa, colour: Base yellowish green, Catalyst white 1:1



| Art.-No. | Packing | Unit |
|----------|--|---------------------------------|
| 02139 | Standard packing base + catalyst mixing cannulas, rose contouring tips | 2 x 50 ml 8 pcs. 8 pcs. |
| 02020 | Eco packing base + catalyst mixing cannulas, rose contouring tips | 8 x 50 ml 32 pcs. 32 pcs. |
| 02611 | contouring tips | 12 pcs. |

greenbite apple



Bite registration material, A-silicone based, medium bodied consistency, thixotropic. Direct application onto teeth, does not flow into interdental spaces, easily removed from mouth. Short time in the mouth (60 sec.), maximum final hardness, approx. 45 Shore D, easy processing. Medical device Cl. IIa, colour: Base green, Catalyst yellow 1:1. With green apple flavour!



| Art.-No. | Packing | Unit |
|----------|---|---------------------------------|
| 02641 | Standard packing base + catalyst mixing cannulas, green contouring tips | 2 x 50 ml 8 pcs. 8 pcs. |
| 02647 | Eco packing base + catalyst mixing cannulas, green contouring tips | 8 x 50 ml 32 pcs. 32 pcs. |
| 02705 | mixing cannulas, green | 48 pcs. |
| 02648 | contouring tips | 12 pcs. |

greenbite colour



Bite registration material, A-silicone based, with visual setting check. The colour change from green to yellow indicates that the setting process is completely finished. Medium bodied consistency, thixotropic. Direct application onto teeth. Short time in the mouth (50 sec.). Maximum final hardness, approx. 45 Shore D. Easy processing (cutting, grinding). Medical device Cl. IIa, colour: Base blue, Catalyst yellow 1:1. With apple flavour!



| Art.-No. | Packing | Unit |
|----------|---|---------------------------------|
| 02207 | Standard packing base + catalyst mixing cannulas, green contouring tips | 2 x 50 ml 8 pcs. 8 pcs. |
| 02726 | Eco packing base + catalyst mixing cannulas, green contouring tips | 8 x 50 ml 32 pcs. 32 pcs. |
| 02705 | mixing cannulas, green | 48 pcs. |
| 02648 | contouring tips | 12 pcs. |



milkbite



VPS bite registration, resistance-free bite registration with particularly light-creamy, airy initial consistency. Thixotropic formulation, stable, extra low extrusion pressure. Makes identification and fixation of the natural bite position without interference significantly easier. Time-saving direct application to the dentition, does not flow away into the interdental spaces, easily removed from mouth. Convenient processing time and extra short oral dwelling time (approx. 50 sec.), patient-friendly. Fracture-resistant, can be milled and cut, perfectly scannable. Final hardness: approx. 90 Shore A, medical device Cl. IIa, colour: creamy white 1:1, aroma: ice-cream



| Art.-No. | Packing | Unit |
|----------|--|-------------------------------|
| 03380 | Standard packing base + catalyst mixing cannulas, rose contouring tips | 2 x 50 ml 8 pcs. 8 pcs. |
| 03381 | Eco packing base + catalyst mixing cannulas, rose | 8 x 50 ml 32 pcs. |
| 02707 | mixing cannulas, rose | 48 pcs. |
| 02648 | contouring tips | 12 pcs. |



glassbite



Bite registration material, A-silicone based. Medium bodied consistency, thixotropic. Direct application onto teeth. Short time in the mouth (60 sec.). High final hardness, approx. 80 Shore A. Also suitable as matrix material during light curing of composites or during positioning of auxiliary parts, x-ray-transparent, easy working on (cutting, grinding). Medical device Cl. IIa, colour: transparent 1:1




| Art.-No. | Packing | Unit |
|----------|---|---------------------------------|
| 02571 | Standard packing base + catalyst mixing cannulas, green contouring tips | 50 ml 4 pcs. 4 pcs. |
| 02572 | Multipack base + catalyst mixing cannulas, green contouring tips | 4 x 50 ml 16 pcs. 16 pcs. |
| 02705 | mixing cannulas, green | 48 pcs. |
| 02648 | contouring tips | 12 pcs. |

IMPRESSION TAKING A-SILICONES

monoprint supra



VPS transfer silicone monophase, medium bodied consistency, with activated surface, specially designed for implantology and fixation impressions with combined prosthetics. Free-flowing, thixotropic setting for easy coverage of abutments. Secure fixation of transfer aids by fast hardening at the end of setting, thermosensitive setting behaviour. Easily removed from mouth, high elasticity. Final hardness approx. 70 Shore A, medical device Cl. IIa, colour: Base neonyellow, Catalyst white 1:1



| Art.-No. | Packing | Unit |
|----------|--|---------------------|
| 02525 | Standard packing base + catalyst mixing cannulas, green | 2 x 50 ml 6 pcs. |
| 02705 | mixing cannulas, green | 48 pcs. |

onetime perfect



Two-phase, one-step A-silicone impression system, for the one-step technique. Extremely short intraoral setting time (2 min.) due to thermally controlled setting. Simultaneous application of putty and heavy material together with wash in only one working stage. Medical device Cl. IIa.

Putty Very high consistency, kneadable, fast setting, 1:1

Heavy High consistency, fast set, 5:1

Wash Low consistency, light bodied. Super hydrophilic, thixotropic, very crisp detail, 1:1
Contrasting colours: Putty and Heavy – turquoise, Wash – violet



| Art.-No. | Packing | Unit |
|----------|--|---------------------------------|
| 03288 | Putty – Standard packing jar of base + catalyst measuring scoops | 2 x 250 ml 2 pcs. |
| 03289 | Wash – Standard packing base + catalyst mixing cannulas, rose Intra-Oral-Tips, white | 2 x 50 ml 12 pcs. 12 pcs. |
| 02155 | Putty – Eco packing jar of base + catalyst measuring scoops | 8 x 250 ml 4 pcs. |
| 02152 | Wash – Eco packing base + catalyst mixing cannulas, rose Intra-Oral-Tips, white | 8 x 50 ml 48 pcs. 48 pcs. |
| 03102 | Jumbo packing heavy base + catalyst dynamic mixers, yellow rotating ring, yellow | 380 ml 6 pcs. 1 pc. |
| 02833 | Dynamic mixers, yellow | 40 pcs. |

silaplast Futur A



Kneadable VPS impression material, normal setting, kneadable consistency. With the familiar soft smoothness of silaplast Classic and the ideal final hardness of 63 Shore A. Excellent transfer effect of the initial impression in correction technique or as tray material in double mix technique. Particularly easy and fast homogeneous mixing, high elasticity, easy to demould. Mixing ratio 1:1, contrast colours allow visual mixing control. Medical device Cl. IIa, colours: base white, catalyst blue, spearmint flavour



| Art.-No. | Packing | Unit |
|----------|------------------|--------|
| 02335 | Standard packing | |
| | jar of base | 450 ml |
| | jar of catalyst | 450 ml |
| | measuring scoops | 2 pcs. |

silasoft® Futur A














VPS precision impression material, easy flowing consistency. For correction technique and double mix technique, standard setting. Exact detail reproduction, highly elastic, extremely tear-resistant and with long-term dimensional stability. Thixotropic, perfect flow properties during positioning of impression. Highly elastic during demoulding, safe control, well readable. Easy to dose and mix, mixing ratio 1:1. Medical device Cl. IIa, colours: base pink, catalyst blue



| Art.-No. | Packing | Unit |
|----------|-------------------|-------|
| 02416 | Standard packing | |
| | base, tube of | 90 ml |
| | catalyst, tube of | 90 ml |

MATRIX C-SILICONES

| TRAY | | Viscosity | | Flowability | | Setting characteristics | | Working time min. approx. | Time in the mouth min. approx. | Impression technique | combinable with | Features |
|-------------------------|---|------------------------|------|-------------|------|-------------------------|------|------------------------------|-----------------------------------|---|--|--|
| | | low | high | low | high | slow | fast | | | | | |
| silaplast Futur |  | very heavy bodied ▼ | | putty ▼ | | regular set ▼ | | 1:15 | 2:30 | Correction & double mix technique | silasoft® N, silasoft® S, silasoft® direct | Catalyst paste or liquid |
| hydro C putty |  | very heavy bodied ▼ | | putty ▼ | | regular set ▼ | | 1:15 | 3:30 | Correction & double mix technique | hydro C wash | Smooth & supple, especially for the mix technique |
| exaplast putty |  | very heavy bodied ▼ | | putty ▼ | | regular set ▼ | | 1:15 | 3:30 | Correction & double mix technique | exasoft wash | Universally applicable for all impression techniques, soft smooth initial consistency |
| silone® |  | medium bodied ▼ | | medium ▼ | | regular set ▼ | | 1:15 | 2:30 | Monophase | silasoft® N, silasoft® S, silasoft® direct | Universal, for dentulous impressions, partial conditions |
| sta-seal |  | medium bodied ▼ | | medium ▼ | | regular set ▼ | | 1:15 | 2:30 | Monophase | | For edentulous impressions, flabby ridges, atrophic conditions |
| CORRECTION / WASH | | | | | | | | | | | | |
| silasoft® N |  | light bodied ▼ | | high ▼ | | regular set ▼ | | 1:15 | 2:30 | Correction & double mix technique | silaplast Futur | Excellent flow properties, flows under pressure into the smallest details, with catalyst paste or liquid |
| silasoft® S |  | light bodied ▼ | | high ▼ | | regular set ▼ | | 1:30 | 2:30 | Correction & double mix technique | silaplast Futur | For manual dosing |
| silasoft® direct |  | light bodied ▼ | | high ▼ | | regular set ▼ | | 1:15 | 3:00 | Correction & double mix technique | silaplast Futur | Automix2 system 10:1, manual dosing and mixing is omitted |
| hydro C wash |  | light bodied ▼ | | high ▼ | | regular set ▼ | | 1:15 | 2:30 | Correction & double mix technique | hydro C putty | Hydrophilic, high elasticity |
| exasoft wash |  | light bodied ▼ | | high ▼ | | regular set ▼ | | 1:15 | 2:30 | Correction & double mix technique | exaplast putty | Excellent „easy flow“ flow properties, hydrophilic |
| FUNCTIONAL IMPRESSION | | | | | | | | | | | | |
| sta-seal f |  | heavy bodied ▼ | | stable ▼ | | slow set ▼ | | 2:30 | 4:00 | Functional impressions for the complete prosthetics, functional border moulding | silasoft® N | With defined plastic phase (2 min. 30 sec.), special catalyst liquid, smooth & supple consistency |

silaplast Futur



Putty – kneadable consistency, smooth initial consistency – ideal final hardness and elastic strength. Multi-purpose use, e.g. preliminary impression material for corrections of preparations, overimpression of copperband impressions. New type of formula, less shrinkage than common C-silicones. Medical device Cl. IIa, colour: white



| Art.-No. | Packing | Unit |
|----------|---|--------------|
| 02001 | Standard packing jar | 900 ml |
| 02002 | 4-pack jars | 4 x 900 ml |
| 02003 | Clinic packing bucket | 1 x 5.400 ml |
| 02681 | Refill bag for jar | 900 ml |
| 02015 | cat p universal Futur, paste suitable for silaplast, silasoft® Normal, silasoft® Special, silone®, sta-seal | 5 x 35 ml |
| 02016 | silaplast Futur cat f, liquid suitable for silaplast | 5 x 50 ml |

Catalysts



Quality catalysts for C-silicones, with new hardener matrix, hypoallergenic, pleasant odour and taste. Blue coloured for exact measuring and mixing control. Available in packs of 5 pieces. Medical device Cl. IIa.



| Art.-No. | Packing | Unit |
|----------|---|-----------|
| 02015 | cat p universal Futur, paste suitable for silaplast, silasoft® Normal, silasoft® Special, silone®, sta-seal | 5 x 35 ml |
| 02016 | silaplast Futur cat f, liquid suitable for silaplast | 5 x 50 ml |
| 02017 | silasoft® cat f, liquid suitable for silasoft® Normal, silasoft® Special | 5 x 10 ml |
| 02018 | sta-seal cat f, liquid suitable for silone®, sta-seal, sta-seal f | 5 x 10 ml |

IMPRESSION TAKING C-SILICONES

silasoft® direct



Precision impression material for correction of preliminary impressions. Light bodied consistency, normal setting. Direct detail reproduction due to its excellent flow properties, hydrophilic. It flows under pressure into the smallest details. Optimal elastic recovery, exceptional tear strength and dimensional stability. Available with the Automix2 cartridge system 10:1, for direct application. Medical device Cl. IIa, colour: Base violet, Catalyst blue



| Art.-No. | Packing | Unit |
|----------|--|---------------------------------|
| 02809 | Standard packing base + catalyst 10:1 mixing cannulas, yellow 10:1 Intra-oral-tips, yellow | 2 x 50 ml 12 pcs. 12 pcs. |
| 02810 | Eco packing base + catalyst 10:1 mixing cannulas, yellow 10:1 Intra-oral-tips, yellow | 8 x 50 ml 48 pcs. 48 pcs. |
| 02807 | Mixing cannulas 10:1, yellow | 48 pcs. |
| 02712 | Intra-Oral-Tips, yellow | 96 pcs. |
| 02808 | Mixing gun for standard cartridges 10:1 | 1 pz. |

silasoft® Special



Precision impression material, with optimised formula, light bodied consistency, for correction of preliminary impressions, for the double impression and double mix technique. Optimum flow properties – the material flows freely when spatulated or syringed. It stays on the impression or preparation. High elasticity, precise detail reproduction and exceptional tear strength. Medical device Cl. IIa, colour: Base pink-transparent, Catalyst blue



| Art.-No. | Packing | Unit |
|----------|--|------------|
| 02275 | Standard packing tube | 160 ml |
| 02276 | 4-pack tubes | 4 x 160 ml |
| 02015 | cat p universal, paste suitable for silaplast, silasoft® Normal, silasoft® Special, silone®, sta-seal | 5 x 35 ml |
| 02017 | silasoft® cat f, liquid suitable for silasoft® Normal, silasoft® Special | 5 x 10 ml |

silasoft® Normal



Precision impression material, light bodied consistency, new type of formula, for correction of preliminary impressions, for the double-impression and double-mix technique. Exact detail reproduction, high elasticity, dimensional stability and tear-strength. Also available in a practical plastic bottle for long-term storage, continuous flow of material and exact dosage. Medical device Cl. IIa, colour: rose-opaque



| Art.-No. | Packing | Unit |
|----------|---|------------|
| 02004 | Standard packing tube | 160 ml |
| 02005 | 4-pack tubes | 4 x 160 ml |
| 02320 | Standard packing bottle | 160 ml |
| 02321 | 4-pack bottles | 4 x 160 ml |
| 02015 | cat p universal, paste suitable for silaplast, silasoft® Normal, silasoft® Special, silone®, sta-seal | 5 x 35 ml |
| 02017 | silasoft® cat f, liquid suitable for silasoft® Normal, silasoft® Special | 5 x 10 ml |

silone®



Universal precision impression material, medium bodied consistency, for all dentulous impressions. Especially recommended for partial-, copperband- and inlay impressions. Medical device Cl. IIa, colour: red-brown

sta-seal



Precision impression material, medium bodied consistency, for edentulous impressions, especially recommended for atrophic conditions, flabby ridges and relin impressions. Medical device Cl. IIa, colour: lime green



| Art.-No. | Packing | Unit |
|----------|---|------------|
| 02008 | silone® Standard packing tube | 160 ml |
| 02009 | silone® 3-pack tubes | 3 x 160 ml |
| 02010 | sta-seal Standard packing tube | 160 ml |
| 02011 | sta-seal 3-pack tubes | 3 x 160 ml |
| 02015 | cat p universal, paste suitable for silaplast, silasoft® Normal, silasoft® Special, silone®, sta-seal | 5 x 35 ml |
| 02018 | sta-seal cat f, liquid suitable for silone®, sta-seal, sta-seal f | 5 x 10 ml |

IMPRESSION TAKING C-SILICONES

sta-seal f



Special impression material, heavy bodied consistency, for functional impressions and functional border moulding. Excellent smooth initial consistency, easy to mould, ideal final hardness. The prolonged setting time in the mouth facilitates a convenient muscle trimming.

Medical device Cl. IIa, colour: Base green, Catalyst blue

| Art.-No. | Packing | Unit |
|----------|---|------------------------|
| 02012 | Standard packing measuring syringes sta-seal cat f bottle | 3 x 80 ml 1 x 10 ml |
| 02013 | Large packing measuring syringes sta-seal cat f bottles | 9 x 80 ml 3 x 10 ml |
| 03069 | easy application syringe, short | 1 pc. |
| 02802 | moulding syringes | 10 pcs. |
| 02018 | sta-seal cat f, liquid suitable for silone®, sta-seal, sta-seal f | 5 x 10 ml |



hydro C



Precision impression material based on C-silicone for two-step impression technique and two-step putty-wash technique, hydrophilic. Visual check of measuring and mixing with Colour-Safety system. Also suitable for anatomical and relined impressions, bite wafers. Pleasant apple aroma. Medical device Cl. IIa

Putty Kneadable consistency, dimensionally stable. Soft, smooth initial consistency, ideal final hardness.

Wash Light bodied consistency. Tearresistant, highly elastic, very crisp detail, hydrophilic.

Activator Hypoallergenic.

Contrast colours: Putty pastel-green, Wash rose-opaque, Activator blue



| Art.-No. | Packing | Unit |
|----------|-----------------------------------|----------|
| 02869 | Set | |
| | Putty, jar | 1.000 ml |
| | Wash, tube | 150 ml |
| | Activator, tube | 60 ml |
| 02242 | Standard packing Putty | 1.000 ml |
| 02471 | Eco packing Putty | 5.400 ml |
| 02243 | Standard packing Wash | 150 ml |
| 02245 | Standard packing Activator | 60 ml |

exaplast



High-precision premium impression material based on C-silicone, universally applicable for all impression techniques. Excellent readability and mixing check with clear contrast colours, patient-friendly „fresh-mint“ aroma. Medical device Cl. IIa

exaplast putty Kneadable consistency, soft smooth initial consistency, ideal final hardness, highly elastic, very dimensionally stable.

exasoft wash Excellent „easy flow“ flow properties, light bodied consistency, hydrophilic. High elastic strength, extremely tear-resistant, exact detail reproduction.

exagel cat Catalyst gel, hypoallergenic.

Colours: Putty yellow, Wash pink, Catalyst blue



| Art.-No. | Packing | Unit |
|----------|-------------------------------|--------|
| 03133 | Set | |
| | Putty, jar | 910 ml |
| | Wash, tube | 140 ml |
| | Cat, tube | 60 ml |
| 02907 | Standard packing Putty | 910 ml |
| 02910 | Standard packing Wash | 140 ml |
| 03125 | Standard packing Cat | 60 ml |

ACCESSORIES IMPRESSION TAKING

AUTOMIX

| | |
|---------------------------|-------|
| mixing gun Automix 1 & 2 | p. 27 |
| mixing tips Automix 1 & 2 | p. 27 |
| Contouring tips | p. 27 |
| Intra-Oral-Tips | p. 27 |

ACCESSORIES

| | |
|-----------------------------|-------|
| Anmischblock, Anmischspatel | S. 28 |
| Deta-Cut | S. 28 |
| Inlay-Löffel | S. 29 |

TRAY ADHESIVE

| | |
|--------------------------|-------|
| sili Adhesive | p. 29 |
| sili Solvent | p. 29 |
| trayloc® A | p. 30 |
| Detaseal® Adhesive rapid | p. 30 |





02699
Automix2 Mixing gun for standard
cartridges 1:1 / 2:1



02808
Automix2 Mixing gun for standard
cartridges 10:1



02343
Automix1 Mixing gun for standard
cartridges 1:1



Automix2

| Art.-No. | Packing | Unit |
|----------|--|---------|
| 02699 | Mixing gun for standard cartridges 1:1 / 2:1 | 1 pc. |
| 03617 | Mixing gun for standard cartridges „T“ | 1 pc. |
| 02808 | Mixing gun for standard cartridges 10:1 | 1 pc. |
| 02706 | Mixing cannulas, yellow | 48 pcs. |
| 03615 | T-Mixer, yellow | 25 pcs. |
| 02712 | Intra-Oral-Tips, yellow | 96 pcs. |
| 02707 | Mixing cannulas, rose | 48 pcs. |
| 03614 | T-Mixer, rose | 25 pcs. |
| 02345 | Intra-Oral-Tips, white | 96 pcs. |
| 02611 | Contouring tips | 12 pcs. |
| 02705 | Mixing cannulas, green | 48 pcs. |
| 03616 | T-Mixer, green | 25 pcs. |
| 02648 | Contouring tips | 12 pcs. |
| 02704 | Mixing cannulas, violet | 48 pcs. |
| 02711 | Intra-Oral-Tips, violet | 96 pcs. |
| 02708 | Mixing cannulas, blue, long | 48 pcs. |
| 02709 | Mixing cannulas, blue, short | 48 pcs. |
| 03617 | Mixing gun for standard cartridges „T“ | 1 pc. |
| 02808 | Mixing gun for standard cartridges 10:1 | 1 pc. |
| 02807 | Mixing cannulas 10:1, yellow | 48 pcs. |
| 02183 | Mixing cannulas 10:1, blue | 25 pcs. |
| 04055 | T-Mixer 10:1, blue | 25 pcs. |
| 02591 | Mixing cannulas, brown 4:1, mini-mix | 25 pcs. |



Automix1

| Art.-No. | Packing | Unit |
|----------|--|---------|
| 02343 | Mixing gun for standard cartridges 1:1 | 1 pc. |
| 02344 | Mixing cannulas, 6 mm | 48 pcs. |
| 02442 | Mixing cannulas, 9 mm | 48 pcs. |
| 02345 | Intra-Oral-Tips, white | 96 pcs. |
| 02611 | Contouring tips | 12 pcs. |



ACCESSORIES IMPRESSION TAKING

Mixing pad

Mixing pad vellum paper (P) or plastic coated (K), 30 sheets, with scaling, in two sizes, sheets separable. Medical device Cl. I



| Art.-No. | Packing | Unit |
|----------|---|----------------|
| 02031 | Mixing pad (P) large | 21,0 x 14,8 cm |
| 02032 | Mixing pad (P) small | 14,5 x 9,5 cm |
| 02033 | Mixing pad (K) large | 21,0 x 14,8 cm |
| 02034 | Mixing pad (K) small | 14,5 x 9,5 cm |
| 02036 | Mixing spatula, stainless steel, with wooden handle | 1 pc. |
| 02037 | Mixing spatula, blue | 1 pc. |
| 02830 | Mixing spatula, white | 1 pc. |

Deta-Cut

Sterilizable special instrument, curved, for time saving cutting out of thin and equal flow off grooves in preliminary silicone based impressions.



| Art.-No. | Packing | Unit |
|----------|----------|-------|
| 02323 | Deta-Cut | 1 pc. |

Inlay Trays



Inlay trays for making safer impressions, mechanically improved tray retention. Suitable for partial impressions, preliminary impressions with provisionals, inlays, crown- and bridge work. Disposable item, made of plastic material, perforated. Medical device Cl. I, colour: light blue



| Art.-No. | Packing | Unit |
|----------|-------------------|---------|
| 02042 | Trays, large | 12 pcs. |
| 02043 | Trays, small | 12 pcs. |
| 02044 | Trays, mixed each | 6 pcs. |

sili Adhesive



Adhesive for adhesion of condensation curing impression materials to plastic, metal and individual trays, extremely fast drying (approx. 1 min.). Medical device Cl. I, colour: blue-transparent



| Art.-No. | Packing | Unit |
|----------|--|--------|
| 02051 | sili Adhesive, spray can | 250 ml |
| 02053 | sili Adhesive Fluid, bottle with brush | 25 ml |

sili Solvent



Solvent for removal of the adhesive and cleaning of trays after impression taking.

| Art.-No. | Packing | Unit |
|----------|-------------------------|--------|
| 02052 | sili Solvent, spray can | 250 ml |

ACCESSORIES IMPRESSION TAKING

trayloc® A



Adhesive for A-silicones, very high adhesion with thinly coated surface conditioning, cementfree. Can be applied more easily and economically than conventional adhesives due to its aqueous-like, thin-flowing consistency. For all plastic, metal and individual trays. With colour indicator for visual application control, fast drying (approx. 3 min.). Also available as a spray, for simple, targeted handling. Medical device Cl. I, colour: blue-transparent



| Art.-No. | Packing | Unit |
|----------|-------------------|-------|
| 03098 | Bottle with brush | 17 ml |
| 03392 | Vaporizer of | 15 ml |

Detaseal® Adhesive rapid



For adhesion of addition curing impression materials to plastic, metal and individual trays, fast drying (approx. 1-2 min.). Medical device Cl. I, colour: orange-transparent



| Art.-No. | Packing | Unit |
|----------|-------------------|-------|
| 02324 | Bottle with brush | 15 ml |

for all trays
plastic, metal and
individual trays



SOFT RELINING

COLD CURING

| | |
|----------------|-------|
| mollosil® plus | p. 32 |
| mollosil® | p. 33 |

GLOSS VARNISH

| | |
|---------|-------|
| Lustrol | p. 33 |
|---------|-------|

ACCESSORIES

| | |
|--------------------|-------|
| mollosil® Remover | p. 32 |
| good morning Spray | p. 34 |

HEAT CURING

| | |
|---------------|-------|
| PERMALINER | p. 70 |
| Molloplast® B | p. 70 |

SOFT RELINING COLD CURING

mollosil® plus



Long-term soft relining material for chairside use, A-silicone based, cold curing, biocompatible. The material is directly applied from the cartridge onto the denture base. Strong bonding to the denture acrylic. Also suitable for application with the indirect technique. Primer with colour indicator for better application check. Booklets for patients and grinding set inclusive. Medical device Cl. IIa, colour: gingiva



Automix2



Automix1

| Art.-No. | Packing | Unit |
|----------|--|----------|
| 02353 | Standard packing cartridge Automix2 | 50 ml |
| | mixing cannulas, pink | 7 pcs. |
| | mollosil® plus Primer, bottle of | 5 ml |
| | mollosil® plus Polish, bottle of | 2 x 7 ml |
| | grinding sleeve, pointed | 1 pc. |
| | grinding sleeve holder, pointed | 1 pc. |
| | pipettes | 2 pcs. |
| 02354 | Refill packing cartridge Automix2 | 50 ml |
| 02436 | Standard packing cartridge Automix1 | 50 ml |
| | mixing cannulas, 6 mm | 7 pcs. |
| | mollosil® plus Primer, bottle of | 5 ml |
| | mollosil® plus Polish, bottle of | 2 x 7 ml |
| | grinding sleeve, pointed | 1 pc. |
| | grinding sleeve holder, pointed | 1 pc. |
| | pipettes | 2 pcs. |
| 02439 | Refill packing cartridge Automix1 | 50 ml |
| 02440 | mollosil® plus Primer | 5 ml |
| 02441 | mollosil® plus Polish | 2 x 7 ml |
| 02631 | mollosil® Remover, dropping bottle | 10 ml |
| 02615 | grinding sleeve, pointed | 10 pcs. |
| 02616 | grinding sleeve holder, pointed | 1 pc. |



Primer with
colour indicator

mollosil®



Long-term relining material for chairside use, A-silicone based, permanently soft, cold curing. Secure bonding to the denture acrylic. For manual mixing, mixing ratio 1:1. Also suitable for immediate dentures. Biocompatible, taste- and odourless. Also suitable for application with the indirect technique. Medical device Cl. IIa, colour: gingiva



| Art.-No. | Packing | Unit |
|----------|--|--|
| 03005 | Standard packing base, tube of catalyst, tube of adhesive, bottle of Lustrol gloss varnish, bottles of pipettes | 50 ml 50 ml 5 ml 2 x 6 ml 2 pcs. |
| 02274 | Starterkit base, tube of catalyst, tube of adhesive, bottle of Lustrol gloss varnish, bottles of pipettes | 30 ml 30 ml 5 ml 2 x 6 ml 2 pcs. |
| 03007 | mollosil® Adhesive | 5 ml |
| 03008 | Lustrol gloss varnish | 2 x 6 ml |
| 02631 | mollosil® Remover, dropping bottle | 10 ml |

Lustrol



Silicone based gloss varnish for sealing and smoothing of permanent soft mollosil® and Molloplast® B silicone based relinings, odour- and tasteless, transparent. Mixing ratio 1:1. Medical device Cl. IIa



| Art.-No. | Packing | Unit |
|----------|---|------------------------|
| 03008 | Standard packing gloss varnish, bottle of catalyst, bottle of pipettes | 6 ml 6 ml 2 pcs. |



for permanent use
biocompatible
also suitable for application
with the indirect technique

SOFT RELINING ACCESSORIES

good morning spray

Antibacterial spray solution for A-silicone based relinings, inhibits plaque and delays accumulation of oral microorganisms for 24 hours, without any adverse effects on the natural oral flora. Bacteriostatic effect (inhibits bacteria growth) due to the special combination of active substances, with chlorhexidine and ethereal oils.

| Art.-No. | Packing | Unit |
|----------|-------------------------------|------------|
| 02284 | Standard packing spray bottle | 100 ml |
| 02811 | 5-pack spray bottle | 5 x 100 ml |





TEMPORARY TREATMENT

C & B COMPOSITE

| | |
|------------------------|-------|
| tempofit® semi 10:1 | p. 36 |
| tempofit® premium 10:1 | p. 36 |
| tempofit® regular 1:1 | p. 37 |
| tempofit® duomix | p. 37 |
| tempofit® finish | p. 39 |

REPAIR

| | |
|---------------------|-------|
| smartrepair® system | p. 38 |
| smartrepair® rosa | p. 38 |

FIXATION / IMPLANT CEMENT

| | |
|---------------------------|-------|
| tempolink® single | p. 39 |
| tempolink® clear | p. 40 |
| implantlink® semi Classic | p. 40 |
| implantlink® semi Forte | p. 41 |
| implantlink® semi Xray | p. 41 |
| tempofill® 2 | p. 42 |

SEALING

| | |
|-----------|-------|
| Fermin | p. 42 |
| Detaferm® | p. 42 |

TEMPORARY TREATMENT C & B COMPOSITE

tempofit® semi 10:1



High-end bis-acrylic-composite for high quality long-term temporary prosthetics, 10:1, self-curing. For fast chairside fabrication of temporary crowns & bridges of highest mechanical strength and stability. Developed for restorations with a wear period > 30 days. With precisely defined elastic phase and minimal polymerization temperature. Esthetic, natural translucence, long-term colour-stable. A particularly thin inhibition layer and smooth surface accelerate finishing. Significant material savings due to use of shortened mixing tips. Neutral taste, direct application with the Automix2 system. Medical device Class IIa, colours: A2, A3, A3.5



| Art.-No. | Packing | Unit |
|----------|--|------------------|
| 04052 | Standard packing A2 cartridge (75 g) T-Mixer, blue, 10:1 | 50 ml 10 pcs. |
| 04053 | Standard packing A3 cartridge (75 g) T-Mixer, blue, 10:1 | 50 ml 10 pcs. |
| 04054 | Standard packing A3.5 cartridge (75 g) T-Mixer, blue, 10:1 | 50 ml 10 pcs. |
| 04055 | T-Mixer, blue, 10:1 | 25 pcs. |
| 02808 | Mixing gun for standard cartridges 10:1 | 1 pc. |

Long-term temporary restorations > 30 days
can be shade matched
using composites

tempofit® premium 10:1



High-end bis-acrylic composite 10:1, self curing, for high-quality temporary crowns and bridges and long-term temporaries. Natural dental aesthetics, brilliant translucent shades, even in variable light conditions. Extremely fracture-resistant, highest flexural and abrasion resistance. With a defined elastic phase, safe removal of temporaries, minimal curing temperature. Grinds and polishes extremely well. Neutral tasting and absolutely colour-stable throughout the wearing time. Direct application in the Automix2 system. Medical device Cl. IIa, colours: A1, A2, A3, A3.5, B1, BL



| Art.-No. | Packing | Unit |
|----------|--|------------------|
| 03248 | Standard packing A1 cartridge (75 g) mixing cannulas, blue, 10:1 | 50 ml 10 pcs. |
| 02201 | Standard packing A2 cartridge (75 g) mixing cannulas, blue, 10:1 | 50 ml 10 pcs. |
| 02209 | Standard packing A3 cartridge (75 g) mixing cannulas, blue, 10:1 | 50 ml 10 pcs. |
| 02217 | Standard packing A3.5 cartridge (75 g) mixing cannulas, blue, 10:1 | 50 ml 10 pcs. |
| 04315 | Standard packing B1 cartridge (75 g) Mischkanülen, blau 10:1 | 50 ml 10 St. |
| 04316 | Standard packing BL Kartusche (75 g) Mischkanülen, blau 10:1 | 50 ml 10 St. |

tempofit® regular 1:1



Standard bis-acrylic-composite for accurate short and long-term temporary prosthetics, 1:1, self-curing,. Natural translucent tooth aesthetics, with long-term colour stabiliser. Highly resistant to bending and breaking, very abrasion-resistant. High initial hardness, short curing time. Easy removal due to elastic phase, low polymerisation temperature. Optimal processing (grinding, polishing), neutral taste. Direct application with the Automix2-System. Medical device Cl. IIa, colours: A1, A2, A3, A3.5



| Art.-No. | Packing | Unit |
|----------|--|------------------|
| 03531 | Standard packing A1 cartridge (75 g) mixing cannulas, blue, short | 50 ml 10 pcs. |
| 03526 | Standard packing A2 cartridge (75 g) mixing cannulas, blue, short | 50 ml 10 pcs. |
| 03532 | Standard packing A3 cartridge (75 g) mixing cannulas, blue, short | 50 ml 10 pcs. |
| 03533 | Standard packing A3.5 cartridge (75 g) mixing cannulas, blue, short | 50 ml 10 pcs. |
| 02709 | mixing cannulas, blue, short | 48 pcs. |
| 02699 | Mixing gun for standard cartridges 1:1 / 2:1 | 1 pc. |

tempofit® duomix 5:1



Bis-acrylic composite for short- & long-term provisionals 5:1, self curing. Pressure and abrasion resistant, normal setting, highest occlusal fitting accuracy. Handmix in a synchron-dispenser, specially suitable for manual mixing, additional dispenser devices are not required. Exact dosage, without any loss of material in the mixing cannula, not dependent on type of mixing gun. Medical device Cl. IIa, colours: A2, A3



| Art.-No. | Packing | Unit |
|----------|--|-------------------------|
| 02763 | Introkit dispenser 2 x A2, 1 x A3 application syringe 1 mixing spatula, 1 mixing pad | each 25 g each 1 pc. |
| 02106 | Single-bag A2 application syringe 1 mixing spatula, 1 small mixing pad | 1 x 25 g each 1 pc. |
| 02270 | Single-bag A3 application syringe 1 mixing spatula, 1 small mixing pad | 1 x 25 g each 1 pc. |
| 02795 | Standard packing A2 | 4 x 25 g |
| 02762 | Standard packing A3 | 4 x 25 g |
| 02800 | Eco packing A2 | 8 x 25 g |
| 02872 | Eco packing A3 | 8 x 25 g |
| 02802 | moulding syringes/application syringes | 10 pcs. |

smartrepair® system



System-set of repair material for tempofit®, light curing, ideal consistency. For fast and precise repairs, delivered in economic direct application syringes with interchangeable application needles. Matched with the colours of tempofit®. With smartbond®, the one-step bonding for absolutely secure adhesion of the repair material. Medical device Cl. IIa, colours: A2, A3, A3.5



| Art.-No. | Packing | Unit |
|----------|--|---|
| 02694 | smartrepair® set application syringe A2, A3.5 smartbond®, bottle smartbrush application needles flowchart | each 1,5 g 5 ml 5 pcs. 6 pcs. 1 pc. |
| 02772 | smartrepair® A2 application needles | 1,5 g 5 pcs. |
| 02862 | smartrepair® A3 application needles | 1,5 g 5 pcs. |
| 02695 | smartrepair® A3.5 application needles | 1,5 g 5 pcs. |
| 02696 | smartbond®, bottle smartbrush | 5 ml 5 pcs. |
| 02715 | smartbrush | 100 pcs. |
| 02714 | application needles | 10 pcs. |

smartrepair® rosa

Modelling and repair resin, light-curing. For aesthetic adjustments and corrections of the gingival region of temporary crowns and bridge restorations or chairside repairs of acrylic prostheses. Universally applicable in practice and laboratory. Saves time with short polymerisation times. Fast, precise application directly from the syringe, economical to use. Ideal modelling properties, thixotropic, MMA-free. Can be polymerised with all standard dental polymerisation lights. Medical device Cl. IIa, colour: pink-transparent



| Art.-No. | Packing | Unit |
|----------|--|----------------|
| 03364 | smartrepair® rose application needles | 3 g 8 pcs. |
| 02696 | smartbond®, bottle smartbrush | 5 ml 5 pcs. |
| 02715 | smartbrush | 100 pcs. |
| 02714 | application needles | 10 pcs. |

smartrepair® System
also for individual shade matching of all
tempofit® C&B composites/FREEPRINT® temp p. 51

tempofit® finish



Light-curing single component varnish for surface sealing of temporary crowns and bridges. "Paint on" polish for chairside use by blue light polymerisation. Brilliant natural finish, without mechanical polishing. Creates a very hard and smooth surface. Also protects against plaque deposits and discolouration. High gloss, colourless, strong bonding. Medical device Cl. IIa



| Art.-No. | Packing | Unit |
|----------|---------------------------------|----------|
| 03536 | bottle | 4 ml |
| | disposable brush, blue metallic | 10 pcs. |
| 03542 | disposable brush, blue metallic | 100 pcs. |

tempolink® single



Temporary fixing cement, polycarboxylate based 1:1, for crowns and bridges, inlays and onlays, zirconium works, self-curing, eugenol-free, thixotropic. Optimum flowability due to its hydrophilic formulation. Secure sealing without any marginal leakage, maximum bonding. No dehydration of the preparation, easy to remove. Just 60 sec. for the fixation of single crowns and smaller applications. Medical device Cl. IIa, colour: natural opaque



| Art.-No. | Packing | Unit |
|----------|---------------------------------------|---------|
| 02598 | Standard packing | |
| | cartridge mini-mix | 5 ml |
| | mixing cannulas, brown, 1:1 | 8 pcs. |
| | flowchart | 1 pc. |
| 02761 | mixing cannulas, brown, 1:1, mini-mix | 25 pcs. |

tempolink® clear



Transparent, temporary fixing cement, resin based 4:1, dual curing. Film-forming, therefore highest marginal tightness and easy removal. Distinctly reduced displacement resistance than with common cements. Flexible working times due to optional light curing. Especially suitable for veneers, crowns and bridges in the anterior tooth region, without any aesthetic impairment. Eugenol-free, compatible with all permanent cements. Antibacterial with Triclosan, of indifferent colour and odour. Medical device Cl. IIa, colour: transparent



| Art.-No. | Packing | Unit |
|----------|---------------------------------------|---------|
| 02593 | Standard packing | |
| | cartridge mini-mix | 5 ml |
| | mixing cannulas, brown, 4:1 | 8 pcs. |
| | flowchart | 1 pc. |
| 02591 | mixing cannulas, brown, 4:1, mini-mix | 25 pcs. |


implantlink® semi Classic



Cement for fixing long-term temporary dentures on a resin base 4:1, for frictionlocked, chewing-stable cementing of superstructures. Undamaged removal by perfectly matched adhesive properties, significantly lower displacement resistance than conventional cements. Extremely low film thickness (8 µm) for maximum fitting accuracy, excess can be removed easily in large pieces during the plastic phase, does not adhere to soft tissue. Compatible with all material combinations, dual-curing, odour and taste-neutral. Medical device Cl. IIa, colour: natural opaque



| Art.-No. | Packing | Unit |
|----------|---------------------------------------|---------|
| 03092 | Standard packing | |
| | cartridge mini-mix | 5 ml |
| | mixing cannulas, brown, 4:1 | 10 pcs. |
| | flowchart | 1 pc. |
| 02591 | mixing cannulas, brown, 4:1, mini-mix | 25 pcs. |



implantlink® semi
Application-Clip!




implantlink® semi Forte



Implant cement for superstructures on a resin base 4:1, specially designed for semipermanent fixing on customised abutments, abutments with reduced adhesion surface or particularly small abutments. With reinforced adhesion and compressive strength. The restoration can still be removed without damage even after extended time in position. Accurate positioning with low displacement resistance and particularly low film thickness. Excess easily removed during the plastic phase. Compatible with all material combinations, dual curing system. Medical device Cl. IIa, colour: natural opaque




| Art.-No. | Packing | Unit |
|----------|--|--------------------------|
| 03371 | Standard packing cartridge mini-mix mixing cannulas, brown, 4:1 flowchart | 5 ml 10 pcs. 1 pc. |
| 02591 | mixing cannulas, brown, 4:1, mini-mix | 25 pcs. |



implantlink® semi Xray



Semi-permanent implant cement based on plastic 4:1, for securing suprastructures. Radiopaque with enhanced adhesion and compressive strength, even with extended wear. The restoration can be taken out without damaging it. Exact positioning and highest precision fit due to very thin film thickness. Easy, removal of large pieces of excess in the elastic phase, does not stick to soft tissue. Compatible with all material combinations, dual-hardening, neutral odour and taste. Medical device Cl. IIa, colour: natural opaque



| Art.-No. | Packing | Unit |
|----------|--|--------------------------|
| 02195 | Standard packing cartridge mini-mix mixing cannulas, brown, 4:1 flowchart | 5 ml 10 pcs. 1 pc. |
| 02591 | mixing cannulas, brown, 4:1, mini-mix | 25 pcs. |



TEMPORARY TREATMENT SEALING

tempofill® 2



One-part composite material, for temporary treatment of cavities and fillings (inlay/onlay technique and veneers), light curing. Direct, cost-effective dispensing from a syringe. Can be very easily condensed and contoured due to its soft, elastic consistency. Does not stick to the instrument. Cured tempofill® 2 fillings are easily removed without damaging the preparation margins. Medical device Cl. IIa, colour: light-opaque



| Art.-No. | Packing | Unit |
|----------|--|--------------------------|
| 02117 | Standard packing tempofill® 2 inlay syringes | 3 x 3 g |
| 02159 | Standard packing smartprotect® bottle of smartbrush flowchart | 7 ml 20 pcs. 1 pc. |
| 02715 | smartbrush | 100 pcs. |



Fermin



Ready-to-use filling material for provisional sealing. Setting in the mouth under the effect of saliva, great adhesive strength in the cavity, excellent marginal sealing, easy to remove. Medical device Cl. IIa, colour: white opaque

Detaferm®



Ready-to-use filling material for temporary cavity care, with calcium fluoride for the remineralisation of teeth. Setting in the mouth under the effect of saliva, great adhesive strength in the cavity. Simple processing, excellent marginal sealing, easy to remove. Medical device Cl. IIa, colour: white opaque



| Art.-No. | Packing | Unit |
|----------|---|------|
| 02054 | Fermin, Standard packing jar of 15 ml | 40 g |
| 02000 | Detaferm®, Standard packing jar of 15 ml | 40 g |



PREVENTION/ REGENERATION

SEALING OF ROOT DENTINE

smartseal® cervi p. 44

FISSURE SEALING

smartseal® & loc p. 45

smartseal® & loc F p. 45

MOUTH RINSING SOLUTION

smartclean p. 46

smartseal® cervi



Sealing varnish to protect exposed, hypersensitive root dentine. Increases mechanical stability, forms an effective erosion barrier and reduces the cervical sensitivity of root dentine. High abrasion resistance due to high degree of linking of highly multifunctional (meth)acrylates. Forms a permanent protective layer against mechanical abrasion and the erosive effects of acids. Structure-reinforcing, desensitising, light-curing. Precise, hygienic application with flexible, fibre-free, plastic micro-applicators. Medical device Cl. IIa



| Art.-No. | Packing | Unit |
|----------|---|----------------------------|
| 03537 | Standard packing dropper bottle applicators, pink disposable dual dishes | 5 ml 50 pcs. 24 pcs. |
| 03543 | Applicators, pink | 100 pcs. |



forms an effective
erosion barrier
reduces cervical
sensitivity



smartseal® & loc / F



Fissure sealant with wet-bonding technology, composite-based for permanent sealing of fissures and enamel defects. Moisture tolerant, hydrophilic. Application on moist or slightly dry enamel surfaces, very effective wetting of tooth enamel. Increased strong bonding due to additional enamel integration of reactive coated fillers. Low viscosity and easy flow for margin-free, tight transitions without cracking or splintering. Absolutely precise direct application with extra long disposable needles. Light-curing, abrasion-proof. Medical device Cl. IIa

smartseal® & loc

Colour: transparent, without fluoride, especially for aesthetic sealing

smartseal® & loc F

Colour: natural opaque, with fluoride



Wet surface welcome!
with or without fluoride

| Art.-No. | Packing | Unit |
|----------|--|---------------------------|
| 02612 | smartseal® & loc Combi Pack smartseal® & loc, syringe smartseal® & loc F, syringe application needles, black | 1,5 g 1,5 g 12 pcs. |
| 02573 | Standard packing smartseal® & loc application syringes application needles, black | 2 x 1,5 g 12 pcs. |
| 02581 | Standard packing smartseal® & loc F application syringes application needles, black | 2 x 1,5 g 12 pcs. |
| 02456 | application needles, black | 25 pcs. |

smartclean

Mouth rinsing solution for use after professional cleaning and pretreatment, delays formation of plaque. Contains chlorhexidine. Easy, safe use in 5 ml syringes with flexing rinsing cannulas. Non-alcoholic, sugar-free. Refill maxi syringe 50 ml, with integrated adapter for effortless, clean refilling of the 5 ml syringes.

| Art.-No. | Packing | Unit |
|----------|---|--------------------|
| 02554 | Standard packing syringes rinsing cannulas | 4 x 5 ml 4 pcs. |
| 02557 | Refill-Container maxi-syringe with integrated adapter | 50 ml |
| 02629 | rinsing cannulas | 25 pcs. |



LIQUIDS, SOLUTIONS

SPRAY SOLUTION

good morning spray

p. 51





LIQUIDS, SOLUTIONS

good morning spray

Antibacterial spray solution for A-silicone based relinings, inhibits plaque and delays accumulation of oral microorganisms for 24 hours, without any adverse effects on the natural oral flora. Bacteriostatic effect (inhibits bacteria growth) due to the special combination of active substances, with chlorhexidine and ethereal oils.

| Art.-No. | Packing | Unit |
|----------|--|------------|
| 02284 | good moring spray, spray bottle | 100 ml |
| 02811 | good moring spray, 5-pack spray bottle | 5 x 100 ml |



3D PRINTING

CROWNS & BRIDGES

FREEPRINT® TEMP p. 51

FREEPRINT® CROWN p. 52

DENTURE BASES

FREEPRINT® DENTURE p. 53

FREEPRINT® DENTURE IMPACT p. 54

FREEPRINT® TRYIN p. 54

ORTHODONTIC BASE COMPONENTS

FREEPRINT® ORTHO p. 55

SPLINTS

FREEPRINT® SPLINT 2.0 p. 55

FREEPRINT® SPLINTMASTER p. 56

BRACKET TRANSFER TEMPLATES

FREEPRINT® IBT p. 56

TRAYS & PLATES

FREEPRINT® TRAY 2.0 p. 57

MODEL PRODUCTION

FREEPRINT® MODEL p. 57

FREEPRINT® MODEL 2.0 p. 58

FREEPRINT® MODEL PRO p. 58

FREEPRINT® MODEL T p. 59

FREEPRINT® MODEL WW p. 59

GINGIVAL MASKS

FREEPRINT® GINGIVA p. 60

CASTINGS

FREEPRINT® CAST 2.0 p. 60



MATRIX 3D RESINS

| Material type | Application | Color | Characteristics | Medical Devices Class MDR | Medical Devices Class FDA | Medical Devices Class NMPA |
|--|---|--------------------------------|--|---------------------------|---------------------------|----------------------------|
| TEMP  | <ul style="list-style-type: none"> Temporary crowns & bridges Temporary anterior and posterior tooth restorations | A1, A2, A3 | <ul style="list-style-type: none"> Natural transparency and tooth esthetics Extremely high construction precision High mechanical stability Biocompatible | IIa | II | – |
| CROWN  | <ul style="list-style-type: none"> Permanent crowns, denture teeth Long-term temporary bridges | A1, A2, A3, B1, B3, C2, D3, BL | <ul style="list-style-type: none"> Natural transparency and tooth esthetics Highest abrasion resistance | IIa | II | – |
| DENTURE  | <ul style="list-style-type: none"> Removable denture bases Total prosthesis | Pink-transparent, pink | <ul style="list-style-type: none"> Long-term stable and biocompatible dentures Fast printing Perfect fit | IIa | II | – |
| DENTURE IMPACT  | <ul style="list-style-type: none"> Removable denture bases Total prosthesis | Pink-transparent, pink | <ul style="list-style-type: none"> Long-term stable and biocompatible dentures High impact resistance Easy post-processing | in process | in process | – |
| TRYIN  | <ul style="list-style-type: none"> Individual functional try-ins | A2 | <ul style="list-style-type: none"> Fast, material-saving production High mechanical stability | IIa | I | – |
| ORTHO  | <ul style="list-style-type: none"> Surgical guides for implant dentistry Orthodontic base components | Clear-transparent | <ul style="list-style-type: none"> Very high mechanical stability & construction precision High printing speed Sterilizable Biocompatible | IIa | I | TEC resin |
| SPLINT 2.0  | <ul style="list-style-type: none"> Hard splints | Clear-transparent | <ul style="list-style-type: none"> High mechanical flexural strength and stability High initial final hardness Biocompatible | IIa | I | TEC resin |
| SPLINTMASTER  | <ul style="list-style-type: none"> Flexible splints Repositioners Mouthguards Nightguards | Clear-transparent | <ul style="list-style-type: none"> Flexible High tensile strength High tension-free wearing comfort | IIa | II | – |
| IBT  | <ul style="list-style-type: none"> Orthodontic transfer trays for positioning brackets | Transparent | <ul style="list-style-type: none"> Elastic and tear-resistant Reliable fixing of brackets Biocompatible | I | I | – |
| TRAY 2.0  | <ul style="list-style-type: none"> Individual impression trays Functional impression trays Base resin plates | Green | <ul style="list-style-type: none"> High dimensional stability, torsional rigidity Max. construction speed Compatible with all impression materials Biocompatible | I | I | MED resin |

| Material type | Application | Color | Characteristics | Medical Devices Class MDR | Medical Devices Class FDA | Medical Devices Class NMPA |
|--|---|---------------------------------|--|---------------------------|---------------------------|----------------------------|
| MODEL  | <ul style="list-style-type: none"> Dental working and situation models Control models | Ivory, grey, sand | <ul style="list-style-type: none"> Maximum surface hardness Dimensional stability Comfortable haptic Very good production precision | TEC resin | TEC resin | I |
| MODEL 2.0  | <ul style="list-style-type: none"> Dental master and working models Control models | Caramel, grey, light grey, sand | <ul style="list-style-type: none"> High detail reproduction Maximum surface hardness and dimensional stability Plaster-like appearance and haptic Very good construction precision | TEC resin | TEC resin | I |
| MODEL PRO  | <ul style="list-style-type: none"> Dental working and situation models Control models | Caramel, grey, sand | <ul style="list-style-type: none"> Plaster-like haptic Good dimensional stability High flexural strengths Easy post-processing | TEC resin | TEC resin | in process |
| MODEL KFO  | <ul style="list-style-type: none"> Laboratory model preparation Orthodontic models | White | <ul style="list-style-type: none"> Plaster-like haptic Distinctive edge and dimension stability Highest surface quality | TEC resin | TEC resin | – |
| MODEL T  | <ul style="list-style-type: none"> Working models for thermoforming technique and aligner technology | Light blue | <ul style="list-style-type: none"> High temperature resistance to process-related temperature stress High edge strength | TEC resin | TEC resin | I |
| MODEL WW  | <ul style="list-style-type: none"> Working models for thermoforming technique and aligner technology | Blue-transparent | <ul style="list-style-type: none"> Water-washable High temperature resistance | TEC resin | TEC resin | – |
| GINGIVA  | <ul style="list-style-type: none"> Flexible gingival masks for dental models | Gingiva | <ul style="list-style-type: none"> 3D reproduction of functional gingival model segments Excellent elasticity and tear resistance Natural gingiva esthetics | TEC resin | TEC resin | I |
| CAST 2.0  | <ul style="list-style-type: none"> Dental casting objects for precision casting | Red-transparent | <ul style="list-style-type: none"> Residue-free burning out High dimensional stability after printing Precise and distortion-free results, even for delicate constructions | TEC resin | TEC resin | I |

MDR Medical Device Regulation EU FDA Food and Drug Administration USA NMPA National Medical Products Administration China

FREEPRINT® TEMP



Light-curing, biocompatible resin for the 3D printing of temporary crowns & bridges. Exceptional breaking and flexural resistance. The natural-looking translucent shades (according to VITA classical A1–D4 shade guide) can be esthetically modified for single-crown and bridge restorations. High level of oral stability and in conjunction with tempolink®, enable excellent marginal seal during the period of wear. Easy polishing results in very high surface quality with exceptional abrasion resistance. No odour or taste irritations, MMA- & THF-MA-free. Medical device Class IIa, colours: A1, A2, A3

| Art.-No. | Packing | Unit |
|----------|--------------------|---------|
| 04058 | FREEPRINT® TEMP A1 | 500 g |
| 04059 | FREEPRINT® TEMP A2 | 500 g |
| 04060 | FREEPRINT® TEMP A3 | 500 g |
| 04062 | FREEPRINT® TEMP A1 | 1.000 g |
| 04063 | FREEPRINT® TEMP A2 | 1.000 g |
| 04064 | FREEPRINT® TEMP A3 | 1.000 g |



FREEPRINT® TEMP
brilliant, translucent colours!

can be shade matched using the
smartrepair® system and composites p. 38



FREEPRINT®
Application-Clip!



3D PRINTING

FREEPRINT® CROWN



Light-curing formulation for 3D printing of permanent crowns, denture teeth and long-term temporary bridges. Brilliant colors (according to VITA classical A1-D4 shade guide) thanks to perfectly matched transparency and opacity. Optimum dimensional stability of the restorations due to maximum bending and abrasion resistance. No tendency to discolour thanks to low water absorption. Fast, uncomplicated cleaning process. Easy grinding and polishing due to minimal surface chalking. MMA- & THF-MA-free. Medical device Cl. IIa, colour: A1, A2, A3, B1, B3, C2, D3, BL

| Art.-No. | Packing | Unit |
|----------|---------------------|---------|
| 02372 | FREEPRINT® CROWN A1 | 500 g |
| 02378 | FREEPRINT® CROWN A2 | 500 g |
| 02417 | FREEPRINT® CROWN A3 | 500 g |
| 02481 | FREEPRINT® CROWN B1 | 500 g |
| 02645 | FREEPRINT® CROWN B3 | 500 g |
| 02766 | FREEPRINT® CROWN C2 | 500 g |
| 02783 | FREEPRINT® CROWN D3 | 500 g |
| 02845 | FREEPRINT® CROWN BL | 500 g |
| 02376 | FREEPRINT® CROWN A1 | 1.000 g |
| 02415 | FREEPRINT® CROWN A2 | 1.000 g |
| 02446 | FREEPRINT® CROWN A3 | 1.000 g |
| 02519 | FREEPRINT® CROWN B1 | 1.000 g |
| 02758 | FREEPRINT® CROWN B3 | 1.000 g |
| 02782 | FREEPRINT® CROWN C2 | 1.000 g |
| 02825 | FREEPRINT® CROWN D3 | 1.000 g |
| 02884 | FREEPRINT® CROWN BL | 1.000 g |



NEW:
FREEPRINT® CROWN
SCIENTIFIC FACTS





FREEPRINT® DENTURE

Light-curing, biocompatible resin for the 3D printing of removable denture bases. High wearing comfort. Process reliability thanks to high initial hardness following printing. Low viscosity for reduced material consumption and faster cleaning. Very high surface quality, excellent to polish. Displays maximum mechanical flexural and tensile strength without becoming brittle. The rigid denture base withstands high loads in the oral environment. Extremely low shrinkage values compared to PMMA materials. Fully compatible with prefabricated, 3D printed (FREEPRINT® CROWN) or milled teeth. Validated with the VITA VIONIC VIGO System. Natural aesthetics and a light transparency enable alignment with natural gingival color. Odorless and tasteless, MMA- and THF-MA-free, Medical device Cl. IIa, colour: pink-transparent, pink

| Art.-No. | Packing | Unit |
|----------|-------------------------------------|--------|
| 02060 | FREEPRINT® DENTURE pink-transparent | 500 g |
| 02040 | FREEPRINT® DENTURE pink-transparent | 1000 g |
| 04092 | FREEPRINT® DENTURE pink | 1000 g |



NEW:
FREEPRINT®
3D CATALOGUE



FREEPRINT® DENTURE IMPACT

Light-curing formulation for 3D printing of impact resistant denture bases. The increased impact strength provides extremely high fracture stability and ensures lasting functionality over the entire wearing time. Colours and transparency optimally adapted to classic denture resins. cation of fit function and occlusion. Easy to process. Compatible with fabricated, printed (FREEPRINT®CROWN) or milled teeth. MMA-, THF-MA- and TPO-free, tasteless. Medical device Cl. IIa, colors: pink-transparent, pink

| Art.-No. | Packing | Unit |
|----------|---|--------|
| 04436 | FREEPRINT® DENTURE IMPACT pink-transparent* | 1000 g |
| 04437 | FREEPRINT® DENTURE IMPACT pink* | 1000 g |

*approx. available 2024



FREEPRINT® TRYIN



Light-curing formulation for the 3D printing of individual functional try-ins of digitally manufactured restorations for complete and partial dentures in esthetically pleasing tooth color. Easy verification of fit function and occlusion. Fast and easy additive manufacturing of functional try-ins of individual tooth setups. Easy to process. MMA- & THF-MA-free. Medical device Cl. IIa, colour: A2

| Art.-No. | Packing | Unit |
|----------|----------------------|---------|
| 04427 | FREEPRINT® TRYIN A2* | 1.000 g |



3D PRINTING

FREEPRINT® ORTHO



Light-curing, biocompatible resin for the 3D printing of base parts for orthodontic apparatuses, drilling and X-ray templates, occlusal and fixation splints. Benefits: maximum reliability of construction process, mechanical strength, easy to polish. Breaking strength, elasticity and influence of moisture following the orthodontic standard. Validated processes for autoclave sterilization according to EN ISO 17664! Adaptation of the occlusion field, easy fixation of orthodontic wire elements and auxiliary components with FREEFORM® fixgel. Resistant to oral conditions, MMA-free, neutral odour and taste. Medical device Cl. IIa, colour: clear-transparent

| Art.-No. | Packing | Unit |
|----------|----------------------|---------|
| 03989 | FREEPRINT® ORTHO | 1.000 g |
| 03988 | FREEPRINT® ORTHO 405 | 1.000 g |



FREEPRINT® SPLINT 2.0



Light-curing, biocompatible resin for the 3D printing of splints, fixation and transfer keys. Hard occlusal splint, rigid type, with high efficiency. Highest mechanical bending and breaking strength, without being brittle. Clear-transparent formulation with high accuracy of fit for visual control in the working area. Process reliability through high initial hardness and final strength. Low-viscosity approach for low material consumption and more rapid cleaning. Compatible with FREEFORM® plast for additional occlusal design in practice. MMA- and THF-MA-free, neutral odour and taste. Medical device Cl. IIa, colour: clear-transparent

| Art.-No. | Packing | Unit |
|----------|-----------------------|---------|
| 02080 | FREEPRINT® SPLINT 2.0 | 500 g |
| 02076 | FREEPRINT® SPLINT 2.0 | 1.000 g |





FREEPRINT® SPLINTMASTER

Light-curing formulation for 3D printing of flexible high-end splints. Available in two levels of flexibility: Taff & Flex. The biocompatible material is characterized by extreme impact strength. Splints made of FREEPRINT® SPLINTMASTER are long-term stable and robust, without becoming brittle. Thanks to its high elastic recovery, the splints, always return to their original shape. The perfectly adjusted properties guarantee highest, tension free wearing comfort and easy insertion and removal. The flexible splints are odour- and tasteless and free from MMA, THF-MA and TPO. Medical device Class IIa. Colour: clear-transparent


| Art.-No. | Packing | Unit |
|----------|-------------------------------|---------|
| 04433 | FREEPRINT® SPLINTMASTER TAFF* | 1.000 g |
| 04432 | FREEPRINT® SPLINTMASTER FLEX* | 1.000 g |



FREEPRINT® IBT

Light-curing, biocompatible resin for the 3D printing of orthodontic transfer trays for positioning brackets. Easy, precise positioning and application of the brackets due to the indirect bonding technique. The transparent bracket transfer templates allow easy visual control. High tensile strength and flexibility allow the easy positioning and subsequent removal of the templates in one single work step. Bisphenol A-, MMA- and THF-MA-free, odorless and tasteless. Medical device Cl. I, colour: transparent

| Art.-No. | Packing | Unit |
|----------|----------------|---------|
| 04248 | FREEPRINT® IBT | 500 g |
| 04249 | FREEPRINT® IBT | 1.000 g |





3D PRINTING

FREEPRINT® TRAY 2.0



Light-curing, biocompatible resin for the 3D printing of individual impression and functional trays and base plates. Low viscosity for low material consumption and more rapid cleaning. Highest dimensional stability and torsional rigidity for accurate and distortion-free impression taking. Perfect for implant impression taking within the digital workflow. Printable with 200 µm layer thickness. Compatible for all tray adhesives and impression materials. MMA- and THF-MA-free. Medical device Cl. I, colour: green

| Art.-No. | Packing | Unit |
|----------|---------------------|---------|
| 02505 | FREEPRINT® TRAY 2.0 | 1.000 g |



FREEPRINT® MODEL

Light-curing formulation for 3D printing of dental master and working models. Perfect detail reproduction due to plaster-like colours: grey, ivory, sand. Fast-printing, maximum surface hardness and shape stability. Maximum precision in construction, feel and stability meet the high demands of model production. Opaque colour adjustment for the optical detection of surface structure, preparation margins, etc. Bisphenol A- & MMA-free. Colours: ivory, grey, sand

| Art.-No. | Packing | Unit |
|----------|----------------------------|---------|
| 03778 | FREEPRINT® MODEL sand | 1.000 g |
| 03780 | FREEPRINT® MODEL ivory | 1.000 g |
| 03782 | FREEPRINT® MODEL grey | 1.000 g |
| 03065 | FREEPRINT® MODEL 405 sand | 1.000 g |
| 03779 | FREEPRINT® MODEL 405 ivory | 1.000 g |
| 03781 | FREEPRINT® MODEL 405 grey | 1.000 g |



FREEPRINT® MODEL 2.0

Light-curing formulation for 3D printing of dental models, master models, situation and orthodontic models. High detail reproduction, maximum surface hardness and shape stability provide an excellent mechanical strain. Faster workflow due to shortened post-processing. Maximum precision in construction, feel and stability meet the high demands of model production. The distinct edge stability and abrasion resistance make the models comparable to conventional plaster models in terms of handling. Opaque colour adjustment for the optical detection of surface structure, preparation margins, etc. MMA- & THF-MA-free. Colours: caramel, grey, light grey, sand, white


| Art.-No. | Packing | Unit |
|----------|--------------------------------|---------|
| 02128 | FREEPRINT® MODEL 2.0 sand | 1.000 g |
| 02099 | FREEPRINT® MODEL 2.0 ight-grey | 1.000 g |
| 02177 | FREEPRINT® MODEL 2.0 grey | 1.000 g |
| 02850 | FREEPRINT® MODEL 2.0 caramel | 1.000 g |
| 02148 | FREEPRINT® MODEL 2.0 white | 1.000 g |



FREEPRINT® MODEL PRO

Light-curing formulation for 3D printing of dental master and working models. Fast-printing. Precise reproduction of details, maximum surface hardness and shape stability. Maximum precision in construction, feel and stability meet the high demands of model production. The further increase in mechanical strength contributes to optimum functionality. The formulation has been revised according to the latest REACH Regulation. Colours: caramel, grey, sand

| Art.-Nr. | Packung | Einheit |
|----------|------------------------------|---------|
| 04440 | FREEPRINT® MODEL PRO caramel | 1.000 g |
| 04438 | FREEPRINT® MODEL PRO grey | 1.000 g |
| 04439 | FREEPRINT® MODEL PRO sand | 1.000 g |




3D PRINTING

FREEPRINT® MODEL T

Light-curing, temperature-resistant resin for the 3D printing of dental models for the thermoforming technique. Precise reproduction of details, maximum surface hardness and edge stability of the models. Precision in construction, tactility and stability meet the high demands of model fabrication. The pronounced intrinsic stability also allows the fabrication of hollow thermoformed models. Low viscosity for considerably easier cleaning and low material consumption. High temperature resistance. The stability of the models is preserved even during heating in thermoforming. MMA-free. Colour: light blue


| Art.-No. | Packing | Unit |
|----------|--------------------|---------|
| 02332 | FREEPRINT® MODEL T | 1.000 g |



FREEPRINT® MODEL WW

Light-curing formulation for the 3D printing of dental models for the thermoforming technique. The water-washable material meets all requirements for digitally produced models in thermoforming. No use of solvent necessary. The high edge strength and good intrinsic stability of the material allow production of hollow thermoformed models. High temperature resistance. The stability of the models is not affected by the heat. Cost-efficient model production. MMA- & THF-MA-free. Colour: blue-transparent

| Art.-No. | Packing | Unit |
|----------|---------------------|---------|
| 03105 | FREEPRINT® MODEL WW | 1.000 g |



FREEPRINT® GINGIVA

Light-curing for the 3D printing of flexible gingival masks for dental models. For 3D reproduction of functional gingival model segments in a digital workflow in conjunction with FREEPRINT® MODEL. Excellent elasticity and tear-resistance. Excellent dimensional stability, no shrinkage or aging, ductile even when stored over a long period. Natural gingiva esthetics, no annoying or unpleasant odors from the completed gingival masks. Bisphenol A-, MMA- and THF-MA-free. Colour: gingiva


| Art.-No. | Packing | Unit |
|----------|--------------------|---------|
| 02820 | FREEPRINT® GINGIVA | 500 g |
| 02843 | FREEPRINT® GINGIVA | 1.000 g |



FREEPRINT® CAST 2.0

Light-curing resin for the 3D printing of high-precision casting objects. Excellent dimensional stability after printing ensures safe application with fit check, in particular where cast models and molds are concerned. Printing results are precise and without any distortions, even for extremely delicate and intricate models. Any corrections or repairs after printing are possible with easyform gel LC. Low viscosity and optimized depth of cure help reduce printing time and material consumption and facilitate faster cleaning. Burns without leaving any residues in the mold and is the basis for high-grade modeling results that ensure the perfect fit. Phosphate-bonded investment materials are used as investments. MMA & THF-MA free. Colour: red

| Art.-No. | Packing | Unit |
|----------|---------------------|---------|
| 02548 | FREEPRINT® CAST 2.0 | 500 g |
| 02632 | FREEPRINT® CAST 2.0 | 1.000 g |



LABORATORY

LIGHT CURING POLYMERS

FREEFORM® p. 62

easyform gel / paste p. 63

INSULATION

FLEXISTONE® p. 64

FLEXISTONE® Plus p. 64

PUTTY

blue eco p. 65

blue eco stone p. 65

blue eco scan® p. 66

compact lab putty p. 66

GINGIVAL MASK

Esthetic Mask p. 67

Implant Mask / scan p. 68

GINGIVAMOLL® p. 71

RELINING

PERMALINER p. 69

Molloplast® B p. 69

Primo p. 70

Tools p. 70

Molloplast® Remover p. 70

RESOURCES

Itopol / Oropol p. 71

xynon p. 72

glad Spray p. 72

modell Spray p. 72



FREEFORM®




Ready-to-use modelling resin, light-curing for multiple indications (e.g., drilling and X-ray templates, splints, base plates, transfer keys). Radiotransparent, can be sterilised and disinfected, biocompatible. Medical device Cl. IIa, colour: clear-transparent plast kneadable paste, can be shaped and modelled, dimension-stable, MMA-free. fixgel stable gel for direct application from the dosing syringe, MMA-free. High mechanical stability for rapid fixation and integration of drill guides or X-ray reference spheres. Also ideal for reinforcing, individualising and repairing. bond light-curing bonding agent for PMMA-compatible thermoforming films and acrylates. Also for repairs & revisions of FREEFORM® plast. coat light-curing, single-component finishing gloss for the final surface sealing. Creates a smooth, hard surface. Brilliant finish, totally without polishing.



| Art.-N. | Packing | Unit |
|---------|---|----------|
| 02204 | Starterset | |
| | plast, jar | 60 g |
| | fixgel, syringe | 3 g |
| | bond, bottle | 5 ml |
| | coat, bottle | 10 ml |
| | disposable brush, blue metallic | 10 pcs. |
| | application needles | 10 pcs. |
| 02542 | Standard packing plast, jar | 60 g |
| 02829 | Standard packing fixgel, syringe | 3 g |
| 02702 | Standard packing bond, bottle | 5 ml |
| 02549 | Standard packing coat, bottle | 10 ml |
| 03542 | Disposable brush, blue metallic | 100 pcs. |
| 02714 | Application needles | 10 pcs. |



FREEFORM®
Application-Clip
YouTube



easyform gel / paste




Ready-for-use modelling resin, light-curing, as gel or paste, MMA-free. Very high fitting accuracy for multiple indications, particularly for casting. Burns off without residue, even with thicker layer. Matching viscosities of gel and paste for fast, accurate work. Dimensionally stable, no polymerisation delay, high final hardness and stability. Can be combined with conventional or light-curing wax. Visual layer thickness control with slight translucency. Colour: red

Gel Direct application from the syringe, stable, easily guided for modelling.

Paste Dough-like consistency, easily shaped and modelled.



| Art.-No. | Packing | Unit |
|----------|---|--------------------|
| 03376 | Standard packing gel syringes application needles | 3 x 3 g 10 pcs. |
| 02714 | application needles | 10 pcs. |
| 03375 | Standard packing paste jar of | 30 g |



Burns off without residue, even
with thicker layer!



FLEXISTONE®

Insulating and modelling material, C-silicone based, light bodied, excellent stability, very high final hardness (75 Shore A). For insulating acrylic or porcelain teeth during investing. Can be easily applied in one single procedure onto acrylic teeth and wax, forming a thin, very hard insulating layer for fixing of teeth during pressing. No retention crystals are required. Colour: green



| Art.-No. | Packing | Unit |
|----------|--|-----------------------------------|
| 03009 | Standard packing bottle with stand of (350 g) catalyst, bottle of mixing bowl | 190 ml 30 ml 1 pc. |
| 03010 | Large packing bottles with stand of (5 x 350 g) catalyst, bottle of mixing bowls | 5 x 190 ml 3 x 30 ml 2 pcs. |
| 03165 | Catalyst, bottle of | 30 ml |

FLEXISTONE® Plus

Insulating and modelling material, A-silicone based, light bodied, excellent stability, very high final hardness (85 Shore A). Optimum flow properties during spatulation. For insulating dentures during investing. Protection of acrylics and acrylic teeth against plaster. Easy application - 1:1 mixing system. Especially indicated for use with injection moulding techniques, recommended by IVOCAR for the SR-Ivocap-System. Colour: Base blue, Catalyst white 1:1



| Art.-No. | Packing | Unit |
|----------|--|------------------------|
| 02383 | Standard packing base + catalyst bottles with stands mixing bowl | 2 x 160 ml 1 pc. |
| 02329 | Eco packing base + catalyst plastic bottles mixing bowls | 2 x 1.000 ml 4 pcs. |

LABORATORY PUTTY

blue eco

Putty for the dental laboratory, addition curing, very high, easily kneadable consistency. Precise detail rendition, stable in shape and volume, smooth, heat-resistant up to 200 °C, standard final hardness 86 Shore A. Excellent mechanical properties, easy to cut. For versatile use (secure fixation of artificial teeth, precasts, checkbites, models, etc.). High economy, optimum mixing control by contrast colours. Colour: Base blue, Catalyst beige 1:1



| Art.-No. | Packing | Unit |
|----------|---|--------------------------|
| 02038 | Standard packing base, jar of catalyst, jar of measuring scoops | 800 g 800 g 2 pcs. |
| 02467 | Mega pack base, bucket of catalyst, bucket of measuring scoops | 5 kg 5 kg 2 pcs. |

blue eco stone


Laboratory putty, A-silicone based, high viscosity, smooth initial consistency. High definition, dimensionally stable, especially smooth, heat-resistance up to 200 °C. Extra high final hardness, approx. 96 Shore A, for versatile use (fabrication of precasts, checkbites, models for denture repairs, blockings-out, etc.), very economical, easy to cut. Colour: Base dark blue, Catalyst white 1:1



| Art.-No. | Packing | Unit |
|----------|---|------------------------------|
| 03117 | Standard packing base, jar of catalyst, jar of measuring scoops | 800 g 800 g 2 pcs. |
| 02517 | Eco pack base, jar of catalyst, jar of measuring scoops | 1.400 g 1.400 g 2 pcs. |
| 02985 | Mega pack base, bucket of catalyst, bucket of measuring scoops | 5 kg 5 kg 2 pcs. |

blue eco scan®

Scannable laboratory putty, addition curing, high, smooth consistency, non sticky. For powderless digital 3D data acquisition of the counter bite. Can be read by all CAD/CAM/CIM systems, optimum exposure quality and artifact-free image reproduction at very high contrast. For all conventional indications of laboratory putty. Low displacement resistance, precise detail rendition, stable in shape and volume, with firm edges, easy to cut. No breaking out of thin contact point areas. Final hardness: 85 Shore A. Colour: Base light grey, Catalyst white 1:1



| Art.-No. | Packing | Unit |
|----------|------------------|--------|
| 03118 | Standard packing | |
| | base, jar of | 800 g |
| | catalyst, jar of | 800 g |
| | measuring scoops | 2 pcs. |

compact lab putty

Putty for laboratories, C-silicone based, smooth initial consistency, high final hardness (90 Shore A), thermoresistant. Multi-purpose use (for fixation of teeth, insulation, key to occlusal recording, blocking out, etc.). Variable working and setting time due to individual dosage of catalyst. Exact detail reproduction, high dimensional stability. Colour: mouse grey





| Art.-No. | Packing | Unit |
|----------|---------------------------------|-----------|
| 02410 | Standard packing | |
| | plastic bucket of (8 kg) | 5.400 ml |
| | measuring scoops | 2 pcs. |
| 02296 | cat compact, liquid, bottle of | 50 ml |
| 02465 | Eco-Set | |
| | plastic bucket of (8 kg) | 5.400 ml |
| | cat compact, liquid, bottles of | 3 x 50 ml |
| | measuring scoops | 2 pcs. |
| 02797 | Eco-Set | |
| | plastic bucket of (8 kg) | 5.400 ml |
| | cat compact, paste, tubes of | 5 x 25 ml |
| | measuring scoops | 2 pcs. |
| 02823 | cat compact, paste, tube of | 25 ml |

Esthetic Mask

Flexible gingival mask, A-silicone based, for the fabrication of master models and working casts in the dental technique, cold curing, high dimensional stability and tear strength. Fast and easy processing, approx. 45 Shore A, mixing ratio of 1:1, with re-usable syringes. For perfect natural aesthetics. Available in tubes and Automix2-cartridges. Colour: gingiva






| Art.-No. | Packing | Unit |
|----------|---|-----------|
| 02340 | Standard packing | |
| | base, tube of | 50 ml |
| | catalyst, tube of | 50 ml |
| | separating liquid for silicones, vaporizer of | 15 ml |
| | disposable syringes | 2 pcs. |
| 03113 | Cartridge packing | |
| | cartridges | 2 x 50 ml |
| | mixing cannulas, yellow | 10 pcs. |
| | separating liquid for silicones, vaporizer of | 15 ml |
| 02706 | mixing cannulas, yellow | 48 pcs. |
| 02712 | Intra-Oral-Tips, yellow | 96 pcs. |
| 02690 | Separating Agent | |
| | separating liquid for silicones, vaporizer of | 15 ml |



Implant Mask / scan

Addition curing special silicone for manufacturing flexible gingival masks for implant procedures. Free-flowing, dimensionally stable, tear-resistant. Flows evenly around impression coping with free-flowing properties. Very high elasticity allows repositioning as required. Very easily milled with high final hardness, approx. 70 Shore A.

Scan Version, for unpowdered digital 3D recording. Readable by all standard digital scanners, high recording quality for interference-free image reproduction and data processing. Colour: gingiva



| Art.-No. | Packing | Unit |
|----------|---|-------------------------------|
| 02529 | Standard packing cartridges mini-mix 1:1 mixing cannulas, light green, 1:1 separating liquid for silicones, vaporizer of | 2 x 10 ml 6 pcs. 15 ml |
| 03374 | Standard packing scan cartridges mini-mix 1:1 mixing cannulas, light green, 1:1 separating liquid for silicones, vaporizer of | 2 x 10 ml 6 pcs. 15 ml |
| 02605 | mixing cannulas, light green 1:1, mini-mix | 25 pcs. |
| 03552 | Cartridge packing scan cartridges mixing cannulas, blue, short separating liquid for silicones, vaporizer of | 2 x 50 ml 10 pcs. 15 ml |
| 02709 | mixing cannulas, blue, short | 48 pcs. |
| 02690 | Separating Agent separating liquid for silicones, vaporizer of | 15 ml |

PERMALINER



Ready-to-use hot cure silicone for the permanent lining of dentures. Permanently soft with reduced polymerisation time. Particularly suitable for the flask technique in the dental laboratory, ideal for the manufacturing of new prostheses. Expanded indication spectrum, even with the most difficult jaw conditions. Holds absolutely securely to all acrylates, without adhesive. Reduces the deposit of oral-specific bacteria through an anti-plaque effect. No softening, no embrittlement, remains elastic even after long periods of wear. Biocompatible, MMA-free, neutral odour and taste. Medical device Cl. IIa, colour: light pink



| Art.-No. | Packing | Unit |
|----------|-------------------------------------|------|
| 02092 | Standard packing jar of | 30 g |
| 02146 | Large packing jar of | 70 g |
| 03153 | Pre-polishers (10 pcs. + 1 mandres) | set |
| 03120 | Cutters, 2 pcs. | set |
| 03101 | Cutters, 3 pcs. | set |

50 % shorter
polymerisation time



Molloplast® B
Application-Clip
YouTube

Molloplast® B



Permanent soft relining material for processing in dental laboratories. One-component silicone, heat curing, microwave-polymerising, wide indications (e.g. for obturators, mouth guards, etc.). Improved bonding to new and old denture base acrylics. Optimised consistency due to new vacuum filling technique. Resisting the influences of the oral environment, with anti-plaque effect. It's outstanding properties are clinically tested and documented. Medical device Cl. IIa



| Art.-No. | Packing | Unit |
|----------|--|---------------|
| 03000 | Standard packing jar of | 45 g |
| 03001 | Large packing jar of | 170 g |
| 03002 | Combi packing jar of Primo Adhesive, bottle of | 45 g 15 ml |
| 03004 | Primo Adhesive bottle with brush of | 15 ml |



Primo



Particularly for acrylic dentures, increases and ensures the bonding of Molloplast® B to the denture base. Also especially suitable in case of repeated repairs, extensions and relinings of acrylates. Application by means of the integrated brush. Medical device Cl. IIa



| Art.-No. | Packing | Unit |
|----------|--|-------|
| 03004 | Primo Adhesive bottle with brush of | 15 ml |

Molloplast® Remover

Dropping solution for fast and secure dissolution of the bonding between elastic relining materials and the denture acrylic. Simple application, without any timeconsuming removing by grinding. No residues, no interaction during repeated relinings.

Bearbeitungswerkzeuge

Grinding sleeves

For smoothing and removing of uneven spots and for final finishing of permanent soft relining materials and silicones.

Cutters & pre-polishers

For optimal trimming of permanent soft, silicone based relinings. No smear effect during cutting, negligent heat development thus no separation of acrylic and relining material. Final finishing with pre-polishers (smoothing and removal of uneven spots).

| Art.-No. | Packing | Unit |
|----------|--------------------------------------|---------|
| 02206 | Molloplast® Remover, dropping bottle | 10 ml |
| 03101 | Cutters, 3 pcs. | set |
| 03120 | Cutters, 2 pcs. | set |
| 03153 | Pre-polishers (10 pcs. + 1 mandres) | set |
| 02616 | Grinding sleeve holder, pointed | 1 pc. |
| 02615 | Grinding sleeves, pointed | 10 pcs. |



LABORATORY

GINGIVAMOLL®

Silicone based epithetic material, heat-curing, for the fabrication of flexible gingival epithesis as a temporary replacement for the receded gingiva – after periodontal treatment. For the correction of aesthetic and phonetic insufficiencies and after implantations, clinically tested. Also available as Starterkit. Medical device Cl. IIa



| Art.-No. | Packing | Unit |
|----------|--------------------------------|------------|
| 03053 | Starterkit | |
| | base material, pink/light pink | each 2,5 g |
| | base material, opaque | 4 g |
| | protective lacquer, bottle of | 5 ml |
| | paint liquid, bottle of | 5 ml |
| | paint brush | 1 pc. |
| | paint colours, blue/red | each 1 pc. |

Itopol / Oropol

High gloss polishing pastes, non-greasy, water-soluble and non-smearing.
ITOPOL for cobalt-chromium alloys
OROPOL for precious metals, non-precious metals and acrylic resins



| Art.-No. | Packing | Unit |
|----------|--------------------|-------|
| 04032 | ITOPOL, green, bar | 650 g |
| 04033 | OROPOL, white, bar | 650 g |

xynon

Liquid insulating agent. Insulates plaster, resin, metal and porcelain from wax. For application with brushes or for directed spray coating with a cannula. Does not form an oily film on the plaster surface.

glad spray


Wetting agent for the casting technique, for smooth and bubble-free castings. Easy handling by spraying, thus no damage to the cast object.



| Art.-No. | Packing | Unit |
|----------|------------------------------|----------|
| 02077 | xynon Spray can with cannula | 250 ml |
| 02078 | xynon Plastic bottle | 100 ml |
| 02079 | xynon Plastic bottle | 1.000 ml |
| 02074 | glad spray Spray can | 250 ml |

modell spray

Special lacquer for the model cast technique. For smoothing, hardening and sealing of model surfaces. Easy application by spraying instead of dipping.



| Art.-No. | Packing | Unit |
|----------|------------------------|--------|
| 02075 | modell spray Spray can | 250 ml |

INDEX

B

| | |
|----------------------|----|
| blue eco | 65 |
| blue eco scan® | 66 |
| blue eco stone | 65 |

C

| | |
|-------------------------|----|
| catalysts | 21 |
| compact lab putty | 66 |

D

| | |
|---|----|
| Deta-Cut | 28 |
| Detaferm® | 42 |
| Detaseal® Adhesive rapid | 30 |
| Detaseal® antilock | 14 |
| Detaseal® bite | 15 |
| Detaseal® function | 15 |
| Detaseal® hydroflow heavy | 13 |
| Detaseal® hydroflow lite | 11 |
| Detaseal® hydroflow mono | 13 |
| Detaseal® hydroflow putty | 12 |
| Detaseal® hydroflow soft ² putty | 12 |
| Detaseal® hydroflow Xlite | 11 |

E

| | |
|----------------------------|----|
| easyform gel / paste | 63 |
| Esthetic Mask | 67 |
| exaplast | 25 |

F

| | |
|---------------------------------|----|
| Fermin | 42 |
| FLEXISTONE® | 64 |
| FLEXISTONE® Plus | 64 |
| FREEFORM® | 62 |
| FREEPRINT® CAST 2.0 | 60 |
| FREEPRINT® DENTURE | 53 |
| FREEPRINT® DENTURE IMPACT | 54 |
| FREEPRINT® GINGIVA | 60 |
| FREEPRINT® IBT | 56 |
| FREEPRINT® MODEL | 57 |
| FREEPRINT® MODEL 2.0 | 58 |
| FREEPRINT® MODEL PRO | 58 |
| FREEPRINT® MODEL T | 59 |
| FREEPRINT® MODEL WW | 59 |

| | |
|-------------------------------|----|
| FREEPRINT® ORTHO | 55 |
| FREEPRINT® SPLINT 2.0 | 55 |
| FREEPRINT® SPLINTMASTER | 56 |
| FREEPRINT® TEMP | 51 |
| FREEPRINT® TRAY 2.0 | 57 |
| FREEPRINT® TRYIN | 54 |

G

| | |
|--------------------------|----|
| GINGIVAMOLL® | 71 |
| glad Spray | 72 |
| glassbite | 17 |
| good morning spray | 34 |
| greenbite apple | 16 |
| greenbite colour | 16 |

H

| | |
|---------------|----|
| hydro C | 25 |
|---------------|----|

I

| | |
|---------------------------------|----|
| Implant Mask / scan | 68 |
| implantlink® semi Classic | 40 |
| implantlink® semi Forte | 41 |
| implantlink® semi Xray | 41 |
| Inlay Trays | 29 |
| Intra-Oral-Tips | 27 |
| Itopol / Oropol | 71 |

L

| | |
|---------------|----|
| Lustrol | 33 |
|---------------|----|

M

| | |
|----------------------------------|----|
| milkbite | 17 |
| mixing gun Automix 1 & 2 | 27 |
| mixing pad | 28 |
| mixing tips Automix 1 & 2 | 27 |
| modell Spray | 72 |
| Molloplast® B | 69 |
| Molloplast® Remover | 70 |
| Molloplast® Rotating Tools | 70 |
| mollosil® | 33 |
| mollosil® plus | 32 |
| mollosil® Remover | 32 |
| monoprint supra | 18 |

O

| | |
|-----------------------|----|
| onetime perfect | 18 |
|-----------------------|----|

P

| | |
|------------------|----|
| PERMALINER | 69 |
| Primo | 70 |

S

| | |
|---------------------------|----|
| silaplast Futur | 21 |
| silaplast Futur A | 19 |
| silasoft® Futur A | 19 |
| silasoft® direct | 22 |
| silasoft® Normal | 23 |
| silasoft® Special | 22 |
| sili Haftlack | 29 |
| sili Haftlöser | 29 |
| silone® | 23 |
| smartclean | 46 |
| smartrepair® rosa | 38 |
| smartrepair® system | 38 |
| smartseal® & loc | 45 |
| smartseal® & loc F | 45 |
| smartseal® cervi | 44 |
| sta-seal | 23 |
| sta-seal f | 24 |

T

| | |
|------------------------------|----|
| tempofill® 2 | 42 |
| tempofit® duomix 5:1 | 37 |
| tempofit® finish | 39 |
| tempofit® premium 10:1 | 36 |
| tempofit® regular 1:1 | 37 |
| tempofit® semi 10:1 | 36 |
| tempolink® clear | 40 |
| tempolink® single | 39 |
| trayloc® A | 30 |

X

| | |
|-------------|----|
| xynon | 72 |
|-------------|----|

HIGH IMPRESSION TAKING
SOFT RELINING
PREVENTION & REGENERATION
PERFORMANCE PRODUCT VALUE
INNOVATIVE MATERIALS
UNLIMITED POSSIBILITIES
POLYMERS LABORATORY
TEMPORARY TREATMENT
LIQUIDS & SOLUTIONS

DETAX GmbH · Carl-Zeiss-Straße 4 · 76275 Ettlingen/Germany
Tel. +49 7243 / 510-0 · Fax +49 7243 / 510-100 · post@detax.com · www.detax.com

3D RESINS