

# *simply the best*

## *cost/performance ratio*



## green eco

### Premium Ear Impression Material

Putty-soft consistency; 1:1 mixing ratio. Multipurpose impression material, premium quality, unbeatable cost/performance ratio.

**3-year product guarantee!**

**Features:** High elastic recovery due to an innovating polymer, exceptional tear strength, lasting dimensional stability.  
Ideal final hardness: 35 Shore A.

- Highest detail reproduction, thixotropic (does not drip)
- Green and yellow contrasting colours enable optimised mixing
- Non-sticky with "snap set effect" (fast setting at body temperature)
- Quick and easy processing (trimming), excellent scanning results

