

bioplast



Addition curing silicone for **indirect fitting** (in the lab) of soft-flexible earmoulds, eg. BTE & noise protection. If required, the earmould can be finished with a special silicone lacquer for permanently surface sealing.

Features:



- multi-purpose
- highest wearing comfort due to soft fitting
- robust and long-term durable
- extremely tear resistant
- guaranteed shelf life: 3 years
- final hardness: high elastic 25, 40 or 60 Shore A
- available in clear-transparent, rose-transparent & white
- easy processing (grinding, cutting)
- skin-friendly & biocompatible