

ultra


for perfect impressions



multi purpose,
non-sticky, thixotropic
37 Shore A

supra

the CIC expert material



low viscosity,
particularly for CIC
and narrow canals,
35 Shore A

addition

FirstClass putty

Premium quality
super soft
highest precision



basic A

All-Inclusive silicone



Best cost/performance ratio

medical **silicones**. Made in Germany by
DETAX premium ear impression materials
biocompatible nano fillers. Creative ideas
platinum catalyst high performance chemistry www.detax.de

