

...on your marks, get set...white



That was the countdown on the Detax booth at the EUHA in Frankfurt as their latest product innovations were being presented: addition chroma® the first ear impression with colour change – from blue to white in three minutes. The colour change indicates that the setting process is completely finished and the impression can be removed. This visual control of setting makes work easier, and is thus the logical further development of conventional ear impression compounds.

Also presented was bioplast – the allrounder in white. Addition curing silicone for indirect and direct fitting of soft-flexible earmoulds, e.g. BTE and noise protection. Now available in white, 40 or 60 shore A. Bioplast earmoulds are robust, long-term durable, easy processing, skin-friendly,

biocompatible, and extremely tear resistant. And yet another silicone with colour effect: aquaplus® magic, for direct fitting of long-term noise and water plugs – changes colour dependent on temperature.



Detax would like to thank all visitors to their booth for the lively dialogue and look forward to seeing them all again in 2007.

For further information visit: www.detax.de