

No more – long drying times or inconvenient skin formation.

## The new DETAX silicone lacquers!

Faster drying, better adhesion | lasting brilliance.



### **express coat**

*wiping-proof after 8 min.*

*fast-drying one component lacquer,  
air-drying, for smoothing & sealing,  
ideal flow properties, brilliant finish*

### **bio coat**

*solvent-free*

*special lacquer for silicone-  
otoplastics, odourless,  
biocompatible, heat curing*

### **dip coat**

*no annoying skin formation*

*special dip coat for silicone-earmoulds  
& ear impressions,  
high-gloss, low viscosity,  
strong bonding, cold curing*

### **supercoat nano**

*with nano-silver*

*antibacterial sealing lacquer for  
silicone based earmoulds &  
otoplastics, air-drying, lasting  
protection against bacteria*