



NEW: addition mini

„the first ear impression material for kids“

Precision impression series,
perfectly adjusted to the
needs of little patients and the
requirements for pediatric adaptations

- Maximum safety: does not contain conventional fillers, extremely tear resistant, highest elastic recovery from deformation (99.8%)
- light bodied for entirely pressure-free impression taking
- UNIQUE: Handy minimix cartridges, hassle-free use (even without gun)
- Child friendly: extra short setting time, (just 2 min.!), fun colour, yummy smell
- Skin-friendly, biocompatible



Junior

Children & teenagers, light bodied,
neon green with scent of gummi bear

Sensitive

Babies, infants, patients with allergies, extra
light bodied, no dye or flavouring agents

exclusively distributed in the UK by:



info@puretone.net
www.puretone.net