

# **addition<sup>®</sup> chroma**

the first impression material with colour change



The first addition curing silicone impression material with colour change from blue to white in 3 minutes. The colour change from blue to white indicates that the setting process is completely finished and the impression can be removed. This visual control of setting saves a lot of work and is thus the logical further development of conventional ear impression compounds.

## **Features:**

- maximum detail reproduction, low viscosity, particularly for CIC
- just 3 min. in the ear
- visual setting check due to
- colour change from blue to white
- thixotropic, does not drip
- ideal final hardness: 35 Shore A
- excellently scannable