

**NEW: 60 Shore – now available in white!**

## **bioplast white**



Addition curing silicone for indirect fitting (in the lab) of soft-flexible earmoulds, eg. BTE and noise protection. If required, the earmould can be finished with a special silicone lacquer, e.g. supercoat nano, for permanently antibacterial surface sealing.

### **Features:**

- highest wearing comfort due to soft fitting
- robust and long-term durable
- extremely tear resistant
- final hardness: high elastic 40 or 60 Shore A
- easy processing (grinding, cutting)
- skin-friendly and biocompatible