

# **protect mold**

*colourful protection against noise & water*



The new hybrid silicone for versatile uses eg. for the fabrication of individual noise and water protection. Earmoulds made from **protect mold** are elastic, sturdy and floatable. Short setting time, 1:1 mixing, it can also be used with the indirect method, in the lab.

## **Advantages:**

- direct fitting, previous impression taking is not required
- long-term durable, tear resistant and floatable
- ideal final hardness: approx. 40 Shore A
- new 80 ml dosing syringes, easy processing
- fluorescence colours can be mixed individually

