

Super fast set, maximum hardness & fresh apple flavour

the cutting edge



greenbite *apple*

bite registration material, A-silicone based

- Super fast set: extremely short time in the mouth, **only 60 seconds.**
- **Maximum final hardness**, approx. 45 Shore D, allows easy processing (cut and shape).
- **12 contouring tips free** in each standard packing for precise and easy application directly onto the teeth. Thixotropic, does not flow into interdental spaces.
- With **apple flavour** for stress-free bite taking

GRATIS
hand & nail cream
"green apple"
with the purchase
of each standard pack
until stocks last



EXTRA:
12 contouring tips

