

start

TRIMMING



no
time delay

without DELAY

Immediately after removal only tempofit® temporaries have the necessary hardness for mechanical processing (cutting, grinding, polishing) by a new bis-composite.

NEW tempofit®

Material advantage

Self-curing material for the fabrication of crowns and bridges (2:1), for high-quality temporaries



- colour-stabilized for long-term natural look; polishable to a high sheen.
- minimal shrinkage guarantees accurate occlusion.

- time-saving features: fast setting material and **smartrepair®** UV-light-polymerizing-repair kit.
- patient-friendly features: low-temperature polymerization, highly resistant to pressure and abrasive damage, tasteless and tissue-compatible.

System advantage

complete system of intercompatible components and AUTOMIX2-dispenser gun



- **smartrepair®** set for quick and precise repairs, features UV-light-polymerizable material, economical applicators and exchangeable tips.



- one size fits all: **tempofit®** 2:1 cartridges are compatible with the new standard 1:1 dispenser guns, no additional dispensing tool required.

Yes, I want to know more about the advantages of tempofit®.
Information-Phone +49 / 72 43 / 5 10 - 146
e-mail info@detax.de

+++ DETAX +++
+++ press release +++
+++ new developments +++

tempofit®

Complete system for making and repairing temporary crowns and bridges

tempofit®, the complete new system for aesthetically pleasing crown and bridge temporaries, is now available as a standard package or as our Introkit, complete with a novel repair set (**smartrepair®** + **smartbond®**).

tempofit® is long-term color-stabilized, self-curing, quick-setting and polymerizes at comfortable temperatures. **tempofit®** temporaries provide an accurate fit, resist pressure or abrasive damage, are tasteless and biocompatible. The new Automix2 dispenser system prevents cross-contamination, minimizes material waste and allows easier extrusion.

smartrepair® features improved applicators with exchangeable tips that allow exact application. This makes even complicated repairs, additions and contouring easy and quick to complete. **smartrepair®** is light curing, saving the user time as well as material.

In addition its excellent material properties, **tempofit®** in the new 2:1 cartridge offers a clear practical advantage through its system of intercompatible components and the Automix2-system. **tempofit®** 2:1 cartridges can be used with the new, standard 1:1 dispenser guns and no longer require a separate dispenser.



DETAX

Postbox 10 02 25 76256 Ettlingen / Germany
Phone +49 / 72 43 / 5 10-0
Fax +49 / 72 43 / 5 10-100

Chemistry for health