

smartprotect® soft

100% risk free



Extending its prevention range, DETAX introduces **smartprotect® soft** – new generation of desensitizer. **smartprotect® soft** is based on oxalate and designed for treatment of dental hypersensitivities, for application on enamel and exposed cervical areas. Its mineralic active substances provide effective sealing of the dentinal tubules. The increased sensitivity to a stimulus (mechanical, chemical or osmotic) is effectively reduced, the sensitivities decreased perceptibly. Next, natriumfluoride further supports remineralisation.

smartprotect® soft is tasteless, free of alcohol, solvents, HEMA and glutaraldehyde. Applying the solution does cause neither decoloration of teeth nor any irritation of the gingiva.

The standard packing contains a patients kit for completed treatment at home, holiday etc.

DETAX GmbH & Co. KG
www.detax.de