

tempofit® N colours ***small but powerful***



tempofit®, the complete system for aesthetically pleasing crown and bridge temporaries has been extended by an additional innovative application form:

tempofit® N colours bis-acryl-composite (4:1) self curing is easy and quickly to use. Due to modified nano particles, temporaries stand up excellently to functional stress and are highly resistant to pressure or abrasive damage. **tempofit® N colours** cures quickly in the mouth (approx. 2 min) and can thus be removed quickly, trimmed, finished and polished at a minimal time cost. Brilliant and translucent shades ensure a natural aesthetic.

Small quantities for rarely used shades: the new 16gr-mini-mix-cartridges are ideal for small and single temporaries. Particular mixing cannulas are minimizing waste of material; no mixing dispenser is required.

tempofit® N colours is available in single packing with 1 mini-mix cartridge (plus mixing cannulas); shades A2, B1, B3 and BL.